


User Manual

3900 ControlLog

By Thunder Scientific Corporation



Contents

Getting Started	1
About ControLog	1
Notice.....	1
Technical Support.....	1
License Agreement	2
Warranty	3
Copyright	3
Trademark.....	3
ControLog Minimum System Requirements	4
Installing ControLog.....	5
Data Backup	5
Interface Connections	6
Starting ControLog	7
Crashes and Data Recovery	12
 ControLog Interface	 15
Menu Bar	16
File Menu	17
Mode Menu	21
Run Menu	23
Profile Menu.....	25
Graph Menu.....	27
Data Menu	30
Units Menu	35
Connections Menu.....	38
Help Menu.....	42
Parameters Tab Group	43
3900 Parameter Tab.....	49
Device Parameter Tabs.....	57
Data and Graph Tab Group.....	59
Data Tabs.....	61
Graph Tabs	61
Profile Tab.....	61
Information Tab Group.....	62
Status Log.....	66
3900 Reported Errors	66
System timing	66
3900 Interface Console.....	67
Device Interface Console	68
Status Bar.....	70
Connection and Run Status	70
Auto Profile Controls and Status	72
Current Point Time	73
 Graphing	 75

Creating a New Graph	76
Customizing a Graph	80
Pan	80
Zoom	81
Zoom Graph's X Axis	82
Zoom Graph's Y Axis	83
Auto Scale	84
Scale	84
Scale X Axis	84
Scale Y Axis	84
Graph Properties	85
Saving a Graph	88
Printing a Graph	90
 Data and Data Summary	 91
Device Data Tabs	94
File Data Tabs	95
Data Summary Tabs	98
Creating a Data Summary	99
Recorded Data Points Tab	107
How to Record a Data Point	108
 Auto Profiling	 113
Creating a New Profile	114
Saving a Profile	124
Opening a Profile	125
Printing a Profile	126
Running an Auto Profile	127
Understanding Profile Phases	129
Manual Override of Profile	131
 Connections	 133
Serial Connection	134
Serial Connection Example 1	147
Serial Connection Example 2	158
Analog Connection	174
Analog Connection Example	182
GPIB Connection	190
GPIB Connection Example	203
Manual Connection	212
Manual Connection Example	216
 Glossary of Terms	 225
 Index	 231

Getting Started

This section will provide the user with general information on the 3900 ControLog® application, information on where to obtain technical support, the required hardware specifications to run ControLog, instructions on how to install ControLog and instructions on how to start ControLog for the first time. Following sections will provide further details on how to use and operate ControLog.

About ControLog

ControLog is a Windows based software application that fully automates the operation of a Thunder Scientific 3900 Humidity Generator and allows various device connections through a number of different interfaces. Data from the generator and attached devices is automatically retrieved and stored for viewing in either numerical or graphical format in real time or post process.

Notice

The information contained in this document is subject to change without notice.

Thunder Scientific Corporation makes no warranties, either express or implied, regarding the keystroke procedures and program, or the fitness of these procedures or program for a particular purpose. The keystroke procedures are made available solely on an "as is" basis and the entire risk as to their quality and performance rests with the user. Thunder Scientific Corporation shall not be liable for any incidental or consequential damages in connection with or arising out of the furnishing, use, or performance of the keystroke procedures or program.

Technical Support

If the user requires assistance with any aspect of the ControLog application, technical support can be obtained by contacting Thunder Scientific Corporation by any of the following means:

www.ThunderScientific.com

Tel: 505.265.8701

FAX: 505.266.6203

support@thunderscientific.com

License Agreement

THIS IS A LEGAL AGREEMENT BETWEEN YOU, THE END USER, AND THUNDER SCIENTIFIC CORPORATION ("TSC"). THE ENCLOSED THUNDER SCIENTIFIC CORPORATION SOFTWARE PROGRAM (THE "SOFTWARE ") IS LICENSED BY THUNDER SCIENTIFIC CORPORATION TO THE ORIGINAL CUSTOMER AND ANY SUBSEQUENT TRANSFEREE OF THE PRODUCT FOR USE ONLY ON THE TERMS SET FORTH HERE. PLEASE READ THIS LICENSE AGREEMENT. IF YOU DO NOT AGREE TO THESE TERMS, RETURN THE FULL PRODUCT WITH PROOF OF PURCHASE WITHIN 30 DAYS FOR A FULL REFUND.

- **GRANT OF LICENSE.** TSC grants to you the right to use one copy of the enclosed SOFTWARE on a single terminal connected to a single computer (i.e. single CPU) or to a network server. If you install the SOFTWARE on a network server, you must purchase a separate copy of the SOFTWARE for each computer terminal that will be used to operate the SOFTWARE. If the anticipated number of users of the SOFTWARE will exceed the number of applicable Licenses, then you must have a reasonable mechanism in place to ensure that the number of persons using the SOFTWARE concurrently does not exceed the number of Licenses. If the SOFTWARE is permanently installed on the hard disk or other storage device of a computer (other than a network server) and one person uses that computer more than 80% of the time, then that person may also use the SOFTWARE on a portable or home computer.

- **COPYRIGHT.** The SOFTWARE is owned by TSC and is protected by United States copyright laws and international treaty provisions. You may either (a) make two copies of the SOFTWARE solely for backup or archival purposes provided that you reproduce all copyright and other proprietary notices that are on the original copy of the SOFTWARE provided to you, or (b) transfer the SOFTWARE to a single hard disk provided you keep the original solely for backup or archival purpose. You may not copy the written materials accompanying the SOFTWARE.

- **OTHER RESTRICTIONS.** You may not rent or lease the SOFTWARE, but you may transfer the SOFTWARE and accompanying written materials on a permanent basis provided you retain no copies and the recipient agrees to the terms of this Agreement. You may not reverse engineer, decompile, disassemble, or create derivative works from the SOFTWARE.

SOFTWARE and written materials

- **GOVERNMENT LICENSEE.** If you are acquiring the SOFTWARE on behalf of any unit or agency of the United States Government, the following provisions apply:

The Government acknowledges TSC's representation that the SOFTWARE and its documentation were developed at private expense and no part of them is in the public domain.

The Government acknowledges TSC's representation that the SOFTWARE is "Restricted Computer Software" as that term is defined in Clause 52.227-19 of the Federal Acquisition Regulations (FAR) and is "Commercial Computer Software" as that term is defined in Subpart 227.471 of the Department of Defense Federal Acquisition Regulations Supplement (DFARS).

The Government agrees that:

- (i) If the SOFTWARE is supplied to the Department of Defense (DoD), the SOFTWARE is classified as "Commercial Computer Software" and the Government is acquiring only "restricted rights" in the SOFTWARE and its documentation as that term is defined in Clause 252.227-7013(c) (1) of the DFARS, and
- (ii) If the SOFTWARE is supplied to any unit or agency of the United States Government other than DoD, the Government's rights in the SOFTWARE and its documentation will be as defined in Clause 52.227-19 (c) (2) of the FAR.

- **RESTRICTED RIGHTS LEGEND.** Use, duplication, or disclosure by the Government is subject to restrictions as set forth in subparagraph (c) (1) (ii) of the Rights in Technical Data and Computer Software clause at DFARS 252.227-7013. Thunder Scientific Corporation, 623 Wyoming SE, Albuquerque, NM 87123.

- **EXPORT LAW ASSURANCES.** You acknowledge and agree that the SOFTWARE is subject to restrictions and controls imposed by the United States Export Administration Act (the "Act") and the regulations thereunder. You agree and certify that neither the SOFTWARE nor any direct product thereof is being or will be acquired, shipped, transferred

or reexported, directly or indirectly, into any country prohibited by the Act and the regulations thereunder or will be used for any purpose prohibited by the same.

- **GENERAL.** This Agreement will be governed by the laws of the State of New Mexico, except for that body of law dealing with conflicts of law.

Should you have any questions concerning this Agreement, or if you desire to contact TSC for any reason, please write: Thunder Scientific Corporation, 623 Wyoming NE, Albuquerque, NM 87123

Warranty

Thunder Scientific Corporation warrants that the media on which the Software is distributed will be free from defects for a period of twelve (12) months from date of original purchase.

During the warranty period, we will at no charge, repair or replace, at our option, any product that proves to be defective, provided you return the product to Thunder Scientific Corporation. This warranty does not apply to any product that has been repaired, serviced, altered, or modified by anyone other than an authorized representative of Thunder Scientific Corporation.

Thunder Scientific Corporation shall not be liable for any consequential, special, incidental, direct or other damage, loss, or expense arising out of or connected with the manufacture, programming, packaging, storage, delivery (including any delay in or failure to deliver), or use of any products, whether or not resulting from negligence or breach of contract, except to the extent that those items sold and warranted by Thunder Scientific Corporation shall be repaired or replaced in accordance with the previously stated terms of this warranty contract. The repair or replacement of a product is your exclusive remedy.

THUNDER SCIENTIFIC CORPORATION MAKES NO OTHER WARRANTY OR REPRESENTATION OF ANY KIND, EXPRESS OR IMPLIED, AND NO WARRANTY SHALL BE IMPLIED BY LAW.

Copyright

©2009-2011 Thunder Scientific Corporation

623 Wyoming Blvd. SE

Albuquerque, NM 87123, USA

This document contains proprietary information which is protected by copyright. All rights are reserved. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording, for any purpose without the prior written consent of Thunder Scientific Corporation.

Trademark

ControLog is a trademark of Thunder Scientific Corporation.

HumiCalc is a trademark of Thunder Scientific Corporation.

ControLog Minimum System Requirements

The following specifications are the required PC minimum system requirements to run ControLog.

- 1.6 GHz or greater Intel® Pentium® or equivalent processor (2.66 GHz Multi Core or better is recommended when connecting multiple devices)
- 512 MB or greater of RAM (4 GB or more is recommended when connecting multiple devices)
- Minimum 800 x 600 screen resolution
- Microsoft® Windows 7 (x86 or x64), Microsoft Windows Vista® SP2 (x86 or x64), Microsoft Windows XP® SP3 (x86 only)
- Microsoft .NET Framework version 4.0
- Internet Browser

If the user desires to connect multiple RS-232 devices or RS-485 devices to ControLog, then the following interface is recommended:

- MOXA® UPort® 1450 USB-to-serial converter

If the user desires to connect a GPIB device to ControLog, then the following interface is required:

- National Instruments® GPIB-USB-HS USB-to-488.2 converter


If the user desires to connect Analog devices to ControLog, then the following data acquisition equipment is required:

- Agilent® 34970A Data Acquisition/Switch Unit
- Agilent® 34901A 20-Channel Multiplexer


Installing ControLog

The following instructions give a step by step process to install the ControLog application on a PC.

1. Insert the ControLog CD or locate the downloaded ControLog installation package. The ControLog installation package will consist of the following files:



2. Click the setup.exe file to begin the installation.



Note: The user must have administrative rights to install the application

3. Follow the on screen installer directions.
4. If using a MOXA® UPort 1450 USB-to-serial converter then install the appropriate drivers according to the manufactures instructions.
5. If using a National Instruments® GPIB-USB-HS USB-to-488.2 converter then install the appropriate drivers according to the manufacturer's instructions.

Warning: PC anti-virus software can interfere with ControLog operation. It is recommended to avoid anti-virus scanning during ControLog operation.

Data Backup

ControLog does not require any direct backup of its operating files. If for any reason a need to restore the software occurs, simply reinstall the application as described in the [Installing ControLog](#) section.

It is recommended to perform a periodic backup on any user generated files such as uncertainty solutions for the generator, profiles, device setups and/or data files.

The default location for user data is under the follow directory: ...Documents\Thunder Scientific\3900 ControLog\

Warning: The user has full control as to where to save data and the above directory may not be the location they choose.

Interface Connections

The following instructions provide details to connect the Thunder Scientific Model 3900, MOXA® UPort® USB-to-serial converter, National Instruments® 488.2 and the Agilent® 34970A to the ControLog PC.

- **Connecting to the 3900**

ControLog communicates with the 3900 through an RS-232 communication connection on any available COM port. To connect the 3900 to the ControLog PC use a 9 pin cable and connect the computer's 9 pin male serial port connector to the 9 pin female console port connector of the 3900.

- **Connecting the MOXA® UPort® 1450 USB-to-serial converter**

To connect the USB-to-serial converter, plug the USB cable that is supplied with the converter into an open USB port on the ControLog PC. Please refer to the manufacturer's documentation for further information.

- **Connecting the National Instruments® GPIB-USB-HS USB-to-488.2 converter**

To connect the National Instruments® USB-to-488.2 converter, plug the USB end of the converter into an open USB port on the ControLog PC. Please refer to the manufacturer's documentation for further information.

- **Connecting to the Agilent® 34970A**

ControLog can communicate with an Agilent® 34970A to record analog devices through an RS-232 communication connection on any available COM port. To connect a 34970A to the ControLog PC use a 9 pin cable and connect the computer's 9 pin male serial port connector to the 9 pin female on the back of the 34970A. Please refer to the manufacturer's documentation for further information and on how to connect and select the RS-232 interface for communication.


Starting ControLog

This section is intended to quickly familiarize the user with the operation of starting ControLog.


1. Ensure that the computer's serial port is properly connected to the 3900 and all desired interfaces are connected.
2. Locate and double-click the ControLog icon.




3. The splash screen will appear once the program is initialized. Note: It may take some time for the program to initialize the first time after installing the application.
4. The splash screen will display a loading status near the bottom, once loading is complete the splash screen will close and the main window will be displayed.
5. A ControLog Product Key Dialog will appear the first time ControLog starts if no valid license is present.



6. Enter your "Name", "Organization" and the "Product Key" from the back of the ControLog CD case or as received in the email receipt when purchasing ControLog online.
7. Click Ok to save the license information.
8. Upon completion of loading, ControLog will automatically attempt to establish communication with the 3900 using the last known or default communication settings.



9. If ControLog is unable to establish communication with the 3900 generator a message will be displayed to inform the user. The message will also ask the user if they would like to check the communication settings.



10. Selecting yes, will open the 3900 communication settings.

3900 Settings

Select the port the 3900 is connected to

Connect using **COM1**

Select the port settings for the 3900

Baud Rate **2400**

Data Bits **8**

Parity **None**

Stop Bits **1**

Handshake **None**

☒ Enable RTS ☐ Enable DTR

< Back **Next >** **Finish** **Cancel**

11. In the 3900 settings, select the COM port that is connected to the 3900 along with the port settings for the 3900. Refer to your 3900 manual to check the current port settings for your 3900.

Default 3900 Port Settings:

Baud Rate: **2400**

Data Bits: **8**

Parity: **None**

Stop Bits: **1**

Handshake: **None**


Enable RTS: **Yes**

Enable DTR: **No**


12. Selecting the “Next” button will display the access interval for the 3900. The access interval is the rate at which ControLog sends and receives commands to and from the 3900. It is recommended to leave the access interval at the default 1.5 seconds.



13. Selecting the “Finish” button will result in ControLog attempting to establish communication with the 3900 using the new communication settings.




14. Once Communication is established, the 3900 parameter tab will appear along with a 3900 data tab. The status log will also show that communication has been established.




Crashes and Data Recovery

If for any reason the PC crashes or power is interrupted, ControLog will automatically load all unsaved data on the next successful startup. ControLog actively stores all device data to Tab delimited text files using a *.backup file extension during normal operation. These files reside in a “Backup” folder and upon a successful closing of the application this folder is removed. If on startup the “Backup” folder exists, then all residing backup files will be automatically loaded to allow the user to save or append the data.


ControLog will display a message telling the user that data has been recovered.



If the generator is running, ControLog will ask the user if they would like to append the data for the generator. If “Yes” is selected, all new data recordings will be appended to the current data. If “No” is selected, ControLog will prompt the user to save the data before clearing it.



If a device is connected, ControLog will ask the user if they would like to append the data for the device. If “Yes” is selected, all new data recordings will be appended to the current data. If “No” is selected, ControLog will prompt the user to save the data before clearing it.



Individual Data Tabs for each backup file will be automatically loaded and the status log will record the event.

The screenshot displays the 3900 ControLog software interface. The main window is titled "3900 ControLog" and features a menu bar with options: File, Mode, Run, Profile, Graph, Data, Units, Connections, and Help. The interface is divided into several sections:

- 3900 Parameters:** A table with columns for Setpoint, Actual, and U (Uncertainty). It lists various parameters such as Frost Point, Saturation Pressure, Saturation Temperature, Test Pressure, Test Temperature, Mass Flow Rate, %RH, Dew Point, PPMv, PPMw, Grains/lb, Enthalpy, SVP@Tt, SVP@Td, and SVP@Ts.
- 3900 Data:** A table showing data for the date 11/9/2009. It includes columns for Date/Time, Setpoint [kPa], and Saturation Pressure [kPa]. A yellow callout box points to the "Recovered Data is load into a Data Tab" message.
- Status Log:** A log of system events. A yellow callout box points to the message "Status Log shows the recovered data has been loaded".
- 3900 Reported Errors:** A section at the bottom showing "3900 Reported Flow Valve Not Closing".


The Status Log contains the following entries:

- 11/9/2009 9:06:23 AM: The system has recovered from a crash and the recovered data has been load
- 11/9/2009 9:06:38 AM: Communication established with the 3900 generator
- 11/9/2009 9:06:38 AM: System i
- 11/9/2009 9:06:41 AM: 3900 Te
- 11/9/2009 9:06:43 AM: 3900 Te
- 11/9/2009 9:06:44 AM: 3900 Reported Expansion Valve Not Closing
- 11/9/2009 9:06:44 AM: 3900 Reported Flow Valve Not Closing

ControLog Interface


This section will provide the user with a detailed overview of ControLog's layout and design. It is intended to allow the user to gain familiarity with ControLog's user interface. The sections following will provide a deeper operational view of the functionality that ControLog offers. However, fundamental operation of the Model 3900 is assumed, and may be found in the Operation and Maintenance Manual provided with the 3900 generator.

The ControLog program is composed of five basic features: Menu Bar, Parameters Tab Group, Data and Graph Tab Group, Information Tab Group and the Status Bar. Each feature is designed to be intuitive to use and to provide the user with detail information on the operation of the generator and/or connected device(s).




Menu Bar

The **Menu Bar** is located at the top of the application and contains various dropdown menus, which provide access to the different ControLog functions and controls. The Profile, Graph, Data menus are dynamic and can change depending on which tab is selected. For example, the Data menu has operations that are specific to the selected data tab and will hide those specific operations when another non-data tab is selected.



File Menu

The **File Menu** allows the user to perform file specific commands. It allows the user to open previous data files, save all open data, close all open data or graph tabs, clear the status log, save the status log, print the status log, import uncertainty solutions and exit the application.



Open Data File

The **Open** file menu command allows the user to open previous data files for further review and analysis. Selecting this command will open a file dialog that will allow the user to browse to the desired location for the file to open. ControLog can open data saved in the follow type and format:

- Text File (Comma Delimited) (*.csv;*.txt)
- Text File (Tab Delimited) (*.dat;*.txt)
- Excel Workbook (*.xlsx;*.xls)
- Backup ControLog File (*.backup)

Save All Data

The **Save All** file menu command allows the user to save all current open data tabs to individual files using a common name. This feature allows the user to quickly save multiple data tabs in a single operation. Selecting this command will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. ControLog can save data in the following type and format:

- Comma-Separated Values (*.csv)
- Text File (Comma Delimited) (*.txt)
- Text File (Tab Delimited) (*.txt)
- Excel Workbook (*.xlsx)
- Excel 97-2003 Workbook (*.xls)

Example: If the user had two data tabs open, one called “3900 Data” and the other called “Agilent 34970A Data” and the user wanted to save the files as Excel Workbooks using the name “Test Data 10Aug09”. ControLog would save two files to the user specified location with the follow names:

Test Data 10Aug09 (3900).xls
Test Data 10Aug09 (Agilent 34970A).xls

Close All Data Tabs

The **Close All** Tabs file menu command allows the user to close all open data and graph tabs. ControLog will ask the user to confirm before closing the tabs and any unsaved data tab will also prompt the user to save the data before the tab is closed.

Note: The 3900 must be shutdown or disconnected before this command becomes available.

Clear Status Log

The **Clear Status Log** file menu command allows the user to clear all current entries in the status log. The user will be asked to save the status log data before the log is cleared.

Save Status Log

The **Save Status Log** file menu command allows the user to save the current entries in the status log. Selecting this command will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. All status log files are saved in HTML format (*.html).

Print Status Log

The **Print Status Log** file menu command allows the user to print the current entries in the status log. Selecting this command will open a print dialog that will allow the user to select the desired print options.

Import Uncertainty Solution



The **Import Uncertainty Solution** file menu command allows the user to import a HumiCalc with Uncertainty solution into ControLog to define the uncertainty for the 3900 generator. Selecting this command will open an “Import Uncertainty” dialog that will step the user through a two step import process. The 3900 uses two pressure ranges, one pressure transducer for saturation pressures below 50 psia and one pressure transducer for saturation pressures above 50 psia. These pressure transducers have different uncertainties and therefore require different HumiCalc solutions. Refer to your HumiCalc with Uncertainty reference manual for more information on creating uncertainty solutions.

Import Uncertainty

Below 50 psia range

Select the name and location of the HumiCalc solution file for the given generator pressure range.

The first step in the import process is to select the uncertainty solution for the 3900 when it is operating below 50 psia. Clicking the “Browse” button will open a file dialog that will allow the user to browse to the desired location for the solution file to open. Once a file has been selected, clicking the “Next” button will display a status dialog as the first solution is imported into ControLog.




The second step in the import process is to select the uncertainty solution for the 3900 when it is operating above 50 psia. Clicking the “Browse” button will open a file dialog that will allow the user to browse to the desired location for the solution file to open. Once a file has been selected, clicking the “Finish” button will complete the import process and will display a status dialog as the second solution is imported into ControLog.

Exit ControLog

The **Exit** file menu command allows the user to exit the application. If there is current data that has yet to be saved the user will be asked to save the data before the application closes.

Mode Menu

The **Mode Menu** allows the user to change the operating mode of the 3900. It allows the user to select between Frost Point, Dew Point, PPMv, PPMw, Percent Relative Humidity and Saturation Pressure. The current mode is indicated by a checkmark.



Saturation Pressure Control Mode

The Saturation Pressure Control Mode, Ps, is controlled at a constant value independent of any other pressure, temperature, or humidity value. While saturation pressure is held constant, all humidity parameters may vary.

Percent Relative Humidity Control Mode

The Percent Relative Humidity Control Mode or %RH Control Mode, %RH, is controlled at a constant value by varying saturation pressure, Ps, to compensate for any changes in saturation temperature, Ts, test temperature, Tt, or test pressure, Pt. While %RH is held constant, all other humidity parameters may vary. While in %RH control mode, the saturation temperature setpoint is automatically determined.

Frost Point Control Mode

The Frost Point Control Mode, Tf, is controlled at a constant value by varying the saturation pressure, Ps, to compensate for changes in either saturation temperature, Ts, or test pressure, Pt. While Frost Point is held constant other humidity parameters may vary. While in Frost Point control mode, the saturation temperature setpoint is automatically determined. Frost Point is independent of test temperature.

Dew Point Control Mode

The Dew Point Control Mode, Td, is controlled at a constant value by varying saturation pressure, Ps, to compensate for any changes in either saturation temperature, Ts, or test pressure, Pt. While Dew Point is held constant, other humidity parameters may vary. While in Dew Point control mode, the saturation temperature setpoint is automatically determined. Dew Point control mode is valid both above and below 0 °C, and Dew Point is independent of test temperature.

PPMv Control Mode


The PPMv Control Mode is controlled at a constant value by varying saturation pressure, Ps, to compensate for any changes in saturation temperature, Ts. While PPMv is held constant, other humidity parameters may vary. While in PPMv control mode, the saturation temperature setpoint is automatically determined. PPMv is independent of test pressure and test temperature.

PPMw Control Mode

The PPMw Control Mode is controlled at a constant value by varying saturation pressure, P_s , to compensate for any changes in saturation temperature, T_s . While PPMw is held constant, other humidity parameters may vary. While in PPMw control mode, the saturation temperature setpoint is automatically determined. PPMw is independent of test pressure and test temperature.

Run Menu

The **Run Menu** allows the user to run the 3900 manually in purge mode or generate mode. The Run menu also allows the user to manually shutdown the 3900. The current run mode is highlighted.



Purge Mode

Selecting **Purge** from the run menu commands the 3900 into purge mode. The purge mode is generally used to prevent icing within the saturator and dry the saturator outlet after movement (transportation), storage (power off, no gas flow, etc.), after performing the saturator fill procedure, or while transitioning the saturator from higher to lower temperatures. When the system is not being used (power off, no gas flow, etc.), the saturator is closed off and the gas within is static. As thermal equilibrium is reached, water vapor will condense on all inner surfaces between the saturator outlet and the expansion valve inlet. The Purge mode counteracts this condition by allowing the carrier gas to flow in the opposite direction (expansion valve to saturator), drying the affected sections of tubing. This is a necessary preparatory step in any low humidity system. As a general rule, when starting from an ambient condition the system should be purged for 24 hours or more before attempting to operate in the Generate mode. If sufficient purge time is not allowed, condensed or trapped water will remain and system accuracy will suffer. Insufficient purge time is usually indicated by higher than normal (wetter than normal) indications of the device under test. These indications can be as little as a few tenths of a degree to as much as several degrees dew/frost point. Purging should also be performed while transitioning from warmer to colder saturation temperatures, and for approximately 5 hours after each 500 hours of continuous Generate mode operation. During Purge mode, both flow control and saturation temperature control are active, but saturation pressure control is disabled. The generator will attempt to achieve the indicated flow and saturation temperature setpoint values.

Note: When the saturation temperature is lowered, the fluid jacket surrounding the saturator cools in order to reduce the saturation temperature to its new setpoint value. As the saturator cools during this transition period, temperature gradients will exist between the inside of the saturator and the fluid jacket that surrounds it. The saturator outlet passes through this fluid jacket and will also exhibit temperature gradients along its length. If gas is allowed to flow normally through the saturator during this cooling period, the 100% humidified gas of the saturator may condense at the colder saturator outlet. Therefore, Purge mode should be used while cooling the saturator to lower temperatures. For this reason, the lowest humidity of a generation sequence or profile should be performed first. This low to high order requires that a Purge be performed only once prior to the generation sequence when cooling to the lowest saturation temperature. Then, as humidity values are increased, warming the saturator to higher values, further purging is not required.

Note: During Purge, no gas flows to the 3900 Conditioned Gas Outlet, and consequently no gas flows through the device under test if connected.

Generate Mode

Selecting **Generate** from the run menu commands the 3900 into generate mode. Generate mode is used to operate the system when exact humidity points or associated time intervals have not been determined, when data must be viewed and/or verified manually before proceeding to the next humidity point, or when more immediate control over the generated humidity is required. When in the Generate mode of operation, the system will control at the currently entered setpoint indefinitely. Any time a setpoint is changed, the system immediately begins adjusting to that new value, and will

control at the new point indefinitely. Generate mode offers the flexibility to change the setpoint at any time, and does not force you into any set sequence or for any prescribed amount of time. The Generate mode also allows you to change the humidity control mode at any time. For instance, the system may be controlling Frost Point, then immediately switched to PPMv control mode.

Shutdown


Selecting **Shutdown** from the run menu commands the 3900 to shutdown. The 3900 may be shutdown while either generating or purging. When stopped, all system functions shutdown, pressure is vented, and the idle Control/Display screen is shown. During this idle time when the 3900 is stopped, gas is not flowing through the saturator. After extended periods of this idle time, the 3900 must be purged again prior to further use.

Profile Menu


The **Profile Menu** allows the user to Open, Save, and create new Auto Profiles. Auto Profiles give the user the ability to program a set of humidity and temperature test points and dwell times that will automate the 3900 generation process. The profile menu is dynamic and has operations that are specific to the profile tab. Specific operations will be hidden when another non-profile tab is selected.

The Auto Profiling feature is very similar to the Generate mode with the main exception that profiling relies on a predefined list of setpoints referred to as a profile. The user configurable profile is used as ControLog's road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.


For more information, see [Profiling](#).



Profile menu when the Profile tab is not selected



Profile Menu when the Profile tab is selected



Profile Menu when running an Auto Profile

New Profile

The **New** profile menu command allows the user to create a new Auto Profile. This operation is always available.

Open Profile

The **Open** profile menu command allows the user to open a saved Auto Profile. Selecting this command will open a file dialog that will allow the user to browse to the desired location for the file to open. This operation is always available.

Close Profile

The **Close** profile menu command allows the user to close the Profile tab. Selecting this command will close the Profile tab, but if the profile has not been saved, ControLog will ask the user to save the profile before closing the Profile tab.

Save Profile

The **Save** profile menu command allows the user to save the currently opened Profile. Selecting this command will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. ControLog Auto Profiles are saved in XML format with a *.profile extension. This operation is available only when the Profile tab is selected.

Print Profile

The **Print** profile menu command allows the user to print the currently opened Profile. Selecting this command will open a print dialog that will allow the user to select the desired print options.

Run Auto Profile

The **Run Auto Profile** command allows the user to start an Auto Profile. Selecting this command will open the Profile Starting Point dialog which allows the user to select which point in the profile they would like to start the profile on.

Hold Profile Point

The **Hold Point** profile menu command allows the user to hold or pause the current Auto Profile point. Selecting Hold Point pauses the current point, allowing the system to remain indefinitely at the current point. While in a hold mode, the system is prevented from completing the ramp, assurance, or soak phases for a point. To resume the profile point, select the menu item again. This re-enables the point and allows the profile to resume normal operation.

Next Profile Point

The **Next Point** profile menu command allows the user to skip to the next point in the Auto Profile. Selecting Next Point manually advances to the next point, skipping any remaining ramp, assurance, or soak phase.

Next Profile Phase

The **Next Phase** profile menu command allows the user to skip to the next phase in the Auto Profile. Selecting Next Phase manually advances to the next phase. It causes Ramp Phase to proceed to the Assurance or Soak Phase, Assurance to proceed to Soak, or Soak to proceed to Ramp of the next profile point. This allows for early manual termination of any phase within a profile point.


Stop Auto Profile

The **Stop Auto Profile** command allows the user to stop the Auto Profile where it is at. The setpoints will remain where the profile stopped and the generator will continue in its current mode of operation.


Graph Menu

The **Graph Menu** allows the user to create a New graph, Close the selected graph, Close All graphs, Save the selected graph, Print the selected graph, Zoom, Scale and set the Properties for the selected graph. The graph menu is dynamic and has operations that are specific to the selected graph tab. Specific operations will be hidden when another non-graph tab is selected.

For more information, see [Graphing](#).



Graph menu when another non-graph tab is selected



Graph menu when a graph tab is selected.

New Graph

The **New** graph menu command allows the user to create a new graph. Selecting this command will open a New Graph Wizard dialog that will step the user through the selection process of what data the user would like to include in the new graph. This operation is always available.

Close Graph

The **Close** graph menu command allows the user to close the selected graph. Selection will result in a confirmation message to assure the user wants to close the graph. This operation is only available when a graph tab is selected.

Close All Graphs

The **Close All** graphs menu command allows the user to close all open graphs. Selection will result in a confirmation message to assure the user wants to close all graphs. This operation is available whenever a graph tab is open.

Save Graph

The **Save** graph menu command allows the user to save the selected graph. Selecting this command will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. ControLog graphs can be saved in following graphic file types:

- Bitmap (*.bmp)
- Graphics Interchange Format (*.gif)
- Joint Photographic Experts Group (*.jpg)
- W3C Portable Network Graphics (*.png)
- EMF Enhanced Metafile Format (*.emf)

This operation is available only when a graph tab is selected.

Print Graph

The **Print** graph menu command allows the user to print the selected graph. Selecting this command will open a print dialog that will allow the user to select the desired print options. This operation is available only when a Graph tab is selected.

Zoom Graph

The **Zoom** graph menu command allows the user to zoom a rectangular area of the graph. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to create a rectangular area on the graph that will be zoomed. This operation is available only when a Graph tab is selected.

Zoom Graph's X Axis

The **Zoom X Axis** graph menu command allows the user to zoom along the X Axis of the graph. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to create a section area on the graph that will be zoomed along the X Axis. This operation is available only when a Graph tab is selected.

Zoom Graph's Y Axis

The **Zoom Y Axis** graph menu command allows the user to zoom along the Y Axis of the graph. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to create a section area on the graph that will be zoomed along the Y Axis. This operation is available only when a Graph tab is selected.

Auto Scale Graph

The **Auto Scale** graph menu command allows the user to reset the graph view to encompass all data. Selecting this command will automatically reset both axis of the graph so that the entire data set of each item is contained within the boundaries of the graph. This operation is available only when a Graph tab is selected.

Scale Graph

The **Scale** graph menu command allows the user to scale both the X and Y axis. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to scale. Dragging the cursor up scales the

display in (zooms in) and dragging the cursor down scales the display out (zoom out). This operation is available only when a Graph tab is selected.

Scale Graph's X Axis

The **Scale X Axis** graph menu command allows the user to scale the X axis. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to scale. Dragging the cursor up scales the X Axis in (zooms X Axis in) and dragging the cursor down scales the X Axis out (zooms X Axis out). This operation is available only when a Graph tab is selected.

Scale Graph's Y Axis

The **Scale Y Axis** graph menu command allows the user to scale the Y axis. Selecting this command checks the operation in the menu and allows the user to use the left mouse button to scale. Dragging the cursor up scales the Y Axis in (zooms Y Axis in) and dragging the cursor down scales the Y Axis out (zooms Y Axis out). This operation is available only when a Graph tab is selected.


Graph Properties

The graph **Properties** graph menu command allows the user to modify the properties of the selected graph. Selecting this command opens the Graph Properties dialog that allows the user to make changes to what data is graphed, the display properties for each line and the axis values. This operation is available only when a Graph tab is selected.


Data Menu

The **Data Menu** allows the user to Open, Create a Data Summary, Record Data Point, Close, Clear, Save, Print, Copy and change the Data Storage Interval. The data menu is dynamic and has operations that are specific to the selected data tab. Specific operations will be hidden when another non-data tab is selected.


For more information, see [Data and Data Summary](#)



Data menu when another non-data tab is selected



Data menu when a data tab is selected and a device is connected



Data menu when a data tab is selected and no device is connected

Open Data File

The **Open** data menu command allows the user to open previous data files for further review and analysis. This command has the same functionality as the Open command under the file menu. Selecting this command will open a file dialog that will allow the user to browse to the desired location for the file to open. ControLog can open data saved in the following types and formats:

- Text File (Comma Delimited) (*.csv;*.txt)
- Text File (Tab Delimited) (*.dat;*.txt)
- Excel Workbook (*.xlsx;*.xls)
- Backup ControLog File (*.backup)


This operation is always available.

Create Data Summary

The **Create Data Summary** data menu command allows the user to create a summary of any currently opened data. The feature lets the user specify what items from which device and at what intervals to include in the data summary. The data summary can also calculate error between the specified standard and the device under test. This operation is available whenever there is an open Data tab that contains data.


Record Data Point

The **Record Data Point** data menu command allows the user to record certain data items from any currently connected device either manually, with each manual device entry or at the completion of each soak phase in an auto profile. The user can specify the number of prior data points to include and has the option to automatically calculate average and or standard deviation of the prior data points. There are two submenus for this menu command; **Settings** and **Take Point**. This operation is available whenever device or generator data is being logged.




Settings

The **Settings** submenu allows the user to define which data items from which connected device they would like to record when a point is taken. They can also define the number of points prior to include and whether to calculate average and or standard deviation. The user can also configure when to take points, either manually, at the end of a profile's soak phase or when a manual device entry is taken.




Take Point Ctrl-P

The **Take Point** submenu allows the user to take a data point manually. The user can even use the shortcut key “Ctrl-P” to take a point without the need to access the menu.



If the user has not configured the system to record data points then ControlLog will ask the user if they would like to define the settings before taking a point.



For more information on record data point settings, see [Data and Data Summary](#)

Close Data

The **Close** data menu command allows the user to close the selected data tab. ControLog will ask the user to confirm before closing the tab and any unsaved data tab will also prompt the user to save the data before the tab is closed. This operation is available only when a Data tab is selected and the device for the selected data tab is disconnected.

Clear Data

The **Clear** data menu command allows the user to clear the selected data tab. ControLog will ask the user to save any unsaved data tab before the tab is cleared. This operation is available only when the device for the select data tab is connected and there is at least one data point recorded in the selected data tab.

Save Data

The **Save** data menu command allows the user to save the selected data tab. Selecting this command will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. ControLog can save data in the following type and format:

- Comma-Separated Values (*.csv)
- Text File (Comma Delimited) (*.txt)
- Text File (Tab Delimited) (*.txt)
- Excel Workbook (*.xlsx)
- Excel 97-2003 Workbook (*.xls)

This operation is available only when a Data tab with data is selected.

Print Data

The **Print** data menu command allows the user to print the selected data tab. Selecting this command will open a print dialog that will allow the user to select the desired print options. This operation is available only when a Data tab is selected.

Copy All Data

The **Copy All** data menu command allows the user to copy all the data from the selected tab to the clipboard. Selecting this command will copy all data in the selected data tab to the clipboard in a tab delimited format. This data may then be pasted into another program. This operation is available only when a Data tab is selected.


Data Storage Interval

The **Data Storage Interval** data menu command allows the user to change the storage interval that data is recorded at. This is the rate at which data is recorded to the data tabs for all connected devices. Selecting this command will open the “Data Storage Interval” dialog that will allow the user to change the data storage interval during both purge and generate operations. This operation is available only when a device is connected.

Warning: Storing and maintaining data can become a time consuming process. The more data that is stored in the data tab, the slower and less responsive the computer may seem. For this reason, some consideration should be given to the amount of data desired, the overall time span of the data (i.e. hours, days, possibly weeks), and ultimately the data interval.

Note: Data is only recorded while the 3900 is either in generate or purge mode. Data is also stored at the generate rate whenever a device is connected and the 3900 is not connected. This gives the user the ability to

use ControLog as a logging application for any device they can connect without the need of a 3900 generator.




The dialog box is titled "Data Storage Interval" and contains two main sections. The left section is titled "Data storage interval for all devices during generate" and features a numeric input field with the value "5" and a dropdown menu set to "seconds". Below this, a text label reads: "This is the rate at which ControLog will record and graph data from all devices during generate." The right section is titled "Data storage interval for all devices during purge" and has an identical layout with a numeric input field showing "5", a "seconds" dropdown, and the text: "This is the rate at which ControLog will record and graph data from all devices during purge." At the bottom right of the dialog are "Ok" and "Cancel" buttons.

Section	Interval Value	Unit	Description
Data storage interval for all devices during generate	5	seconds	This is the rate at which ControLog will record and graph data from all devices during generate.
Data storage interval for all devices during purge	5	seconds	This is the rate at which ControLog will record and graph data from all devices during purge.

Units Menu

The **Units Menu** allows the user to change ControLog's displayed units. The Temperature, Pressure, Flow Rate, Density and Enthalpy units can be changed.

Note: All parameter tabs and graph tabs will be updated with the selected unit but the data tabs will not change. All data tab values remain in SI units which provide a consistent unit base for saved data.



Temperature Unit

The **Temperature** Unit allows the user to change the displayed units for temperatures. The selected Unit will be checked.



Pressure Unit

The **Pressure** Unit allows the user to change the displayed units for pressure. The selected Unit will be checked.

<input checked="" type="checkbox"/>	psia
<input type="checkbox"/>	atm
<input type="checkbox"/>	Pa
<input type="checkbox"/>	hPa
<input type="checkbox"/>	kPa
<input type="checkbox"/>	MPa
<input type="checkbox"/>	bar
<input type="checkbox"/>	millibar
<input type="checkbox"/>	Torr
<input type="checkbox"/>	in Hg
<input type="checkbox"/>	cm Hg
<input type="checkbox"/>	mm Hg
<input type="checkbox"/>	in H2O
<input type="checkbox"/>	cm H2O
<input type="checkbox"/>	mm H2O

Flow Unit

The **Flow** Unit allows the user to change the displayed units for flow. The selected Unit will be checked.

<input checked="" type="checkbox"/>	l/m
<input type="checkbox"/>	l/h
<input type="checkbox"/>	cfm
<input type="checkbox"/>	cfh

Density Unit

The **Density** Unit allows the user to change the displayed units for density. The selected Unit will be checked.

<input type="checkbox"/>	g/l
<input checked="" type="checkbox"/>	g/m ³
<input type="checkbox"/>	lb/ft ³

Enthalpy Unit


The **Enthalpy** Unit allows the user to change the displayed units for enthalpy. The selected Unit will be checked.

	btu/lb
✓	J/g

Connections Menu

The **Connections Menu** allows the user to create New connections, Open saved connections, Close connections and change the Settings of a connection.

For more information, see [Connections](#)



Connect to 3900

The **Connect to 3900** connection menu command allows the user to connect to the 3900 generator. Selecting this command will cause ControLog to begin establishing communication with the 3900. This operation is available only when the 3900 is disconnected.

New Connection


The **New** connection menu command allows the user to create a new device connection. Selecting this command will open a Connection Wizard dialog that will step the user through the process of creating a new connection to a device.

Open Connection

The **Open** connection menu command allows the user to open previous saved device connections. Device connection files are user saved configurations for a specific device. Selecting this command will open a file dialog that will allow the user to browse to the desired location for the file to open. Device files are saved in an XML format with the (*.device) extension.


Close Connection

The **Close** connection menu command allows the user to close a specific connection. ControLog will ask the user to confirm before closing the connection.



Settings

The **Settings** connection menu command allows the user to change the settings for a specific connection. Selecting this command will open the Connection Wizard dialog for the selected connection, allowing the user to change connection settings as desired.




Open 3900 Interface Console

The **Open 3900 Interface Console** connection menu command under Settings allows the user to open a 3900 console tab. The 3900 console tab allows the user to send and receive commands to and from the 3900.

Warning: Interfacing with the 3900 using the Console should only be attempted if instructed to do so by Thunder Scientific.

Note: This menu command is only available when the 3900 is connected.



Once a valid password is entered the selection will open a 3900 Interface Console tab in the information tab group.




Commands are entered in the bottom text area of the console tab and the 3900 responses are displayed in the upper text area of the console tab.




Open Device Interface Console

The **Open Device Interface Console** connection menu command under Settings allows the user to open a device console tab. The device console tab allows the user to view the commands being sent to and received from any given connected device.


Note: This menu command is only available when at least one device is connected.



The user can select which connected device to view using the drop down selection at the top of the tab.




The data sent to and received from the selected device is displayed in the lower text area of the console tab.



Help Menu

The **Help Menu** allows the user to access the Help Topics, visit Thunder Scientific on the Web and get information About ControLog.



Help Topics

The **Help Topics** help menu command allows the user to open the 3900 ControLog User manual (this document). The user manual is opened in an easy to use help format that allows the user to search and navigate through the entire 3900 ControLog User manual.

Thunder Scientific on the Web

The **Thunder Scientific on the Web** help menu command opens the Thunder Scientific Website using the default internet browser.

About ControLog

The **About ControLog** help menu command opens a dialog giving information on the ControLog application including version number.

Parameters Tab Group

The **Parameter Tab Group** is located on the left side of the application and contains a parameter tab for each connected device. Each parameter tab displays the current data for its particular device. The parameter tabs are a docking style window that can be “pinned” open or allowed to close when not active. A particular parameter tab can be selected by clicking its tab label at the bottom of the group.

By clicking the pin icon on a parameter tab the user can unpin the Parameter tabs.

The screenshot displays the 3900 ControLog application window. The interface includes a menu bar (File, Mode, Run, Profile, Graph, Data, Units, Connections, Help) and a toolbar. The main area is divided into two panes. The left pane, titled '3900 Parameters', lists various parameters with their Setpoint, Actual, and U values. The right pane, titled 'Agilent 34970A Data', displays a table of temperature data. Below the data table is a Status Log pane showing system events. At the bottom, there is a docked area with tabs for '3900 Parameters', 'Agilent 34970A Parameters', 'Status Log', and '3900 Reported Errors'. A status bar at the very bottom indicates 'System is generating' and 'Point Time: 01:32:36'.

Parameter	Setpoint	Actual	U	Units
Frost Point	-10.000	-34.908	±0.099	°C
Saturation Pressure	10.000	63.551	±0.300	psia
Saturation Temperature	-8.0000	-20.730	±0.080	°C
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	21.110	21.110	±0.080	°C
Mass Flow Rate	0.5000	3.9130		l/m
Grains/lb		0.9747		
Enthalpy		21.569		J/g
SVP@Tt		0.3633		psia
SVP@Td		3E-003		psia
SVP@Ts		0.0140		psia
F@Tt,Pt		1.0040		
F@Td,Pt		1.0048		
F@Ts,Ps		1.0189		
Specific Humidity		1E-004		
Absolute Humidity		0.1670		g/m ³

Date/Time Stamp	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]
10/22/2009 11:58:39 AM	26.541	27.158	27.67
10/22/2009 11:58:44 AM	26.569	27.139	27.655
10/22/2009 11:58:49 AM	26.619	27.135	27.655
10/22/2009 11:58:54 AM	26.621	27.132	27.656
10/22/2009 11:58:59 AM	26.593	27.132	27.656
10/22/2009 11:59:04 AM	26.606	27.132	27.657
10/22/2009 11:59:09 AM	26.588	27.124	27.66
10/22/2009 11:59:14 AM	26.606	27.124	27.66

Status Log

- 10/22/2009 9:40:59 AM: System is shutdown
- 10/22/2009 9:41:02 AM: 3900 Test Pressure transducer is disconnected
- 10/22/2009 9:41:04 AM: 3900 Test Temperature probe is disconnected
- 10/22/2009 9:41:05 AM: 3900 Reported Expansion Valve Not Closing
- 10/22/2009 9:41:05 AM: 3900 Reported Flow Valve Not Closing
- 10/22/2009 10:26:41 AM: Generate Mode Requested
- 10/22/2009 10:26:41 AM: Generate Mode Enabled
- 10/22/2009 11:58:39 AM: Communication established with the Agilent 34970A device

3900 Parameters | Agilent 34970A Parameters | Status Log | 3900 Reported Errors

System is generating | Point Time: 01:32:36

Once the parameter tabs are unpinned they will automatically hide.

The screenshot displays the 3900 ControLog software interface. The main window has a menu bar with options: File, Mode, Run, Profile, Graph, Data, Units, Connections, and Help. Below the menu bar, there are tabs for '3900 Data', 'Frost Point', and 'Agilent 34970A Data'. The 'Agilent 34970A Data' tab is active, showing a table of temperature data. The table has columns for 'Date/Time Stamp', 'Temperature 1 [°C]', 'Temperature 2 [°C]', 'Temperature 3 [°C]', 'Temperature 4 [°C]', and 'Temperature 5 [°C]'. The data rows show timestamps from 10/22/2009 11:58:39 AM to 11:59:14 AM with corresponding temperature readings. Below the table is a 'Status Log' section with a list of system events, including shutdowns, disconnections, and mode changes. At the bottom, there is a 'Status Log' button, a warning icon, and the text '3900 Reported Errors'. The bottom status bar shows 'System is generating' and 'Point Time: 01:32:38'.

Date/Time Stamp	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]	Temperature 4 [°C]	Temperature 5 [°C]
10/22/2009 11:58:39 AM	26.541	27.158	27.67	25.376	25.025
10/22/2009 11:58:44 AM	26.569	27.139	27.655	25.35	25.014
10/22/2009 11:58:49 AM	26.619	27.135	27.655	25.334	25.019
10/22/2009 11:58:54 AM	26.621	27.132	27.656	25.328	25.016
10/22/2009 11:58:59 AM	26.593	27.132	27.656	25.328	25.019
10/22/2009 11:59:04 AM	26.606	27.132	27.657	25.333	25.022
10/22/2009 11:59:09 AM	26.588	27.124	27.66	25.353	25.025
10/22/2009 11:59:14 AM	26.606	27.124	27.66	25.349	25.027

Status Log

- 10/22/2009 9:40:59 AM: System is shutdown
- 10/22/2009 9:41:02 AM: 3900 Test Pressure transducer is disconnected
- 10/22/2009 9:41:04 AM: 3900 Test Temperature probe is disconnected
- 10/22/2009 9:41:05 AM: 3900 Reported Expansion Valve Not Closing
- 10/22/2009 9:41:05 AM: 3900 Reported Flow Valve Not Closing
- 10/22/2009 10:26:41 AM: Generate Mode Requested
- 10/22/2009 10:26:41 AM: Generate Mode Enabled
- 10/22/2009 11:58:39 AM: Communication established with the Agilent 34970A device

Status Log ⚠ 3900 Reported Errors

System is generating Point Time: 01:32:38

The user can access the hidden tabs by clicking the desired parameter label on the left.

The screenshot displays the 3900 ControLog software window. The left sidebar shows a list of parameters under the '3900 Parameters' tab. The main window is divided into two sections: a table of parameters with Setpoint, Actual, and U values, and a table of temperatures (Temperature 2 through Temperature 5) with their respective values in °C.

Parameter	Setpoint	Actual	U	Units
Frost Point	-10.000	-34.809	±0.099	°C
Saturation Pressure	10.000	62.859	±0.300	psia
Saturation Temperature	-8.0000	-20.730	±0.080	°C
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	21.110	21.110	±0.080	°C
Mass Flow Rate	0.5000	3.0856		l/m
Grains/lb		0.9853		
Enthalpy		21.573		J/g
SVP@Tt		0.3633		psia
SVP@Td		3E-003		psia
SVP@Ts		0.0140		psia
F@Tt,Pt		1.0040		
F@Td,Pt		1.0048		
F@Ts,Ps		1.0187		
Specific Humidity		1E-004		
Absolute Humidity		0.1688		g/m ³
Dry Air Density		1199.6		g/m ³

Temperature 2 [°C]	Temperature 3 [°C]	Temperature 4 [°C]	Temperature 5 [°C]
27.158	27.67	25.376	25.025
27.139	27.655	25.35	25.014
27.135	27.655	25.334	25.019
27.132	27.656	25.328	25.016
27.132	27.656	25.328	25.019
27.132	27.657	25.333	25.022
27.124	27.66	25.353	25.025
27.124	27.66	25.349	25.027

The bottom status bar indicates 'System is generating' and 'Point Time: 01:32:38'.

Tabs can also be accessed via the drop down selection by clicking on the arrow icon. The currently selected tab is indicated by a check mark in the drop down list.

The screenshot shows the 3900 ControlLog software interface. The main window has a menu bar (File, Mode, Run, Profile, Graph, Data, Units, Connections, Help) and a tab bar at the top with three tabs: 3900 Parameters (selected), 3900 Data, and Frost Point. A dropdown menu is open for the 3900 Parameters tab, showing a checkmark next to '3900 Parameters' and 'Agilent 34970A Parameters' below it. The main display area is divided into three sections: a parameter table on the left, a data table on the right, and a status log at the bottom.

	Setpoint	Actual	U	
Frost Point	-10.000	-34.885	±0.099	°C
Saturation Pressure	10.000	63.387	±0.300	psia
Saturation Temperature	-8.0000	-20.730	±0.080	°C
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	21.110	21.110	±0.080	°C
Mass Flow Rate	0.5000	3.9146		l/m
Grains/lb		0.9772		
Enthalpy		21.570		J/g
SVP@Tt		0.3633		psia
SVP@Td		3E-003		psia
SVP@Ts		0.0140		psia
F@Tt,Pt		1.0040		
F@Td,Pt		1.0048		
F@Ts,Ps		1.0189		
Specific Humidity		1E-004		
Absolute Humidity		0.1675		g/m ³

	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]
10/22/2009 11:58:44 AM	26.541	27.158	27.67
10/22/2009 11:58:49 AM	26.569	27.139	27.655
10/22/2009 11:58:54 AM	26.619	27.135	27.655
10/22/2009 11:58:59 AM	26.621	27.132	27.656
10/22/2009 11:59:04 AM	26.593	27.132	27.656
10/22/2009 11:59:09 AM	26.606	27.132	27.657
10/22/2009 11:59:14 AM	26.588	27.124	27.66
10/22/2009 11:59:19 AM	26.606	27.124	27.66
10/22/2009 11:59:24 AM	26.598	27.126	27.662
10/22/2009 11:59:29 AM	26.598	27.125	27.664
10/22/2009 11:59:34 AM	26.572	27.126	27.668
10/22/2009 11:59:39 AM	26.574	27.128	27.668

Status Log

- 10/22/2009 9:40:59 AM: System is shutdown
- 10/22/2009 9:41:02 AM: 3900 Test Pressure transducer is disconnected
- 10/22/2009 9:41:04 AM: 3900 Test Temperature probe is disconnected
- 10/22/2009 9:41:05 AM: 3900 Reported Expansion Valve Not Closing
- 10/22/2009 9:41:05 AM: 3900 Reported Flow Valve Not Closing
- 10/22/2009 10:26:41 AM: Generate Mode Requested
- 10/22/2009 10:26:41 AM: Generate Mode Enabled
- 10/22/2009 11:58:39 AM: Communication established with the Agilent 34970A device

System is generating Point Time: 01:32:54

Each Parameter tab can also be torn out by dragging the tab label upward. This allows multiple parameter tabs to be visible at once. To place the tab back into the list, drag the tabs label back to the bottom.

The screenshot displays the 3900 ControLog software interface. The main window has a menu bar (File, Mode, Run, Profile, Graph, Data, Units, Connections, Help) and a toolbar. The interface is divided into several sections:

- 3900 Parameters:** A table with columns Setpoint, Actual, U, and units. It lists parameters like Frost Point, Saturation Pressure, Saturation Temperature, Test Pressure, Test Temperature, and Mass Flow Rate.
- Agilent 34970A Parameters:** A section showing actual values for Temperature 1 through Temperature 7.
- Agilent 34970A Data:** A table with columns Date/Time Stamp, Temperature 1 [°C], Temperature 2 [°C], and Temperature 3 [°C]. It shows a series of temperature readings over time.
- Status Log:** A list of system events and errors, including shutdowns, disconnections, and mode changes.
- System is generating:** A status indicator at the bottom left.
- Point Time:** A timestamp at the bottom right showing 01:33:01.

The interface also features a tabbed view at the bottom where the '3900 Parameters' tab is currently selected, and other tabs like 'Agilent 34970A Parameters' and 'Status Log' are visible.

All Parameter Tabs have a context menu that can be displayed by right clicking in the tab. The context menu allows quick access to the device related functions that are available in the Connections menu.

The screenshot displays the 3900 ControlLog software interface. The main window is titled "3900 ControlLog" and features a menu bar with options: File, Mode, Run, Profile, Graph, Data, Units, Connections, and Help. Below the menu bar, there are several tabs: "3900 Parameters", "3900 Data", "Frost Point", "Agilent 34970A Data" (which is the active tab), and "Agilent 34970A backup".

The "3900 Parameters" tab is visible on the left, showing a list of parameters with their Setpoint, Actual, and U values. The "Agilent 34970A Data" tab on the right displays a table of data points. A context menu is open over the "Agilent 34970A Data" tab, showing options: "Connection Settings" and "Close Connection".

The "Agilent 34970A Data" table has the following columns: Date/Time Stamp, Temperature 1 [°C], Temperature 2 [°C], and Temperature 3 [°C]. The data points are as follows:

Date/Time Stamp	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]
10/22/2009 1:25:11 PM	27.154	27.84	28.393
10/22/2009 1:25:16 PM	27.199	27.842	28.391
10/22/2009 1:25:21 PM	27.188	27.841	28.392
10/22/2009 1:25:26 PM	27.081	27.844	28.392
10/22/2009 1:25:31 PM	27.123	27.852	28.39
10/22/2009 1:25:36 PM	27.149	27.857	28.389
10/22/2009 1:25:41 PM	27.136	27.862	28.387
10/22/2009 1:25:46 PM	27.128	27.867	28.386
10/22/2009 1:25:51 PM	27.16	27.873	28.387
10/22/2009 1:25:56 PM	27.191	27.873	28.387
	27.225	27.874	28.387
	27.136	27.875	28.388
	27.086	27.882	28.389
10/22/2009 1:26:16 PM	27.144	27.885	28.389
10/22/2009 1:26:21 PM	27.193	27.885	28.391
10/22/2009 1:26:26 PM	27.24	27.877	28.394

Below the data table is a "Status Log" section with a list of events:

- 10/22/2009 9:41:04 AM: 3900 Test Temperature probe is disconnected
- 10/22/2009 9:41:05 AM: 3900 Reported Expansion Valve Not Closing
- 10/22/2009 9:41:05 AM: 3900 Reported Flow Valve Not Closing
- 10/22/2009 10:26:41 AM: Generate Mode Requested
- 10/22/2009 10:26:41 AM: Generate Mode Enabled
- 10/22/2009 11:58:39 AM: Communication established with the Agilent 34970A device
- 10/22/2009 11:59:53 AM: Communication terminated with the Agilent 34970A device
- 10/22/2009 1:25:06 PM: Communication established with the Agilent 34970A device

The bottom status bar shows "System is generating" and "Point Time: 02:59:50".

3900 Parameter Tab

The **3900 Parameter Tab** is the primary 3900 interface for the user and is visible whenever the 3900 is connected. From this tab the user can set setpoints, view actual values and view real-time uncertainty values. The 3900 Parameters tab is divided into two different areas; the Control Parameters and the Calculated Humidity Parameters.

3900 Parameters

	Setpoint	Actual	U	
Frost Point	-10.000	-34.922	±0.099	°C
Saturation Pressure	10.000	63.646	±0.300	psia
Saturation Temperature	-8.0000	-20.730	±0.080	°C
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	21.110	21.110	±0.080	°C
Mass Flow Rate	0.5000	3.3305		l/m

Control Parameters

Calculated Humidity Parameters

3900 Generator Parameter Tab

Agilent 34970A Data

Date/Time Stamp	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]
10/21/2009 3:06:07 PM	26.729	27.63	27.925
10/21/2009 3:06:10 PM	26.745	27.604	27.933
10/21/2009 3:06:13 PM	26.742	27.604	27.942
10/21/2009 3:06:16 PM	26.768	27.529	27.95
10/21/2009 3:06:19 PM	26.747	27.365	27.957
10/21/2009 3:17:05 PM	26.859	27.805	28.257
10/21/2009 3:17:10 PM	26.802	27.802	28.252
10/21/2009 3:17:15 PM	26.765	27.803	28.246
10/21/2009 3:17:20 PM	26.682	27.804	28.243
10/21/2009 3:17:25 PM	26.671	27.806	28.242

Status Log


- Test Temperature probe is disconnected
- 10/21/2009 1:50:10 PM: 3900 Reported Expansion Valve Not Closing
- 10/21/2009 1:50:10 PM: 3900 Reported Flow Valve Not Closing
- 10/21/2009 3:05:45 PM: Communication established with the Agilent 34970A device
- 10/21/2009 3:06:03 PM: Generate Mode Requested
- 10/21/2009 3:06:07 PM: Generate Mode Enabled
- 10/21/2009 3:06:29 PM: Communication terminated with the Agilent 34970A device
- 10/21/2009 3:17:03 PM: Communication established with the Agilent 34970A device

3900 Parameters Agilent 34970A Parameters **Status Log** 3900 Reported Errors

System is generating Point Time: 00:11:22

Control Parameters

The Control Parameters contain all the control and measurement parameters critical to the operation of the humidity generator. Notice that each parameter consists of a brief title, unit of measurement, and the data values for Setpoint, Actual, and Uncertainty as applicable.




	Setpoint	Actual	U	
Frost Point	-10.000	-34.816	±0.099	°C
Saturation Pressure	10.000	62.909	±0.300	psia
Saturation Temperature	-8.0000	-20.730	±0.080	°C
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	21.110	21.110	±0.080	°C
Mass Flow Rate	0.5000	2.8685		l/m

All parameters are listed in an order similar to the generator's display with the exception of the first listed parameter. The first parameter on the Control Parameters is the active humidity control parameter. For instance, if ControLog were set to control the generator on dew point, then Dew Point would be listed as the first parameter rather than Frost Point.


In the Setpoint column is the list of the system controlled parameters. The Actual column lists the current values for each of the measured temperatures, pressures, and flow values. The Uncertainty column is a representation of the current uncertainty calculated in real-time by HumiCalc with Uncertainty using the uncertainties solutions imported into ControLog.

Changing Setpoints

To change a setpoint, click on the setpoint field that you would like to change. A Setpoint Entry box will appear and the title of the setpoint being changed will be underlined. For example, to change the Frost Point setpoint click on the Frost Point setpoint field.



Enter the new value into the Setpoint Entry box and select Ok.



-15.0

Cancel Ok

Notice that the Frost Point setpoint field updates to the new value and the values in the Actual column begin moving toward the new setpoint.

3900 Parameters			
	Setpoint	Actual	U
Frost Point	-15.000	-34.710	±0.099 °C
Saturation Pressure	10.000	62.178	±0.300 psia
Saturation Temperature	-8.0000	-20.730	±0.080 °C
Test Pressure	14.700	14.700	±0.050 psia
Test Temperature	21.110	21.110	±0.080 °C
Mass Flow Rate	0.5000	1.9939	l/m

The user can also drop down a Mini version of HumiCalc to help calculate the desired setpoint by clicking the drop down arrow on the Setpoint Entry box. For example, let's say the user wants to calculate the required Saturation Temperature needed to generate a -50.0° Frost Point with a limited supply pressure of only 100 psia. Start by clicking the Saturation Temperature setpoint field and then click the dropdown at the right of the Setpoint Entry box.

HumiCalc

Configuration

Known: Frost Point

Known Values

Frost Point: -15.0

Saturation Pressure: 100.0

Test Temperature: 21.11

Test Pressure: 14.7

Calculate

HumiCalc Drop Down

Next, select the known to be Frost Point and enter the know values for both Frost Point and Saturation Pressure.

The screenshot shows the 'HumiCalc' window. A callout 'Select the known' points to the 'Known' dropdown menu, which is set to 'Frost Point'. Another callout 'Enter the known values' points to the input fields for 'Frost Point' (containing -50.0) and 'Saturation Pressure' (containing 100.0). A third callout 'Select Calculate' points to the 'Calculate' button. Below the input fields are 'Test Temperature' (21.11) and 'Test Pressure' (14.7). A green checkmark icon is visible next to the 'Calculate' button. The background shows a 'Status Log' window with a list of timestamps.

Clicking the “Calculate” button will result in the calculated Saturation Temperature being placed in the Setpoint Entry Box and will close the HumiCalc dropdown.

The screenshot shows a 'Setpoint Entry Box' with a text field containing '-33.5803' and a dropdown arrow. Below the text field are 'Cancel' and 'Ok' buttons.

Clicking “Ok” will close the setpoint entry and will send the setpoint to the 3900 generator.

The screenshot shows a 'Status Log' window with a single entry: '10/21/2009 5:06:43 PM: Saturation Temperature Setpoint set to -33.58 °C'. The 'Status Log' button is highlighted.

Calculated Humidity Parameters

The Calculated Humidity Parameters contain values of the currently generated humidity calculated from current system temperatures and pressures. The calculated humidity parameters will also display a scrollbar as needed to allow for varying display sizes without loss of data.

%RH	0.9011	
Frost Point	-34.920	°C
Dew Point	-38.354	°C
PPMv	223.57	
PPMw	139.09	
Grains/lb	0.9735	
Enthalpy	21.569	J/g
SVP@Tt	0.3633	psia
SVP@Td	3E-003	psia
SVP@Ts	0.0140	psia

3900 Parameters

%RH

Percent Relative Humidity (%RH) is the ratio of the amount of water vapor in a given sample to the maximum amount possible at the same temperature and pressure.

Frost Point

Frost Point Temperature is the temperature to which a gas must be cooled in order to just begin condensing water vapor in the form of frost or ice, and therefore only exists at values below 0.01 °C. When operating the system with indicated Frost Points above 0.01 °C, the values indicated are to be interpreted as Dew Points. However, Frost Point is not the same as Dew Point for values below freezing. Frost Point is independent of test chamber temperature.

Dew Point

Dew Point Temperature is the temperature to which a gas must be cooled in order to just begin condensing water vapor in the form of dew. Generally, Dew Point exists at temperatures above freezing. In many instances, Dew Point may actually exist at indicated values below freezing (super-cooled dew). However, it is important to note that Dew Point is not the same as Frost Point. Dew Point is independent of test chamber temperature.

PPMv

Parts per Million by Volume is a ratio of the number of molecules of water vapor to the number of molecules of the other constituents in the gas. Once established, PPMv is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

PPMw

Parts per Million by Weight is a ratio of the weight of the water vapor in a sample to the weight of the remaining constituents in the gas. Once established, PPMw is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

Grains/lb

Grains per pound is a ratio of the weight, in grains, of water vapor to the weight, in pounds, of the other constituents in the gas. (7000 grains = 1 pound). Once established, Grains/lb is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

Enthalpy

Enthalpy is a measure of the amount of energy required to change a gas from one temperature/humidity value to another. In application, enthalpy is not used as an absolute value, but rather it is the difference in enthalpy between two distinct points which are of interest. The datum point which results in zero enthalpy was therefore arbitrarily chosen at a test temperature of 0 °C and 0 %RH. Applying enthalpy is a matter of computing the difference in enthalpy between two or more distinct data points.

SVP@Tt

Saturation Vapor Pressure (SVP) computed at the Test Temperature.

SVP@Td

Saturation Vapor Pressure (SVP) computed at the Dew/Frost Point Temperature.

SVP@Ts

Saturation Vapor Pressure (SVP) computed at the Saturation Temperature.

F@Tt.Pt

Enhancement Factor at Test Temperature and Pressure.

F@Td.Pt

Enhancement Factor at Dew/Frost Point Temperature and Test Pressure.

F@Ts.Ps

Enhancement Factor at Saturation Temperature and Pressure.

Specific Humidity

Specific Humidity is a ratio of the weight of the water vapor to the total weight of the humidified gas. Specific Humidity is independent of test chamber temperature.

Absolute Humidity

Absolute Humidity is the weight of the water vapor per unit volume of humidified gas.

Dry Air Density

Dry Air Density is the *partial* density in weight per unit volume of only the dry air portion of a moist air sample. In other words, if the water vapor were removed from a fixed volume of air, the remaining dry air would exhibit this density.

Moist Air Density

Moist Air Density is the total weight per unit volume of a moist air sample. This density includes both the weight of the air and the weight of the water vapor.

Wet Bulb Temperature

Wet Bulb temperature is used in wet bulb/dry bulb aspirated Psychrometry, and is the temperature measured by a temperature probe whose tip is coated with water (typically by being covered with a wet sock). When aspirated at a constant air velocity, the wet bulb will cool due to evaporation of the water from the probe. The cool temperature, to which it equilibrates, is used in the calculation of humidity parameters.

Mixing Ratio by Volume

Mixing Ratio by Volume is a ratio of the partial pressure of the water vapor to the partial pressure of the remaining constituents in the sample. Mixing Ratio by Volume is independent of test chamber temperature.

Mixing Ratio by Weight

Mixing Ratio by Weight is a ratio of the weight of the water vapor to the weight of the remaining constituents in the sample. Mixing Ratio by Weight is independent of test chamber temperature.

Percent by Volume

Percent by Volume is a ratio (expressed as a percentage) of the partial pressure of the water vapor to the total pressure of the sample. Percent by Volume is independent of test chamber temperature.

Percent by Weight

Percent by Weight is a ratio (expressed as a percentage) of the weight of the water vapor to the total weight of the sample. Percent by Weight is independent of test chamber temperature.

Vapor Mole Fraction

Vapor Mole Fraction is the mole fraction of water vapor in a sample.

Dry Air Mole Fraction

Dry Air Mole Fraction is the mole fraction of the dry air portion of a sample. The dry air portion is considered to be all constituents in a gas exclusive of the water vapor.

Supply Pressure

Supply pressure is the regulated pressure reading of the gas supply.

Note: The supply pressure value is gauge pressure.

Device Parameter Tabs

The **Device Parameter Tabs** show the current actual values for the given device. The tabs are visible whenever the device is connected. The Device Parameters will display a scrollbar as needed to allow for varying display sizes without loss of data.

Agilent 34970A Parameters

Parameter	Actual	Unit
Temperature 1	27.081	°C
Temperature 2	27.876	°C
Temperature 3	28.388	°C
Temperature 4	25.638	°C
Temperature 5	25.345	°C
Temperature 6	26.295	°C
Temperature 7	24.736	°C
Temperature 8	23.456	°C
Voltage	8.6673	
RTD	24.182	°C

Agilent 34970A Data

Date/Time Stamp	Temperature 1 [°C]	Temperature 2 [°C]	Temperature 3 [°C]
10/22/2009 1:25:11 PM	27.154	27.84	28.393
10/22/2009 1:25:16 PM	27.199	27.842	28.391
10/22/2009 1:25:21 PM	27.188	27.841	28.392
10/22/2009 1:25:26 PM	27.081	27.844	28.392
10/22/2009 1:25:31 PM	27.123	27.852	28.39
10/22/2009 1:25:36 PM	27.149	27.857	28.389
10/22/2009 1:25:41 PM	27.136	27.862	28.387
10/22/2009 1:25:46 PM	27.128	27.867	28.386
10/22/2009 1:25:51 PM	27.16	27.873	28.387
10/22/2009 1:25:56 PM	27.191	27.873	28.387
10/22/2009 1:26:01 PM	27.225	27.874	28.387
10/22/2009 1:26:06 PM	27.136	27.875	28.388

Status Log

- 10/22/2009 9:41:04 AM: 3900 Test Temperature probe is disconnected
- 10/22/2009 9:41:05 AM: 3900 Reported Expansion Valve Not Closing
- 10/22/2009 9:41:05 AM: 3900 Reported Flow Valve Not Closing
- 10/22/2009 10:26:41 AM: Generate Mode Requested
- 10/22/2009 10:26:41 AM: Generate Mode Enabled
- 10/22/2009 11:58:39 AM: Communication established with the Agilent 34970A device
- 10/22/2009 11:59:53 AM: Communication terminated with the Agilent 34970A device
- 10/22/2009 1:25:06 PM: Communication established with the Agilent 34970A device

Agilent 34970A Parameters | 3900 Parameters | **Status Log** | 3900 Reported Errors

System is generating | Point Time: 02:59:28


Device Parameters

Device Parameters contain all the most recent actual measurement parameters received from the device. Notice each parameter consists of a brief title, unit of measurement and the Actual data value.


Agilent 34970A Parameters ▼ ⓘ		
Actual		
Temperature 1	27.065	°C
Temperature 2	27.880	°C
Temperature 3	28.389	°C
Temperature 4	25.625	°C
Temperature 5	25.339	°C
Temperature 6	26.290	°C
Temperature 7	24.721	°C
Temperature 8	23.451	°C
Voltage	8.6673	
RTD	24.185	°C

Data and Graph Tab Group


The **Data and Graph Tab Group** is located in the middle right of the application and can contain data, graph, and profile tabs. Data and Graph tabs are fixed tab style windows that can be selected by clicking the desired tab labels at the top.




Data and Graph Tabs will expand to fill the available area as tabs around are collapsed and hidden.



Data and Graph Tabs can also be selected using the drop down arrow in the upper right hand corner. The currently selected tab is indicated by a check mark in the drop down list.



Data and Graph Tabs can be closed depending on the state of the device or generator using the “X” in the upper right hand corner.



Data Tabs

Data Tabs contain a spreadsheet type view of data.

For more information, see [Data and Data Summary](#)

Graph Tabs

Graph Tabs contain a pictorial view of data.

For more information, see [Graphing](#)

Profile Tab


The **Profile Tab** contains the profile point definitions for an auto profile.

For more information, see [Profiling](#)


Information Tab Group

The **Information Tab Group** is located on the bottom right hand side of the application. This is a docking style window that can be “pinned” open or allowed to close when not active. A particular information tab is selected by clicking its tab label at the bottom of the group. The Information Tab Group contains status information about the operation of the generator and its connected devices. The group can consist of a Status Log tab, 3900 Reported Errors tab and Interface Console tabs.


By clicking the pin icon on any information tab the user can unpin the Information tabs.




Once the information tabs are unpinned they will automatically hide.




The user can access the hidden tabs by clicking the desired information tab label at the bottom.




Tabs are also accessed via the drop down selection by clicking on the arrow icon. The currently selected tab is indicated by a check mark in the drop down list.



Each Information tab can be torn out by dragging the tab label upward. This allows multiple information tabs to be visible at once. To place the tab back into the list, drag the tabs label back to the bottom.




All Information Tabs have a context menu that is displayed by right clicking in the tab. The context menu allows quick access to functions that can clear, save and print the information.





Status Log

The **Status Log** tab contains chronological information about the system status, changes in operational modes, changes in setpoints, and runtime errors due to communication or mechanical difficulties encountered by the generator.



3900 Reported Errors

The **3900 Reported Errors** tab only appears when the 3900 reports an error. This is a very important information tab because it reports 3900 system errors to the user. These types of errors can cause the 3900 to shutdown and require immediate attention by the operator. The  icon will be displayed to help draw the attention of the user to the reported 3900 errors.



System timing

The **System Timing** tab shows information about the current timing associated with the current operation such as elapsed run time at current conditions. This window may be shown at any time by clicking the "Point Time" in the status bar and is automatically shown when an Auto Profile is started. The tab gives detailed information on the Auto Profile as it runs. Elapsed and remaining Phase, Point and Total time are listed along with the detailed assurance conditions values and tolerances.

System Timing							
Auto Profile: Generate - Point: 3 of 6 - Phase: Assurance							
	Elapsed (dd:hh:mm:ss)	Remaining (dd:hh:mm:ss)	Assurance Condition	Actual Value	Min Value	Max Value	Std Dev
Phase Time:	00:00:00:49	00:00:00:00	Saturation Temperature (2500)	-87.869	24.900	25.100	5.5845E-004
Point Time:	00:00:00:49	00:01:00:00	Saturation Pressure (2500)	-23.779	0.3797	0.5797	0.0310
Total Time:	00:00:00:49	00:04:00:00					


At the bottom of the window, there is a tab labeled "Status Log" and a yellow warning icon followed by the text "2500 Reported Errors". The "System Timing" tab is currently selected.

3900 Interface Console


The **3900 Interface Console** tab allows the user to send and receive commands to and from the 3900. The console tab is opened by selecting “Open 3900 Interface Console” from the Connections>Settings menu whenever the 3900 is connected.

Warning: Interfacing with the 3900 using the Console should only be attempted if instructed to do so by Thunder Scientific.


Commands are entered in the bottom text area of the console tab and the 3900 responses are displayed in the upper text area of the console tab.




To enter a command, type the desired command in the bottom text area and press return.




Once the command has been entered it will appear in bold in the upper text area. A processing status will be indicated in the bottom text area while ControLog sends and waits for the response from the 3900. The user is not allowed to enter another command until the previous command has been processed.



Once the command has been successfully processed the 3900 response will appear in the upper text area right below the command sent.



If ControLog does not receive a response from the 3900 within a given time frame an error message will be displayed.




Device Interface Console

The **Device Interface Console** tab allows the user to view the commands being sent to and received from any given connected device. The device console tab is opened by selecting "Open Device Interface Console" from the Connections>Settings menu whenever the a device is connected.


The user can select which connected device to view using the drop down selection at the top of the tab.



The data sent to the device is indicated by a bold “**Sent:**” label in the lower text area of the console tab.




The data received from the device is indicated by a bold “**Received:**” label in the lower text area of the console tab.



Status Bar

The **Status Bar** is located on the bottom of the application window. The Status Bar displays the current Connection and Run Status of the generator, Auto Profile Controls and Status, and current Point Time.




Connection and Run Status

The **Connection and Run Status** is shown on the left hand side of the status bar and gives the user a quick visual and textual reference to the current state of the 3900 generator. As new events are recorded into the Status Log, the event will be displayed for a short period of time in the status bar to inform the user of the new event.


Disconnected

The status bar will show a broken link icon and “Disconnected” when the 3900 is not connected to ControlLog;




Purging

The status bar will show a green circular arrow icon and “System is purging” when the 3900 is purging.




If the 3900 is reporting an error, the green circular arrow icon will also appear with a small yellow warning.




Generating

The status bar will show a green forward arrow icon and “System is generating” when the 3900 is generating.




If the 3900 is reporting an error, the green forward arrow icon will also appear with a small yellow warning.




Shutdown

The status bar will show a red stop icon and “System is shutdown” when the 3900 is shutdown.




If the 3900 is reporting an error, the red stop icon will also appear with a small yellow warning.



New Event

During any state the most recent status log event will be shown briefly in the status bar.




Auto Profile Controls and Status

The **Auto Profile Controls and Status** is shown on the right hand side of the status bar and gives the user quick controls and status over a running profile. The Auto Profile consists of three status parts and three control parts.

For more information, see [Auto Profiling](#)


Auto Profile Controls

The Auto Profile controls consist of a “Pause” button, “Next Point” button and a “Next Phase” button.




Auto Profile Status

The Auto Profile status consists of the generator “Run Mode” for the current profile point, the “Profile Point” the system is currently running and the “Phase” of the current point.




Current Point Time

The **Current Point Time** is shown on the right hand side of the status bar and gives the user a quick display to the amount of time the system has been at point. Point Time is not the amount of time at setpoint but simply the amount of time since the last setpoint or mode change.



The user can click the point time to open the system timing tab.




Graphing

Graphing is a powerful tool used to view previously recorded data or to monitor the current data in real time. The graph works hand in hand with the data tabs. While the generator is in operation, data tabs store the most recent data points from the connected devices at the desired interval. A graph can be used to create a visual picture of this stored data.

Graph operations can be accessed by two means, either by selecting the desired command from the graph menu at the top or by right clicking a graph tab and opening a context menu. The functionality of the commands are the same regardless of which method is used, but remember that the menu commands are dynamic and reflect operations that can be performed on the selected graph tab.

Note: *ControLog has a limit of ten graph tabs opened at any given time.*


Each graph tab consists of a two dimension plot across an X and Y Axis. The graph can be customized to display different point symbols, various line colors, a legend and a header.




Creating a New Graph

A new graph can be created using the New Graph Wizard dialog. The wizard will step the user through the selection process of what data the user would like to include in the new graph and how it should look.


To create a new graph, select “New” from the main menu or right click a graph tab and select “New” from its context menu.




A “New Graph Wizard” dialog will appear.



The first page of the New Graph Wizard is where the user selects which device items they would like to include in the graph. On the left hand side is a list of all available devices. Selecting a device will result in the “Device Items” list being updated to reflect the available items for the selected device. To add an item, highlight the desired item or items in the “Device Items” list and click the “Add” button.




Repeat this process until all the desired items are listed on the right side.




Note: The user can invert any selection by right clicking. This will highlight all items that are not currently highlighted and will remove highlight from any items that are currently highlighted.

To remove an item from the list of items to graph the user can either click the “Clear All” button to remove all items or simply select the desired item from the “Items to Graph” list and click the “Remove” button.



The **Graph Properties** dialog box contains three main sections: **Devices:**, **Device Items:**, and **Items to Graph:**. The **Devices:** list shows '3900' and 'Simulated Chilled Mirror'. The **Device Items:** list shows 'Dew Point', 'Frost Point', and 'Test Pressure'. The **Items to Graph:** list shows 'Frost Point (3900)' and 'Frost Point (Simulated Chilled Mirror)'. Between the lists are three buttons: **ADD** (plus icon), **REMOVE** (minus icon), and **CLEAR ALL** (X icon). At the bottom are navigation buttons: **< Back**, **Next >**, **Finish**, and **Cancel**.

Once complete, clicking the “Next” button will bring up the Properties page. From the properties page the user can customize the look of each graphed line along with the chart itself.



The **New Graph Wizard** dialog box is divided into two sections: **Line Properties:** and **Chart Properties:**. The **Line Properties:** section has a list box with 'Frost Point (3900)'. To its right are settings for **Color:** (DarkGoldenrod), **Symbol:** (Box), a **Smoothed** checkbox, and a **Show Data Label** checkbox with a frequency of **Every 10 Points**. The **Chart Properties** section includes a **Tab Name** field (Frost Point), and checkboxes for **Show Header** (3900 Frost Point) and **Show Legend**. Navigation buttons **< Back**, **Next >**, **Finish**, and **Cancel** are at the bottom.

Line Properties define the color and point symbol for the selected line. Select the desired line from the selection box and set the Color, Symbol, whether to smooth the line, and whether to show data labels and how often to show them.

Color defines the color that the line will be drawn in. To change the color, select the desired color from the drop down list. A sample of the selected color is displayed in the small box to the left of the drop down.

Symbol defines the symbol drawn at each point. The symbol will be the same color as the line and will be shown at each data point. To change the symbol, select the desired symbol from the drop down list.

Smoothed defines whether to smooth the line between points. This will remove any sudden angles between points and will draw a smooth flowing line between points.

Show Data Label defines whether to show data labels and how often to show them. Data labels are small labels appearing next to a point displaying the value at the point.


Chart Properties define the name of the graph tab, indicated whether to show a header on the graph and whether to show a legend of the lines plotted.

Tab Name defines the name of the graph tab. This is the name the user will see appear in the Data and Graph Tab Group.

Show Header defines the header that will appear at the top of the graph. To add a header to the graph check the checkbox and enter the desired text description for the header.

Show Legend defines whether a legend will be displayed on the right hand side of the graph indicating which colored line is which data item.

Once complete, clicking the “Finish” button will display the new graph in the Data and Graph Tab Group.




Customizing a Graph

Each graph tab can be customized to display the data in different means. The user can Pan, Zoom and Scale the graph to the desired appearance.

Pan


The graph can be panned up and down as well as left and right by a left click and hold of the mouse button while the user moves the cursor around. Panning is useful when you have zoomed the graph and want to view different parts of the data without changing the scaling.




Zoom

Selecting the **Zoom** command checks the operation in the menu and then allows the user to use the left mouse button to create a rectangular area on the graph that will be zoomed. This operation is available only when a Graph tab is selected.

To perform the zoom, select “Zoom” from the main menu or right click a graph tab and select “Zoom”.




Move the cursor to the graph, left click and drag the cursor to create a box around the portion to zoom and release the left mouse button. The portion of the graph within the drag box will expand to fill the entire graph. The time and Y-axis limits update accordingly.




Zoom Graph's X Axis

Selecting the **Zoom X Axis** command checks the operation in the menu and then allows the user to use the left mouse button to create a section area on the graph that will be zoomed along the X Axis. This operation is available only when a Graph tab is selected.

To perform the zoom, select “Zoom X Axis” from the main menu or right click a graph tab and select “Zoom X Axis”.




Left click on the graph and drag the cursor so that the portion of the X-axis (time axis) of interest is contained within the two vertical dashed lines. The portion contained within this region will expand to fill the entire X-axis. The Y-axis remains unchanged.




Zoom Graph's Y Axis

Selecting the **Zoom Y Axis** command checks the operation in the menu and then allows the user to use the left mouse button to create a section area on the graph that will be zoomed along the Y Axis. This operation is available only when a Graph tab is selected.

To perform the zoom, select “Zoom Y Axis” from the main menu or right click a graph tab and select “Zoom Y Axis”.



Left click on the graph and drag the cursor so that the portion of the Y-axis of interest is contained within the two horizontal dashed lines. The portion contained within this region will expand to fill the entire Y-axis. The X-axis remains unchanged.



Auto Scale

Selecting the **Auto Scale** command will automatically reset both axis of the graph so the entire data set for each selected item is contained within the boundaries of the graph. This operation is available only when a Graph tab is selected.

Scale

Selecting the **Scale** command checks the operation in the menu and allows the user to use the left mouse button to scale. Dragging the cursor up scales the display in (zooms in) and dragging the cursor down scales the display out (zoom out). This operation is available only when a Graph tab is selected.

Scale X Axis


The **Scale X Axis** command allows the user to scale the X axis. Selecting the command checks the operation in the menu and then allows the user to use the left mouse button to scale. Dragging the cursor up scales the X Axis in (zooms X Axis in) and dragging the cursor down scales the X Axis out (zooms X Axis out). This operation is available only when a Graph tab is selected.

Scale Y Axis

The **Scale Y Axis** command allows the user to scale the Y axis. Selecting the command checks the operation in the menu and then allows the user to use the left mouse button to scale. Dragging the cursor up scales the Y Axis in (zooms Y Axis in) and dragging the cursor down scales the Y Axis out (zooms Y Axis out). This operation is available only when a Graph tab is selected.

Graph Properties

The graph **Properties** command allows the user to modify the properties of the selected graph. Selecting the command opens the Graph Properties dialog that allows the user to make changes to what data is graphed, line properties and axis values. This operation is available only when a Graph tab is selected.



The **Graph Properties** dialog box is divided into two main sections: **Line Properties** and **Chart Properties**.

Line Properties:


- A list box on the left contains the selected data series: **Frost Point (3900)**.
- To the right, the **Color** is set to **DarkGoldenrod** (represented by a color swatch and a dropdown menu).
- The **Symbol** is set to **None** (represented by a dropdown menu).
- There are two unchecked checkboxes: **Smoothed** and **Show Data Label**.

Chart Properties:

- The **Tab Name** is **Frost Point** (text in a field).
- There are two unchecked checkboxes: **Show Header** and **Show Legend**.

At the bottom of the dialog are four buttons: **< Back**, **Next >** (highlighted in yellow), **Finish**, and **Cancel**.

The first two pages of the Graph Properties dialog are the same as the New Graph Wizard that is used to create new graphs. The Graph Properties dialog starts on the properties page from which the user can customize the look of each graphed line as well as the chart itself. Clicking the “Back” button will move back to the data selection page where the user can select which device items they would like to include in the graph. Clicking the “Next” button will move to the Axis page where the user can specify the starting and ending X and Y axis values.



Y Axis Limits defines the maximum, minimum, middle and span values for the Y Axis.

The **Maximum** value defines the maximum Y value for the Y Axis. No values beyond this maximum will be displayed on the graph.

The **Minimum** value defines the minimum Y value for the Y Axis. No values below this minimum will be displayed on the graph.

The **Middle** value defines the middle Y value for the Y Axis.

The **Span** value defines the amount above and below the middle value where the maximum and minimum Y values lie.

The **Auto** check box tells ControlLog to automatically calculate the best Y Axis limits to encompass the current data.

Note: The maximum, minimum, middle and span values are interrelated and changing any one value may result in another value automatically changing to assure all values mathematically equate.

X Axis Limits defines the maximum and minimum values for the X Axis.


The **Minimum** value defines the minimum date and time for the X Axis. No values below this minimum date and time will be displayed on the graph.

The **Auto Min** check box tells ControlLog to automatically use the starting date and time for the current data as the X Axis minimum.

The **Maximum** value defines the maximum date and time for the X Axis. No values above this maximum date and time will be displayed on the graph.

The **Auto Max** check box tells ControlLog to automatically use the last date and time for the current data as the X Axis maximum. The graph will constantly expand as new data points are recorded.

Once complete, clicking the “Finish” button will display the graph in the same tab with the new property settings.




Saving a Graph


Each graph tab can be saved to a file in any of the following graphic file types:

- Bitmap (*.bmp)
- Graphics Interchange Format (*.gif)
- Joint Photographic Experts Group (*.jpg)
- W3C Portable Network Graphics (*.png)
- EMF Enhanced Metafile Format (*.emf)

To perform the save, select “Save” from the main menu or right click a graph tab and select “Save”.




or



Using the “Save” dialog, select the location, name and graphic type you want to save the graph as.

Note: ControLog defaults the file name to the name of the graph tab.




Printing a Graph


Each graph tab can be printed to any of the PC's installed printers.

Note: You must have a printer installed before you can print using ControlLog.

To perform the print, select “Print” from the main menu or right click a graph tab and select “Print”.



Using the “Print” dialog, select the desired printer, range and number of copies you want to print.



Data and Data Summary

ControLog stores data into individual Data Tabs. Each data tab contains a spreadsheet type view that consists of a date/time stamp and the measured data items corresponding to that date/time stamp. Data tabs consist of three similar but different types: Device Data, File Data and Data Summary. Each type has the same spreadsheet type view and operation, but all three have different data sources.

Note: The data tab data is always stored in SI units regardless of the current system unit settings. The only exception is for Data Summary tabs which are created using currently selected system units but will not update on further unit changes.

Date/Time Stamp	Saturation Pressure Setpoint [kPa]	Saturation Pressure [kPa]
10/23/2009 11:36:38 AM	68.947573	435.648
10/23/2009 11:37:08 AM		.71763
10/23/2009 11:37:38 AM		.86837
10/23/2009 11:38:08 AM		.43736
10/23/2009 11:38:38 AM	68.947573	431.9924
10/23/2009 11:39:08 AM	68.947573	429.57234
10/23/2009 11:39:38 AM	68.947573	434.17114
10/23/2009 11:40:08 AM	68.947573	439.35186
10/23/2009 11:40:38 AM		437.19863
10/23/2009 11:41:08 AM		432.04893
10/23/2009 11:41:38 AM		435.03436
10/23/2009 11:42:08 AM	68.947573	433.47822
10/23/2009 11:42:38 AM	68.947573	432.09375
10/23/2009 11:43:08 AM	68.947573	434.48761
10/23/2009 11:43:38 AM	68.947573	432.3261
10/23/2009 11:44:08 AM	68.947573	438.90784

The user can navigate through the data using the scroll bars .

3900 Data | Data Test (3900).txt | Frost Point | Data Summary | Auto Profile

3900 Data

Date/Time Stamp	Frost Point $\pm U$ (k=2) [°C]	Dew Point [°C]	Dew Point $\pm U$ (k=2) [°C]
10/23/2009 11:36:38 AM	0.098796176	-38.2863	0.1006363
10/23/2009 11:37:08 AM	0.098792679	-38.2879	0.1006262
10/23/2009 11:37:38 AM	0.098735141	-38.3133	0.1005555
10/23/2009 11:38:08 AM			0.1004040
10/23/2009 11:38:38 AM			0.1008636
10/23/2009 11:39:08 AM			0.101010
10/23/2009 11:39:38 AM			0.1007272
10/23/2009 11:40:08 AM	0.098612551	-38.3681	0.1004040
10/23/2009 11:40:38 AM	0.098718715	-38.3206	0.1005353
10/23/2009 11:41:08 AM	0.098979286	-38.2061	0.1008686
10/23/2009 11:41:38 AM	0.098827066	-38.2727	0.1006767
10/23/2009 11:42:08 AM	0.098906007	-38.2381	0.1007676
10/23/2009 11:42:38 AM	0.098976977	-38.2071	0.1008585
10/23/2009 11:43:08 AM	0.098854703	-38.2605	0.1007070
10/23/2009 11:43:38 AM	0.098965018	-38.2123	0.1008484
10/23/2009 11:44:08 AM	0.098634311	-38.3583	0.1004343

Date/Time Stamp column stays fixed as the data is scrolled

The user can also select specific data by clicking and dragging the desired cells.


3900 Data | Data Test (3900).txt | Frost Point | Data Summary | Auto Profile

3900 Data

Date/Time Stamp	Saturation Pressure Setpoint [kPa]	Saturation Pressure [kPa]
10/23/2009 12:39:40 PM	68.947573	436.79942
10/23/2009 12:40:10 PM	68.947573	435.07091
10/23/2009 12:40:40 PM	68.947573	426.93992
10/23/2009 12:41:10 PM	68.947573	433.27482
10/23/2009 12:41:40 PM	68.947573	434.62344
10/23/2009 12:42:10 PM	68.947573	435.51148
10/23/2009 12:42:40 PM	68.947573	437.4558
10/23/2009 12:43:10 PM	68.947573	440.44399
10/23/2009 12:43:40 PM	68.947573	432.96938
10/23/2009 12:44:10 PM	68.947573	429.14004
10/23/2009 12:44:40 PM	68.947573	434.60482
10/23/2009 12:45:10 PM	68.947573	437.33446
10/23/2009 12:45:40 PM	68.947573	437.96326
10/23/2009 12:46:10 PM	68.947573	431.87105
10/23/2009 12:46:40 PM	68.947573	439.12089
10/23/2009 12:47:10 PM	68.947573	436.18993

User selected data


The user can then copy selected data to the clipboard by selecting “Copy Selection” from the context menu or by using the keyboard shortcut combination of Ctrl+C.



Device Data Tabs

The **Device Data Tabs** contain stored data values obtained by the connected generator or device. After establishing communication with the generator or device an individual data tab for the device will be created. These tabs store the data readings from the connected device at the specified data storage interval.

Note: Data is only recorded while the 3900 is either in generate or purge mode. Data is also stored at the generate rate whenever a device is connected and the 3900 is not connected. This gives the user the ability to use ControLog as a logging application for any device they connect without the need of a 3900 generator.



The screenshot shows the ControLog software interface with the '3900 Data' tab selected. The table displays data for 'Agilent 34970A Data'. The columns are 'Date/Time', 'Duration Pressure Setpoint [kPa]', and 'Saturation Pressure [kPa]'. The data is recorded at 10-second intervals. A callout points to the '3900 Data' tab, and another points to the 'Generator Data Tab' label. A third callout states 'Data is stored at the selected data storage interval'.

Date/Time	Duration Pressure Setpoint [kPa]	Saturation Pressure [kPa]
11/6/2009	68.947573	438.8513
11/6/2009 1:07:15 PM	68.947573	434.37936
11/6/2009 1:07:25 PM	68.947573	437.09865
11/6/2009 1:07:35 PM	68.947573	432.30059
11/6/2009 1:07:45 PM	68.947573	436.73599
11/6/2009 1:07:55 PM	68.947573	432.8108
11/6/2009 1:08:05 PM	68.947573	431.20364
11/6/2009 1:08:15 PM	68.947573	428.52089
11/6/2009 1:08:25 PM	68.947573	427.09988
11/6/2009 1:08:35 PM	68.947573	434.18769
11/6/2009 1:08:45 PM	68.947573	432.15304
11/6/2009 1:08:55 PM	68.947573	435.25775
11/6/2009 1:09:05 PM	68.947573	437.17312
11/6/2009 1:09:15 PM	68.947573	438.53001
11/6/2009 1:09:25 PM	68.947573	435.72177
11/6/2009 1:09:35 PM	68.947573	437.57922


File Data Tabs

The **File Data Tabs** contain data values loaded from a previously saved Device Data Tab. ControLog can open data saved in the following types and formats:


- Text File (Comma Delimited) (*.csv;*.txt)
- Text File (Tab Delimited) (*.dat;*.txt)
- Excel Workbook (*.xlsx;*.xls)
- Backup ControLog File (*.backup)

The only requirement for loading data from the above mentioned data files is the data must be formatted so date/time values appear in the first column and all other columns contain numeric values.


To open a data file, select “Open” from the main menu or right click a data tab and select “Open”.




or



Using the “Load Data File” dialog, browse and select the data file you want to open.



Once the load is complete the loaded file data tab will be displayed.




The screenshot shows a software window titled '3900 Data' with a menu bar containing 'Units', 'Connections', and 'Help'. A tab labeled 'Data Test (DP 30 with Options).txt' is active. Below the tab is a table with the following data:

Date/Time Stamp	Balanced	Loaded File Data Tab	Dew Point Temperature [°C]
10/7/2009 5:59:18 PM	0		999.9
10/7/2009 5:59:23 PM	0	0	26.2
10/7/2009 5:59:28 PM	0	0	26.2
10/7/2009 5:59:33 PM	0	0	26.2
10/7/2009 5:59:38 PM	0	0	26.2
10/7/2009 5:59:43 PM	0	0	26.2
10/7/2009 5:59:48 PM	0	0	26.2
10/7/2009 5:59:53 PM	0	0	26.2
10/7/2009 5:59:58 PM	0	0	26.2
10/7/2009 6:00:03 PM	0	0	26.2
10/7/2009 6:00:08 PM	0	0	26.2
10/7/2009 6:00:13 PM	0	0	26.2
10/7/2009 6:00:18 PM	0	0	26.2
10/7/2009 6:00:23 PM	0	0	26.2
10/7/2009 6:00:28 PM	0	0	26.2
10/7/2009 6:00:33 PM	0	0	26.2

Data Summary Tabs

The **Data Summary Tabs** contain data values generated from a Data Summary. The data summary allows the user to summarize the available data into a single data tab. The data summary also allows the user to calculate error or differences between a selected standard values and selected device values.

Note: Data Summary tabs are the only types of Data Tabs that are created using the currently selected system units instead of SI units.




The screenshot shows the 3900 ControlLog software interface. The 'Data Summary' tab is selected, displaying a table of data. A yellow callout bubble points to the tab with the text 'Generated Data Summary Tab'.

Date/Time Stamp	Data Test (DP 30 with Op
10/7/2009 5:59:18 PM	
10/7/2009 6:04:18 PM	26.54
10/7/2009 6:09:19 PM	8.33
10/7/2009 6:14:19 PM	-0.84
10/7/2009 6:19:20 PM	-8.53
10/7/2009 6:24:20 PM	-11.95
10/7/2009 6:29:21 PM	-11.18
10/7/2009 6:34:21 PM	-10.96
10/7/2009 6:39:22 PM	-31.64
10/7/2009 6:44:22 PM	-8.75
10/7/2009 6:49:22 PM	-8.72
10/7/2009 6:54:23 PM	-8.73
10/7/2009 6:59:23 PM	-8.78
10/7/2009 7:04:24 PM	-8.77
10/7/2009 7:09:24 PM	-9.85
10/7/2009 7:14:25 PM	-9.86
10/7/2009 7:19:25 PM	-9.88


Creating a Data Summary

To create a Data Summary, select “Create Data Summary” from the main menu or right click a data tab and select “Create Data Summary”. Selection will open the Data Summary Wizard dialog that will step the user through the creation process.



or

The first step in the data summary creation process is to give the data summary a name. This will be the name of the tab that appears in the Data and Graph Tab Group.



The next step in the data summary creation process is to select the standard. The standard will be used as the reference to calculate differences between the device data items if the user chooses to do so.

On the left hand side is the available devices to choose from and on the right hand side is the selected standard to add to the summary. To select a device, highlight it and click the “Add” button. In almost all cases the standard will be the 3900.

Note: Only one device can be selected as the standard for the summary.


The screenshot shows a window titled "Data Summary" with a close button in the top right corner. Inside the window, there is a section titled "Select the standard to include in the summary:" followed by a descriptive paragraph: "Select the single standard to include in the summary. In most cases the standard will be the generator, but can be any device with available data. The standard will be used as the reference to calculate differences between device data items." Below this, there are two lists and two buttons. The left list, titled "Available standards:", contains "3900" and "Simulated Chilled Mirror", with the latter highlighted in blue. The right list, titled "Standard to add to the summary:", is currently empty. Between the lists are two buttons: a green "ADD" button with a plus sign and a grey "REMOVE" button with a minus sign. At the bottom of the window are four buttons: "< Back", "Next >", "Finish", and "Cancel".

To remove an item as the selected standard, select the desired item from the right hand side and click the “Remove” button.

This screenshot shows the same "Data Summary" window. In this state, the "Simulated Chilled Mirror" item has been moved from the right list to the left list, "Available standards:". The "Available standards:" list now contains "3900" and "Simulated Chilled Mirror". The "Standard to add to the summary:" list is empty. The "ADD" button is now grey and the "REMOVE" button is red with a minus sign. The rest of the window, including the descriptive text and bottom navigation buttons, remains the same.

Once complete, selecting the “Next” button will allow the user to select each device they would like to include in the summary. On the left hand side are the available devices to choose from and on the right hand side are the selected devices to add to the summary. To select a device, highlight it and click the “Add” button.

Note: Multiple devices can be selected and added to the summary.



Note: The user can invert any selection by right clicking. This will highlight all items that are not currently highlighted and will remove highlight from any items that are currently highlighted.

To remove an item from the list of devices, select the desired item or items from the right hand side and click the “Remove” button.

The dialog box is titled "Data Summary" and contains the instruction "Select each device to include in the summary:". Below this is a sub-instruction: "Select each device you would like to include in the summary." The interface is divided into two main sections: "Available devices:" on the left and "Devices to add to the summary:" on the right. The "Available devices:" list contains "Test (Manual)" and "Simulated Chilled Mirror". The "Devices to add to the summary:" list also contains "Test (Manual)" and "Simulated Chilled Mirror", with "Test (Manual)" currently selected. Between these two lists are two buttons: a green "ADD" button with a plus icon and a red "REMOVE" button with a minus icon. At the bottom of the dialog are four buttons: "< Back", "Next >", "Finish", and "Cancel".

Once complete, selecting the “Next” button will allow the user to select the standard’s data items they would like to include in the summary. On the left hand side are the available data items to choose from and on the right hand side are the selected data items to add to the summary. To select a data item or items, highlight them and click the “Add” button.

The dialog box is titled "Data Summary" and contains the instruction "Select the 3900 data to use in the summary:". Below this is a sub-instruction: "Select each data item for this device that you would like to include in the summary." The interface is divided into two main sections: "Available data:" on the left and "Data to add to the summary:" on the right. The "Available data:" list contains "Error Code", "F@Td,Pt", "F@Ts,Ps", "F@Tt,Pt", "Frost Point", "Frost Point $\pm U$ (k=2)", "Grains/lb", "Mass Flow Rate", "Mass Flow Rate Setpoint", and "Mixing Ratio by Volume". The "Frost Point $\pm U$ (k=2)" item is currently selected. The "Data to add to the summary:" list is empty. Between these two lists are two buttons: a green "ADD" button with a plus icon and a red "REMOVE" button with a minus icon. At the bottom of the dialog are four buttons: "< Back", "Next >", "Finish", and "Cancel".

To remove an item or items from the list of data items, select the desired item or items from the right hand side and click the “Remove” button.

Data Summary

Select the 3900 data to use in the summary:

Select each data item for this device that you would like to include in the summary.

Available data:

- Error Code
- F@Td,Pt
- F@Ts,Ps
- F@Tt,Pt
- Frost Point
- Frost Point $\pm U$ ($\kappa=2$)
- Grains/lb
- Mass Flow Rate
- Mass Flow Rate Setpoint
- Mixing Ratio by Volume

Data to add to the summary:

- Frost Point
- Frost Point $\pm U$ ($\kappa=2$)
- Grains/lb

Buttons: ADD, REMOVE, < Back, Next >, Finish, Cancel

Selecting the “Next” button will repeat the process of selecting data for the next device in the series. This will continue until the user defines all the device data they wish to include in the data summary.

Once the user has completed selecting data, ControLog will ask the user to pair each device data item with its corresponding standard data item. This tells ControLog which device data item corresponds with which standard data item so that a difference can be calculated and the error between the standard and the device can be included in the summary. If the user does not desire to calculate error for a particular data item they simple need not pair it.

To pair a data item, select the desired device data item and the standard data item you would like to pair it with. Select “Add” to create the pairing.

The screenshot shows a window titled "Data Summary" with a close button in the top right corner. The main heading is "Pair the Simulated Chilled Mirror data with its corresponding 3900 data:". Below this is a paragraph of instructions: "Pair each device data item with its corresponding standard data item. This tells the ControLog which device data item corresponds with which standard data item so that a difference can be calculated to show the error between the standard and the device. If you do not desire to calculate error for a particular data item then simply do not pair it." The interface is divided into three main sections: "Device Data:", "Standard data:", and "Paired Data:". Under "Device Data:", there is a list box containing "Frost Point". Under "Standard data:", there is a list box containing "Frost Point" and "Frost Point $\pm U$ ($k=2$)". Between these two list boxes are two buttons: a green "ADD" button with a plus sign and a grey "REMOVE" button with a minus sign. The "Paired Data:" section is an empty list box. At the bottom of the window are four buttons: "< Back", "Next >", "Finish", and "Cancel".

To remove a data item pairing, select the desired pair or pairings from the right hand side and click the “Remove” button.


This screenshot shows the same "Data Summary" window, but now the "Paired Data:" list box contains the entry "Frost Point = Frost Point". The "ADD" button is now grey and disabled, while the "REMOVE" button is red and active. The other elements of the interface, including the "Device Data:", "Standard data:", and navigation buttons at the bottom, remain the same as in the previous screenshot.

Selecting the “Next” button will repeat the process of pairing data for the next device in the series. This will continue until the user defines all the desired data pairs they wish to include in the data summary.

Once data pairing is complete, the user will be asked to select the time range and interval that the data summary will cover. The selected device’s time range will be used to determine which points to include. The data will begin at the

closest point to the selected start time and will end at the closest point to the selected end time. The summary will include all available data points between the start and end times at the closest points to the specified interval. If a particular device does not have a corresponding time for a given base time, then the value fields will be left blank for that device for that given time.

Note: For manual devices it is recommended to select the manual device as the bases of the time range. It is also recommended to select a small time interval (1 second) to assure all data is encompassed because often manual devices have irregular time intervals between points.



The image shows a software dialog box titled "Data Summary". It contains a section titled "Select the time range the summary will cover:" with a descriptive paragraph. Below this, there are three input fields: "Device Time Range to Use:" with a list box containing "3900" and "Simulated Chilled Mirror"; "Start Time:" with a date/time picker set to "10/29/09 12:09:14 PM"; and "End Time:" with a date/time picker set to "10/29/09 1:10:47 PM". To the right of these fields is a note. Below the time fields is a "Time Interval:" section with a numeric input set to "1" and a dropdown menu set to "minutes". At the bottom are four buttons: "< Back", "Next >", "Finish" (highlighted in yellow), and "Cancel".

Data Summary

Select the time range the summary will cover:

Select the device for which you want to base the data summary's time range on. The data will begin at the closest point to the selected start time and will end at the closest point to the selected end time. The summary will include all available data points between the start and end times at the closest points to the specified interval.

Device Time Range to Use:

- 3900
- Simulated Chilled Mirror

Start Time: 10/29/09 12:09:14 PM

End Time: 10/29/09 1:10:47 PM

Time Interval: 1 minutes

Note: For manual devices it is recommend to select the manual device as the bases of the time range. It is also recommended to select a small time interval (1 second) to assure all data in encompassed because often manual devices have irregular time intervals between points.

< Back Next > **Finish** Cancel

Upon clicking the “Finish” button, ControlLog will open a new data tab with the newly created data summary


The screenshot shows the 'Data Summary' window in ControlLog. The window has three tabs: '3900 Data', 'Simulated Chilled Mirror Data', and 'Data Summary'. The 'Data Summary' tab is active, displaying a table with the following columns: 'Date/Time Stamp', '3900: Frost Point [°C]', '3900: Frost Point ±U (k=2) [°C]', 'Simulated Chilled Mirror: Frost Point [°C]', and 'Simulated Chilled Mirror: Frost Point Error [°C]'. The table contains 18 rows of data. Callouts are present: 'Standard' points to the '3900: Frost Point [°C]' column; 'Device' points to the 'Simulated Chilled Mirror: Frost Point [°C]' column; 'Calculated difference between the standard and the device' points to the 'Simulated Chilled Mirror: Frost Point Error [°C]' column; and 'Date/Time Stamp increments per the selected time interval' points to the 'Date/Time Stamp' column.

Date/Time Stamp	3900: Frost Point [°C]	3900: Frost Point ±U (k=2) [°C]	Simulated Chilled Mirror: Frost Point [°C]	Simulated Chilled Mirror: Frost Point Error [°C]
10/29/2009 12:09:14 PM	-34.7702	0.099005483	-37.033	-2.262
10/29/2009 12:10:14 PM	-34.7702	0.0989675	-37.022	-2.236
10/29/2009 12:11:14 PM	-34.8	0.098922345	-37.16	-2.358
10/29/2009 12:12:14 PM	-34.8737	0.098753913	-37.19	-2.316
10/29/2009 12:13:14 PM	-34.908	0.098671535	-37.248	-2.340
10/29/2009 12:14:14 PM	-34.8649	0.098774875	-37.164	-2.299
10/29/2009 12:15:14 PM	-34.8046	0.098921358	-37.136	-2.331
10/29/2009 12:16:14 PM	-34.7917	0.098952794	-37.128	-2.336
10/29/2009 12:17:14 PM	-34.8689	0.098765386	-37.194	-2.325
10/29/2009 12:18:14 PM	-34.7075	0.099160918	-37.117	-2.409
10/29/2009 12:19:14 PM	-34.8555	0.098797631	-37.087	-2.231
10/29/2009 12:20:14 PM	-34.7558	0.099041336	-37.068	-2.312
10/29/2009 12:21:14 PM	-34.7884	0.098998488	-37.184	-2.414
10/29/2009 12:22:14 PM	-34.7884	0.098998488	-37.035	-2.225
10/29/2009 12:23:14 PM	-34.7884	0.098998488	-37.145	-2.250
10/29/2009 12:24:15 PM	-34.7884	0.098998488	-37.024	-2.199
10/29/2009 12:25:15 PM	-34.843	0.098827831	-37.04	-2.199
10/29/2009 12:26:15 PM	-34.7787	0.098984689	-37.142	-2.364

Recorded Data Points Tab

The **Recorded Data Points Tab** contains the recorded data points that have been taken either manually by the user, after a manual device entry or at the completion of a soak phase during an auto profile. Each data point can also calculate the average and or standard deviation for the defined number of prior points taken with each recorded data point.


Note: Each entry point in the Recorded Data Points tab is created using the currently selected system units.




Recorded Data Points			
Point 1	2500		
Date/Time Stamp	Dry Air Density [g/m ³]		
8/27/2012 1:32:31 PM	37.572888		0.466373
8/27/2012 1:32:36 PM	37.600163		0.466694
8/27/2012 1:32:41 PM	37.897638		0.470621
8/27/2012 1:32:46 PM	37.897638		0.470621
8/27/2012 1:32:51 PM	38.159161		0.473587
8/27/2012 1:32:56 PM	38.42108		0.47691
8/27/2012 1:33:01 PM	38.298111		0.475418
8/27/2012 1:33:06 PM	38.370131		0.476314
8/27/2012 1:33:11 PM	38.370131		0.476314
8/27/2012 1:33:13 PM	38.459779		0.477408
Average	38.104672		0.473026
Std Dev	0.33839308		0.0041882154

How to Record a Data Point


To Record a Data Point, select “Record Data Point” from the main menu or right click a data tab and select “Record Data Point”.




or




The first step in recording data points is to configure what and how to take each point. Use the Settings menu to open the Record Data Point Wizard to define the data points to be taken.




A “Record Data Point Wizard” dialog will appear.



The first page of the Record Data Point Wizard is where the user selects which device items they would like record. On the left hand side is a list of all available devices. Selecting a device will result in the “Device Items” list being updated to reflect the available items for the selected device. To add an item, highlight the desired item or items in the “Device Items” list and click the “Add” button.



Repeat this process until all the desired items are listed on the right side.



Note: The user can invert any selection by right clicking. This will highlight all items that are not currently highlighted and will remove highlight from any items that are currently highlighted.

To remove an item from the list of items to record, the user can either click the “Clear All” button to remove all items or simply select the desired item from the “Items to Record” list and click the “Remove” button.

The "Record Data Point Wizard" dialog box is shown with three main sections: "Devices:", "Device Items:", and "Items to Record:". The "Devices:" list contains "2500" and "Manual Device". The "Device Items:" list contains "Frost Point", "Test Pressure", and "Test Temperature". The "Items to Record:" list contains "%RH@PcTc (2500)", "Frost Point (Manual Device)", "Test Temperature (Manual Device)", and "Test Pressure (Manual Device)". Between the "Device Items:" and "Items to Record:" lists are three buttons: "ADD" (a plus sign), "REMOVE" (a minus sign), and "CLEAR ALL" (a red X). At the bottom are four buttons: "< Back", "Next >", "Finish", and "Cancel".

Once complete, clicking the “Next” button will bring up the Properties page. From the properties page the user can define what will be recorded for each point.

The "Record Data Point Wizard" dialog box is shown with the "Recording Properties" section. It contains a text box "Include the last" followed by a spinner box set to "10" and the text "points". Below this are two checked checkboxes: "Include Average" and "Include Standard Deviation". At the bottom are four buttons: "< Back", "Next >", "Finish", and "Cancel".


Recording Properties define the number of prior (last) points to include and whether to calculate the average and or standard deviation for the defined number of prior points taken with each recorded data point.

Include the last ... points defines the number of prior points to include with each recorded data point. The prior points are determined from the Data Tab for each point being recorded. ControLog will include the number of prior (last) points directly using the entries in the Data Tab for the given point from the time the data point was recorded.

Include Average defines whether to include an average of the defined number of prior points for the recorded data items for each point taken.


Include Standard Deviation defines whether to include the standard deviation of the defined number of prior points for the recorded data items for each point taken.

Once complete, clicking the “Next” button will bring up the frequency page. From the frequency page the user can define when and how to take a data point.



Recording Frequency defines when and how the system will take a data point. The user should select (place check mark) by each method they would like to use.

Manually indicates that the user will manually take points when they want, using either the “Take Point” menu item or by pressing the Ctrl-P key combination on the keyboard.



Automatically with each profile point indicates that a point will be taken automatically at the completion of each profile soak phase.

Note: If the user manually bypasses the soak phase using the manual profile control buttons, an automatic point will NOT be taken regardless of the selection of “Automatically with each profile point” setting.

Automatically after manual device entry indicates that a point will be taken automatically when the user completes a manual device entry.

Once complete, clicking the “Finish” will close the wizard and save the user selection for recording data points.

2500 Data Auto Profile Manual Device Data %RH@PcTc Recorded Data Points						
Recorded Data Points						
Manual Device Entry Point 1	2500	2500	Manual Device	All recorded data points will appear in the Recorded Data Points tab		
Date/Time Stamp	Frost Point [°C]	%RH@PcTc	Frost Point [°C]	Test Pressure [psia]	Test Temperature [°C]	
9/11/2012 1:57:46 PM	-0.46677127	25.000000	9.3924	12.128	20.127	
9/11/2012 1:58:16 PM	-0.47586164	24.995000	9.3924	12.128	20.127	
9/11/2012 1:58:46 PM	-0.47233929	24.995000	9.3924	12.128	20.127	
Average	9.3914635	50.006167	9.3924	12.128	20.127	
Std Dev	0.0017532385	0.0039106692	0	0	0	
Point 1	2500	2500	Manual Device	Manual Device	Manual Device	
Date/Time Stamp	Frost Point [°C]	%RH@PcTc	Frost Point [°C]	Test Pressure [psia]	Test Temperature [°C]	
9/11/2012 1:57:46 PM	-0.46677127	50.0105				
9/11/2012 1:58:16 PM	-0.47586164	50.0051	9.3924	12.128	20.127	
9/11/2012 1:58:46 PM	-0.47233929	50.0051				
Average	9.3906544	50.0069	9.3924	12.128	20.127	
Std Dev	0.001965995	0.0031176915	0	0	0	
Profile Point 1	2500				Manual Device	
Date/Time Stamp	Frost Point [°C]	%RH@PcTc		Test Pressure [psia]	Test Temperature [°C]	
9/11/2012 1:57:46 PM	-0.46677127	25.000000				
9/11/2012 1:58:16 PM	-0.47586164	24.995000				
9/11/2012 1:58:46 PM	-0.47233929	24.995000				
Average	-0.4716574	24.9944	0	0	0	
Std Dev	0.0045833863	0.010877959	0	0	0	
Profile Point 2					Manual Device	
Date/Time Stamp					Test Temperature [°C]	
9/11/2012 2:09:16 PM						
9/11/2012 2:09:46 PM	9.357444	49.988				
9/11/2012 2:10:16 PM	9.3583322	49.9935				

Auto Profiling

The Auto Profiling feature is very similar to the Generate and Purge modes with the main exception that profiling relies on a predefined list of setpoints referred to as a profile. The user configurable profile is used as ControLog's road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next setpoint.

Note: Before attempting to operate the system in an Auto Profile mode, you should become thoroughly familiar with the manual Generate and Purge mode of operation.

The Profile Tab is used to create and modify auto profiles. The tab consist of drop drowns, entry fields and a data grid.

Auto Profile

Control Mode: Dew Point

Test Pressure: 14.700 psia

Test Temperature: 23.000 °C

Mode and Test setpoints

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
1	Warning	-90	0.5	0		Yes
2		-90	0.5	0		No
3	Generate	-80	0.5	0		No
4	Generate	-85	0.5	0	minutes 1 hours	No
*						

Assurance Conditions


Warning

New Point Line


Creating a New Profile

A profile is a list of humidity, temperature, pressure, flow, and time parameters that are used during automated control of the 3900 generator. The profile essentially programs the computer/controller operations.

Open a new profile tab by selecting “New” from the profile menu.




Select the desired **Control Mode** for the entire profile. The Control Mode is what the user would like to generate during the profile. The control modes are the same modes available during manual operation.



For more information, see [Mode Menu](#).

Enter the **Test Pressure** to be used during the entire profile. This value is not required and is ignored if the 3900 test pressure transducer is plugged in. When the test pressure transducer is not plugged in, calculations of Frost Point, Dew Point, and %RH are referenced to this user entered test pressure value rather than a measured value. For accurate humidity generation under these conditions, the absolute pressure at the device under test should be entered as the test pressure.



Enter the **Test Temperature** to be used during the entire profile. This value is not required and is ignored if the 3900 test temperature probe is plugged in. When the test temperature probe is not plugged in, calculations of %RH are referenced to this user entered test temperature value rather than a measured value. For accurate humidity generation under these conditions, the temperature at the device under test should be entered as the test temperature.

The screenshot shows the '3900 Data' window with the 'Auto Profile' tab selected. The 'Auto Profile' sub-window contains the following settings:

- Control Mode:** Frost Point (dropdown menu)
- Test Pressure:** 14.700 psia
- Test Temperature:** 23.0 °C

A text box explains: 'The profile is used as a road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.'

Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
*	[dropdown arrow]					

The first column, next to the point numbers, is the Generator Mode. The generator mode defines the run mode that the 3900 will operate in for this profile point. The generator modes are the same run modes available during manual operation.

Note: Shutdown is only available for the last point.

This screenshot is similar to the previous one, but the 'Generator Mode' dropdown menu for the first point is open, showing the following options:

- Purge
- Generate
- Shutdown

For more information, see [Run Menu](#).

The first column next to the Generator Mode contains the humidity value to generate and is titled at the top according to the currently selected control mode. In the example shown, the profile control mode is set to Dew Point. Whenever the profile control mode is changed, the title on this column changes to reflect the control mode selected.

The screenshot shows the '3900 Data' window with the 'Auto Profile' tab selected. The 'Control Mode' is set to 'Dew Point'. The 'Test Pressure' is 14.700 psia and the 'Test Temperature' is 23.000 °C. A text box explains that the profile is used as a road map during Auto Profile operation. Below this is a table with columns: Point, Generator Mode, Dew Point [°C], Mass Flow Rate [l/m], Ramp Time, Soak Time, and Assurance. The first row is labeled 'Purge' with a value of -90 in the Dew Point column. A '*' symbol is in the first column of the second row.

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
	Purge	-90				
*						

When the Control Mode is set to Saturation Pressure, the first column next to the Generator Mode contains the Saturation Temperature followed by the Saturation Pressure. In this control mode the user can also drop down a Mini version of HumiCalc to help calculate the desired Saturation Temperature or Pressure by clicking the drop down arrow on the desired cell. For example, let's say the user wants to calculate the required Saturation Pressure needed to generate a -80.0 °C Frost Point with a Saturation Temperature of -70.0 °C. Start by entering a -70.0 °C Saturation Temperature and then click the Saturation Pressure HumiCalc drop down at the right of the cell box.

The screenshot shows the '3900 Data' window with the 'Auto Profile' tab selected. The 'Control Mode' is set to 'Saturation Pressure'. The 'Test Pressure' is 14.700 psia and the 'Test Temperature' is 23.000 °C. A text box explains that the profile is used as a road map during Auto Profile operation. Below this is a table with columns: Point, Generator Mode, Saturation Temperature [°C], Saturation Pressure [psia], and Mass Flow Rate [l/m]. The first row is labeled 'Generate' with values -70 and 177.26. The second row is labeled 'Generate' with values -70 and 150. A dropdown arrow is visible next to the 150 value. A callout box labeled 'HumiCalc Drop Down' points to this dropdown arrow.

Point	Generator Mode	Saturation Temperature [°C]	Saturation Pressure [psia]	Mass Flow Rate [l/m]
1	Generate	-70	177.26	1
2	Generate	-70	150	1
*				

Next, select the known to be Frost Point and enter the desired Frost Point temperature of -80.0 °C. Notice that ControLog automatically entered the other known values based on the already entered profile values.

The screenshot shows the HumiCalc dialog box. A callout 'Select the Known' points to the 'Known' dropdown menu which is set to 'Frost Point'. Another callout 'Enter the known values' points to the input fields for 'Frost Point' (-80.00), 'Saturation Temperature' (-70.0), 'Test Pressure' (14.7), and 'Test Temperature' (23.0). A third callout 'Select Calculate' points to the 'Calculate' button.

Saturation Pressure [psia]	Mass Flow Rate [l/m]	Ramp
177.26	1	0
150.00	1	0

HumiCalc

Configuration

Known: Frost Point

Known Values

Frost Point: -80.00

Saturation Temperature: -70.0

Test Pressure: 14.7

Test Temperature: 23.0

Calculate

Clicking the “Calculate” button will result in the calculated Saturation Pressure being placed in the profile cell and the HumiCalc dropdown will close.

The screenshot shows the '3900 Data' window with the 'Auto Profile' tab selected. The 'Control Mode' is 'Saturation Pressure'. The 'Test Pressure' is 14.700 psia and the 'Test Temperature' is 23.000 °C. A text box explains: 'The profile is used as a road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.'

Point	Generator Mode	Saturation Temperature [°C]	Saturation Pressure [psia]	Mass Flow Rate
1	Generate	-70	177.26	1
2	Generate	-70	72.3	1
*				

A callout 'Calculated Saturation Pressure' points to the value 72.3 in the Saturation Pressure column for point 2.

The **Mass Flow Rate** column contains values of the air flow at which the generator will operate. Although not affecting the generated value of humidity, flow rate does affect the air exchange rate, equilibration time of the instruments under test and permeation.

3900 Data Auto Profile

Auto Profile

Control Mode: Dew Point

Test Pressure: 14.700 psia

Test Temperature: 23.000 °C

The profile is used as a road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
	Purge	-90	0.5			
*						

Ramp Time is the desired amount of time the 3900 should take to transition from one profile test point to another. Setting a ramp time of zero instructs the 3900 to make the transition as quickly as possible. Zero is the setting used for most profiles.

Note: The first point should not have a ramp time, because the starting setpoints can vary.

3900 Data Auto Profile

Auto Profile

Control Mode: Dew Point

Test Pressure: 14.700 psia

Test Temperature: 23.000 °C

The profile is used as a road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
	Purge	-90	0.5	0		
*						

Unit dropdown menu options: minutes, milliseconds, seconds, minutes (selected), hours, days.

Soak Time is the desired amount of time to generate at a particular profile point. The soak time required depends on the application, but should be a significant amount of time based upon the humidity measurement devices being calibrated.

The screenshot shows the '3900 Data Auto Profile' window. The 'Control Mode' is set to 'Dew Point'. The 'Test Pressure' is 14.700 psia and the 'Test Temperature' is 23.000 °C. A text box explains: 'The profile is used as a road map during Auto Profile operation. It defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.'

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
	Purge	-90	0.5	0 minute	5	Yes
*						

A dropdown menu is open for the 'Soak Time' column, showing options: milliseconds, seconds, minutes, hours (selected), and days.

Note: Both Ramp Time and Soak Time are limited to a maximum time of 24.8551348032407 days.

Assurance, if set to “Yes”, forces the system to wait until the measured values are within a specified tolerance and stability before ControLog will start the Soak Phase. If “No” is set, the Soak Phase will start immediately upon completion of the Ramp Phase.

This screenshot is identical to the one above, but the dropdown menu for the 'Assurance' column is open, showing 'Yes' (selected) and 'No' options.

Point	Generator Mode	Dew Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
	Purge	-90	0.5	0 minute	5 hours	Yes
*						

If Assurance is set to “Yes”, a small “Add” button will appear on the right hand side. Clicking the “Add” button will open the “Assurance Conditions” dialog for the point. The dialog allows the user to enter various assurance conditions for the profile point.

On the upper left hand side is a list of all available devices. Clicking a device will result in the Device Items list being updated to reflect the available items for the selected device. To add an item, highlight the desired item in the Device Item list, enter the desired Tolerance and/or Standard Deviation and click the “Add” button. For quicker assurance times, increase the Tolerance and/or the Standard Deviation values. Tighter tolerances or standard deviations (smaller values) result in longer assurance times. Setting these values too small could prevent assurance conditions from being met, therefore preventing the system from advancing to the next profile point.

The **Tolerance** is the allowable variation between the setpoint and the actual. This is best thought of as a window based on a minimum and maximum value, the minimum being the setpoint minus the tolerance and the maximum being the setpoint plus the tolerance. Once the actual value is within the window the tolerance portion of the condition is considered met.

Note: Tolerances can only be entered for 3900 setpoints. The field will be grayed out for all other non 3900 setpoint device items. This is because ControLog can only assure a tolerance for an item that it has the ability to control.

Standard Deviation is a statistic used to measure the variation in the actual data and can be thought of as how spread out or stable the data is. ControLog calculates the Standard Deviation from the device data tab for the points within the given **Sample Size**. When the actual standard deviation is less than the defined limit, the standard deviation portion of the condition is considered met.

Note: The sample size should always be carefully considered based on the data storage interval. Too small of a sample size in relation to the data storage interval will result in a small number of points used to calculate the Standard Deviation.

To update an assurance condition, select the desired condition from the list at the bottom, make the desired changes and select the Update button. To delete an assurance condition, select the desired condition from the list at the bottom and select the Delete button.



Point 1 Assurance Conditions

Device:
3900

Device Item:
Mixing Ratio by Volume
Mixing Ratio by Weight
Moist Air Density
Percent by Volume
Percent by Weight
PPMv
PPMv $\pm U$ ($k=2$)
PPMw
PPMw $\pm U$ ($k=2$)
Saturation Pressure
Saturation Pressure $\pm U$ ($k=2$)
Saturation Temperature

3900 Saturation Temperature Assurance Values

Tolerance: ± 0.1
Standard Deviation: 0.1
Sample Size: 15.0 minutes

Device	Item	Tolerance	Standard Deviation	Sample Size
3900	Saturation Temperature	± 0.1	0.1	15 minutes
*				

Ok Cancel

Once all assurance conditions have been completed, select the “Ok” button to close the dialog.


Adding more points to the auto profile is the same process as entering the first point, but the user can let ControLog help fill in values for the new point by simply entering the desired values and then by selecting the new point line (indicated by the * asterisk). ControLog will predict values for any empty field by either copying the values from the point above or by continuing the pattern from the previous points. For example if the previous Frost Points were -90 °C and -80 °C ControLog will automatically use -70 °C for the next point if the user leaves that field empty.

Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
1	Purge	-90	0.5	0 minutes	5 hours	Yes
2	Generate	-90	0.5	0 minutes	1 hours	No
3	Generate	-80	0.5	0 minutes	1 hours	No
4	Generate	-70	0.5	0 minutes	1 hours	No
*						

New points can also be inserted between existing points by right clicking and selecting “Insert Point” from the context menu. ControLog will insert a new point at the selected location and will automatically predict the values.


Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
1	Purge	-90	1	0 minutes	5 hours	Yes
2	Generate	-90	1	0 minutes	1 hours	No
3	Generate	-80	1	0 minutes	1 hours	No
4	Generate	-70	1	0 minutes	1 hours	No
5	Generate	-60	1	0 minutes	1 hours	No
*						

Existing points can also be deleted by selecting the desired point and then by right clicking and selecting “Delete Point” from the context menu.



Profile points that cause operational issues for the 3900 are automatically flagged by ControLog and are indicated by a red circle with an explanation point . An explanation of the issue will display once the user places the mouse cursor over the warning icon.


Warning: The user should address and resolve all issues before attempting to run the Auto Profile.




Saving a Profile

The Profile tab can be saved to file for future recall. ControLog Auto Profiles are saved in XML format with a *.profile extension.

To perform the save, select “Save” from the main menu or right click a profile tab and select “Save”.




Using the “Save Profile” dialog, select the location and name you want to save the profile as.




Opening a Profile

Profiles can be loaded from previously saved profile files.

To perform the open, select “Open” from the main menu or right click a profile tab and select “Open”.



Using the “Open Profile” dialog, browse and select the profile file that you want to open.




Once the load is complete the profile tab will be displayed with the loaded profile points.

Printing a Profile


The profile tab can be printed to any of the PC's installed printers.

Note: You must have a printer installed before you can print using ControlLog.

To perform the print, select "Print" from the main menu or right click a profile tab and select "Print".




Using the "Print" dialog, select the desired printer, range and number of copies you want to print.




Running an Auto Profile


To run an auto profile, select “Run Auto Profile” from the main menu or right click a profile tab and select “Run Auto Profile”.

Note: While operating the system in the Auto Profile mode, manual setpoint and mode changes are not allowed. Also, if changes are made with the 3900 keypad directly, they will quickly be overridden by ControLog.




Selection will open the “Profile Starting Point” dialog which allows the user to select which point in the profile they would like to start the profile on. This feature provides more flexibility by allowing the user to skip ahead to a desired point within the auto profile.




If the profile contained errors, which were indicated by a red circle with an explanation point  on the profile, a warning message will appear when the user tries to run the profile.

Warning: Running a profile with errors may cause adverse behavior during the profile run. The user is strongly encouraged to address and fix all profile issues before attempting to run the profile.




Once the auto profile begins, ControLog will begin sending the commands and setpoints for the starting profile point. The System Timing tab will be displayed in the information tab group and the Auto Profile Controls and Status will appear in the status bar.

The **System Timing** tab gives detailed information on the Auto Profile as it runs. The elapsed and remaining Phase, Point and Total time are listed along with the detailed assurance conditions values and tolerances.



The **Auto Profile Controls and Status** consist of the generator run state for the current profile point, the profile point the system is currently running and the phase of the current point. It also consists of three shortcut buttons to control the operation of the auto profile. There is a “Pause Profile” button, advance to “Next Point” button and advance to the “Next Phase” button. These are the same commands that are available in the profile menus.



Understanding Profile Phases

Each profile point consists of three distinct phases; Ramp, Assurance and Soak. Each phase accomplishes a specified task.


Ramp Phase

The **Ramp Phase** is used to linearly transition from one point to the next point in a given amount of time.

Assurance Phase

The **Assurance Phase** forces the system to wait until measured parameters and setpoint values are within a specified tolerance and/or stability before the computer starts the Soak Phase.

During the assurance phase, assurance conditions that have not been met will be displayed in yellow on the parameter tab of the device containing the condition.



	Setpoint	Actual	U	
Frost Point	-90.000	-34.896	±0.099	°C
Saturation Pressure	11.0	63.465	+0.300	psia
Saturation Temperature				
Test Pressure	14.700	14.700	±0.050	psia
Test Temperature	23.000	23.000	±0.080	°C
Mass Flow Rate	1.0000	2.8601		l/m

The system timing tab provides a detailed view of each assurance condition for the current point. The actual value of each condition will be displayed in red if it has yet to be met and will be displayed in green once the condition has been met. The system timing tab also displays criteria the condition needs to meet for both tolerance and standard deviation.

System Timing

Auto Profile: Generate - Point: 2 of 5 - Phase: Assurance

	<u>Elapsed</u>	<u>Remaining</u>	<u>Assurance Condition</u>	<u>Actual Value</u>	<u>Min Value</u>	<u>Max Value</u>	<u>Std Dev</u>	<u>Std Dev Limit</u>
Phase Time:	00:07:56	00:00:00	Frost Point (3900)	-34.744	-90.100	-89.900	0.0485	0.1000
Point Time:	00:07:56	01:00:00						
Total Time:	00:07:56	04:00:00						

Assurance Condition Status

Status Log

3900 Reported Errors

System Timing

Note: The assurance phase will be active for a minimum of 30 seconds. This delay allows the 3900 to calculate setpoints before ControLog begins to assure each condition.

Soak Phase

The **Soak phase** is the desired amount of time to generate at a particular point before proceeding to the next point.

Example 1

Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
1	Generate	-80	1	0 minutes	1 hours	No

Example 1 causes the Soak phase to begin immediately at the start of the profile point, even though the 3900 may still be adjusting to the point. The next point will start after the 1 hour soak phase.

Example 2

Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
1	Generate	-80	1	0 minutes	1 hours	Yes

Example 2 causes the Assurance phase to begin immediately at the start of the profile point. Measured values are continually compared with the setpoint values until they agree with the set tolerance and/or the measured values are stable to within the specified degree. Once assured (tolerances met) the Soak phase begins. Total time required for this point varies and depends upon the amount of time that is required for assurance (dictated by tolerances). Total time for example 2 is elapsed assurance time plus 1 hour.


Example 3

Point	Generator Mode	Frost Point [°C]	Mass Flow Rate [l/m]	Ramp Time	Soak Time	Assurance
2	Generate	-80	1	0.5 minutes	1 hours	Yes


The Ramp phase (Ramp Timer) begins at the start of the point. The system adjusts slowly toward the setpoint, taking 30 minutes. Once the ramp time has elapsed, assurance starts and waits for tolerances to be met. When tolerances are met, the soak phase begins and lasts 1 hour. Total time for example 3 is 30 minutes, plus elapsed assurance time, plus 1 hour. assurance time is a variable and depends on tolerances.

Manual Override of Profile

Although the system is operating automatically, some manual control is allowed using the Hold Point, Next Phase and Next Point menu items on the Run menu or Status Bar.




or




Holding the Profile


Selecting **Hold Point** from the profile menu or status bar will stop the current "Remaining Time" timers, allowing the system to remain indefinitely at the current point. While in a hold mode, the system is prevented from completing the ramp, assurance, or soak phases of a point.



When holding, a hold indicator appears in the System Timing tab and the hold menu buttons change into resume buttons.




To resume the profile point, select the **Resume Point** from the profile menu or status bar. This re-enables the timing functions and allows the profile to resume normal operation.




Advancing to the Next Point

Selecting **Next Point** from the Run menu manually duplicates the action which automatically occurs when the Point time counter reaches zero.




Advancing to the Next Phase

Selecting **Next Phase** from the profile menu or status bar manually duplicates the action which automatically occurs when a Remaining Ramp or Soak Time counter reaches zero, or when the assurance conditions are met. In other words, it causes Ramp Phase to proceed to the Assurance or Soak Phase, Assurance to proceed to Soak, or Soak to proceed to Ramp of the next profile point. This allows for early manual termination of any phase within a profile point.



Stopping the Auto Profile

Selecting **Stop Auto Profile** from the profile menu or context menu will terminate the profile at the current point and the generator will continue at its current setpoints for Saturation Pressure, Saturation Temperature and Flow. Another way to exit the Auto Profile is to switch from Auto Profile to Generate, Purge or Shutdown.




Connections

ControLog supports a customizable interface that works with most devices. ControLog will allow the user to create a new device connection using the “Connection Wizard” or open previously saved connections. The wizard will open a separate dialog window containing various steps that will guide the user in defining the communication required to receive the desired data items from the device. The user can create as many (up to 60) or as few data items as they see fit for any one device. Each data item can be uniquely named and once connected will be recorded in its own parameter and data tab. ControLog also allows the user to save these interfaces for future use.

***Note:** ControLog has a limit of ten devices connected at any given time.*

The “Connection Wizard” allows the user to step through the connection configuration. Using the “Next” and “Back” buttons the user is allowed to progress through the connection configuration steps. At any time the user may cancel the new connection or opening of a connection by selecting the “Cancel” button. Once the last step has been completed the “Finish” button will be available to complete the new connection.




***Note:** It is always recommended to have the manufactures documentation for the device being connected handy while creating the new connection. It is also recommended to use a terminal based application to test the various commands before creating a new connection.*


Serial Connection

A Serial Connection uses either an RS-232 or RS-485 interface to acquire data from a given ASCII based serial device. The customizable interface provided by ControlLog allows the user to define the ASCII commands that are sent and/or received through the RS-232/485 interface to communicate with the serial device. The system supports both request to receive type of communication as well as receive only type of communication.

To create a new serial connection, select “New” from the Connections menu. This will open a “Connection Wizard” dialog that will step the user through the connection definition process.



Select “Serial” as the type of connection the device uses.



Enter a unique name for the device.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a light gray background. In the center, there is a white rectangular box with the text 'Enter a name for the device' above it. Below the text is a single-line text input field. At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select whether the device requires a setup command or commands. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed.



The image shows the next step of the 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a light gray background. In the center, there is a white rectangular box. At the top of this box is the question 'Do you need to send a setup command or commands to configure this device?'. Below the question are two radio button options: 'Yes' (which is unselected) and 'No' (which is selected with a checkmark). Below these options is a paragraph of text: 'Select 'Yes' if you need to send a setup command or commands to configure the device. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed.' Below this paragraph is another line of text: 'Select 'No' if you don't need to send any setup commands.' At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

If setup commands are required, then enter the ASCII setup command or commands that will be sent at the start of communication. Refer to the legend to enter special characters such as carriage returns and/or line feeds.


Note: All setup commands are case sensitive.

Note: End of Line (EOL) or End of Transmission (EOT) characters such as carriage returns and/or line feeds are very important and are the leading cause to failed communication. Refer to the manufacturer's documentation for the device to verify the required EOL or EOT characters.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area is titled "Setup Command or Commands to send" and contains a large text input field. Below the input field, there is instructional text: "Enter the ASCII setup command or commands that will be sent once at the start of communication. These are commands that are required to configure or setup the device. Refer to the legend to enter special characters such as carriage returns or line feeds." To the right of this text is a "Legend" box containing the following entries: "<CR> = Carriage Return", "<LF> = Line Feed", and "<T> = Tab". At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Enter the number of responses you expect to receive from the device. This is the amount of response messages that you will receive from the device, not necessarily the number of data items. A device response message could contain multiple data items. Later you will select the number of data items that each response message contains.

For example, if you need to send three commands to receive three response messages, then you would enter three. Or if your device automatically sends one message periodically you would only have a single response and you would enter one.



The dialog box is titled "Connection Wizard" and contains the question "How many responses do you expect to receive?". Below the question is a text input field containing the number "1". A detailed instruction block explains that the user should enter the number of response messages expected, providing examples for three commands and periodic messages. At the bottom are four buttons: "< Back", "Next >", "Finish", and "Cancel".


How many responses do you expect to receive?

1

Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains.

< Back Next > Finish Cancel

Select whether a command needs to be sent to request a response. If the device requires a command to be sent to receive a response then select "Yes". If the device automatically outputs data without any request then select "No".



The dialog box is titled "Connection Wizard" and contains the question "Do you need to send a command to receive a response from this device?". Below the question are two radio button options: "Yes" (which is selected) and "No". A detailed instruction block explains the selection criteria. At the bottom are four buttons: "< Back", "Next >", "Finish", and "Cancel".

Do you need to send a command to receive a response from this device?


☒ Yes
☐ No

Select 'Yes' if you need to send a command to request data from the device.
Select 'No' if the device automatically transmits data periodically.

< Back Next > Finish Cancel


If a command was required then enter the ASCII Command or Commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns and/or line feeds.

Note: All commands are case sensitive.




The screenshot shows the 'Connection Wizard' window with the title bar 'Connection Wizard' and a close button. The main area is titled 'Command or Commands to send' and contains a large text input field. Below the input field, there is instructional text: 'Enter the ASCII command or commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns or line feeds.' To the right of the input field is a 'Legend' box containing the following entries: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. At the bottom of the window are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Select the ASCII identification that will be used to identify the response if no command is required but the device automatically sends messages periodically. The identification is the unique characters at the beginning of the message that identifies the response received. This is often used on devices that automatically transmit multiple data messages without the use of a command to request data.



The screenshot shows the 'Connection Wizard' window with the title bar 'Connection Wizard' and a close button. The main area is titled 'Response Identification' and contains a text input field. Below the input field, there is instructional text: 'Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.' At the bottom of the window are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Select whether the response has a delimiter that is separating each data item. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then each item is separated by a comma delimiter.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'Are the response data items separated by a delimiter?' is displayed in bold. Below this, there are two radio button options: 'Yes' (which is unselected) and 'No' (which is selected with a checkmark). A text box below the options provides instructions: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.' and 'Select 'No' if the data items are not separated by a unique character.' At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Are the response data items separated by a delimiter?

☐ Yes


☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

Enter the number of data items the response message contains. A data item is the numeric value or flag portion for a single device parameter within the response message. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'How many data items does the response contain?' is displayed in bold. Below this, there is a numeric spinner control showing the value '1'. A text box below the spinner provides instructions: 'Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.' At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

How many data items does the response contain?


1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

Enter a name or description for the data item and define the syntax, unit and scaling. The data item syntax is defined using the symbols in the Legend. Use the “X” symbol to indicate a character that should be ignored, use the “N” symbol to represent a numeric ASCII character and the “A” symbol to indicate a flag or any ASCII character. This dialog will repeat for each data item in the response.

Note: The syntax cannot contain both a Flag and a Numeric syntax definition. If the user requires both, then create another data item to define them separately.



The Connection Wizard dialog box is used to define a data item. It includes fields for the Data Item Name, Data Item Syntax, and Data Item Unit. A legend on the right defines the symbols used in the syntax: X for Ignore, N for Numeric, and A for Flag. A scaling table allows for defining signal and data values. A large text area is provided for the Response Syntax. Navigation buttons at the bottom include < Back, Next >, Finish, and Cancel.

Connection Wizard

Data Item Name

Data Item Syntax

Data Item Unit
 None ▾

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

X = Ignore
 N = Numeric
 A = Flag

Response Syntax

< Back Next > Finish Cancel

When defining a numeric syntax enter an “N” for each possible digit in the response. For example, if you know the device returns a six digit numeric value you would enter “NNN.NNN”. The decimal point is not required and its location is not important. Decimal points, plus signs and minus signs are treated the same as an “N” and are allowed merely to help make the syntax resemble a number value.

Note: It is important to have sufficient numeric definition to assure all possible numeric responses will be covered, especially when a device responses with scientific notation or varying precession.

Connection Wizard

Data Item Name
Data Point

Data Item Syntax
NNN.NNNN

Data Item Unit
None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax
NNN.NNNN

< Back Next > Finish Cancel

It is possible to scale a numeric data item response. The scaling consists of a two point definition for a linear scaling or a three to seven point definition for polynomial interpolation scaling. The number of points determines the degree of the polynomial used to scale the data item response. Lagrangian Interpolation is used to determine the coefficients of the polynomial. Each point definition consists of a signal value and a data value. The signal value represents the “raw” output signal from the device. The data value represents the actual value or real world value at the given signal value.

Scaling allows the user to scale a numeric data item response into a given humidity value. For example, if you have a numeric data item response that ranges from -1 to 1 and it is known that 0 corresponds to -100 and 1 corresponds to 100. The user can then enter these scaling values and ControlLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Connection Wizard

Data Item Name
Data Point

Data Item Syntax
NNN.NNNN

Data Item Unit
None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value
-1	-100
1	100

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax
NNN.NNNN

< Back Next > Finish Cancel

Selecting the unit for the data item will allow ControlLog to convert the value to the selected system units for display in the parameter tab and record the value in the default SI units in the data tab. Remember this is the unit the device is sending the data item in, not the unit you wish to display the data item as. If “None” is selected then ControlLog will treat the data item as a simple number and will display and record the value exactly as it is received.

Connection Wizard

Data Item Name
Data Point

Data Item Syntax
NNN.NNNN

Data Item Unit
Temperature °C

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value
-1	-100
1	100

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax
NNN.NNNN

< Back Next > Finish Cancel

When defining flag type syntax enter an “A” for each character in the response that represents the flag. The Flag Definitions define what each possible ASCII flag represents. The user must enter a numeric value for each flag definition which will be recorded in the data tab and a description for the flag that will be shown in the parameters tab.

Connection Wizard

Data Item Name: Stable

Data Item Syntax: XXXA

Data Item Unit: None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Flag Definitions

Flag	Numeric Value	Description
*	0	No
S	1	Yes

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax

XXXA

< Back Next > Finish Cancel

Enter the End of Transmit (EOT) character that is sent after the response. This is the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns and/or line feeds.

Connection Wizard

End of Transmit (EOT) character sent after the response

Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend

<CR> = Carriage Return
<LF> = Line Feed
<T> = Tab

< Back Next > Finish Cancel

Select the name and location to save the new serial connection. Selecting the “Browse” button will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. All device connection files are saved in XML format with a (*.device) extension.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area contains a section titled "Save Device Configuration as" with a text input field and a "Browse" button. Below this, a message reads: "Select the name and location to save this device configuration to." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".


Next, the user can select whether to connect to the device now or to exit without connecting.

Note: The user can connect at any time by loading the device from the Connections menu.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area contains a section titled "Would you like to connect to this device now?" with two radio button options: "Yes" (which is selected) and "No". Below these options, a message reads: "Select 'Yes' if you would like to connect to the device now." followed by "Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Select the communication port that the device is connected to and select the port settings for the device.

Warning: *Incorrect RTS and DTR settings can prohibit communication. Refer to the device manufacturer's documentation to verify all device port settings.*



The image shows a 'Connection Wizard' dialog box with a title bar and a close button. It is divided into two main sections. The left section, titled 'Select the port the device is connected to', contains a label 'Connect using' and a dropdown menu currently showing 'COM1'. The right section, titled 'Select the port settings for the device', contains five dropdown menus: 'Baud Rate' (9600), 'Data Bits' (8), 'Parity' (None), 'Stop Bits' (1), and 'Handshake' (None). Below these are two checkboxes, 'Enable RTS' and 'Enable DTR', both of which are unchecked. At the bottom of the dialog are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select the port the device is connected to		Select the port settings for the device	
Connect using	COM1	Baud Rate	9600
		Data Bits	8
		Parity	None
		Stop Bits	1
		Handshake	None
		<input type="checkbox"/> Enable RTS	<input type="checkbox"/> Enable DTR

< Back Next > Finish Cancel

Select the access rate at which ControLog will communicate with the device.

Note: It is always recommended to start with the default 1.5 second access interval and to modify later as needed.

Warning: Do not set the access interval too small. If the device is not capable of communicating at the set interval then ControLog may inadvertently think communication has been lost when the device does not reply within the desired amount of time.



Connection Wizard

Access interval to the device

1.5 seconds

This is the rate at which ControLog will communicate with the device.

☐ Show Interface Console


< Back Next > **Finish** Cancel

Select “Show Interface Console” to automatically open a Device Interface Console tab once communication has been established with the device.


Serial Connection Example 1

This example will demonstrate the creation of a serial connection to an RH Systems® 373 Dew Point Mirror. We will request the Frost Point temperature and Atmospheric Pressure from the mirror as data items.


Start by selecting “New” from the Connections menu.



Select “Serial” as the type of device connection.



Enter “373 Dew Point Mirror” as the name for the device.




The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Enter a name for the device' and a text input field containing '373 Dew Point Mirror'. At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

For the 373 we do not need to send any setup commands.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Do you need to send a setup command or commands to configure this device?'. Below the heading are two radio button options: 'Yes' (unselected) and 'No' (selected). Below the options is a paragraph of text: 'Select 'Yes' if you need to send a setup command or commands to configure the device. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed. Select 'No' if you don't need to send any setup commands.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

The 373 will send two responses; one with the requested Frost Point value and the other with the requested Atmospheric Pressure value.




Connection Wizard

How many responses do you expect to receive?

Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains.

< Back **Next >** Finish Cancel

The 373 requires a request command to be sent to receive either the Frost Point value or the Atmospheric Pressure value.



Connection Wizard

Do you need to send a command to receive a response from this device?

☒ **Yes**
☐ **No**

Select 'Yes' if you need to send a command to request data from the device.
Select 'No' if the device automatically transmits data periodically.

< Back **Next >** Finish Cancel

The 373 documentation tells us we need to send the frost point command to receive the Frost Point value. The documentation also states the command must be terminated with a carriage return or a carriage return line feed combination.

The screenshot shows the 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area is titled 'Command or Commands to send'. Below this title is a text input field containing 'FP?<CR>'. To the left of the input field, there is a small text block: 'Enter the ASCII command or commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns or line feeds.' To the right of the input field is a 'Legend' box containing the following text: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. At the bottom of the dialog box are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

The 373 will respond to the frost point command with a single numeric value and therefore no delimiter is used.

The screenshot shows the 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area is titled 'Are the response data items separated by a delimiter?'. Below this title are two radio button options: 'Yes' (unchecked) and 'No' (checked). Below these options is a text block: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom of the dialog box are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

The response from the frost point command will only contain a single numeric value.

Connection Wizard

How many data items does the response contain?

1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

The first item is given the name Frost Point. The 373 can respond with a high precision numeric value so the Data Item Syntax was defined with the maximum digits possible that the 373 can respond with. The Data Item's Unit was defined as a temperature in degrees Celsius as stated in the documentation for the 373. The response value is the actual value so there is no need to define any type of scaling.

Connection Wizard

Data Item Name
FrostPoint

Data Item Syntax
NNNNNNNNNNNNNNNNNN

Data Item Unit
Temperature °C

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax
NNNNNNNNNNNNNNNNNN

< Back Next > Finish Cancel


According to the 373 documentation, each response is terminated with a carriage return and line feed.

The screenshot shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'End of Transmit (EOT) character sent after the response'. Below this heading is a text input field containing '<CR><LF>'. To the left of the input field is a paragraph: 'Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.' To the right of the input field is a 'Legend' box containing three entries: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. At the bottom of the dialog are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Next we define the command for the second response. Referring to the 373 documentation we enter the command for atmospheric pressure. Again the command must be terminated with a carriage return or carriage return line feed combination.

The screenshot shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Command or Commands to send'. Below this heading is a text input field containing 'Pa?<CR>'. To the left of the input field is a paragraph: 'Enter the ASCII command or commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns or line feeds.' To the right of the input field is a 'Legend' box containing three entries: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. At the bottom of the dialog are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

The 373 will respond to the atmospheric pressure command with a single numeric value and therefore no delimiter is used.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'Are the response data items separated by a delimiter?' is displayed in bold. Below it are two radio button options: 'Yes' (unselected) and 'No' (selected). A text box provides instructions: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Are the response data items separated by a delimiter?

☐ Yes


☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

The response from the atmospheric pressure command will only contain a single numeric value.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'How many data items does the response contain?' is displayed in bold. Below it is a spin box containing the number '1'. A text box provides instructions: 'Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

How many data items does the response contain?

1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

The second item is given the name Pressure. Again the Data Item Syntax was defined with the maximum digits possible that the 373 can respond with. The Data Item's Unit was defined as a pressure in Pascal's as stated in the documentation. The response value is the actual value so there is no need to define any type of scaling.

Connection Wizard

Data Item Name: Pressure

Data Item Syntax: NNNNNNNNNNNNNNNNN

Data Item Unit: Pressure (dropdown), Pa (dropdown)

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax

NNNNNNNNNNNNNNNN

< Back Next > Finish Cancel

The response is terminated with a carriage return and line feed.

Connection Wizard

End of Transmit (EOT) character sent after the response

<CR><LF>

Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend

- <CR> = Carriage Return
- <LF> = Line Feed
- <T> = Tab

< Back Next > Finish Cancel


Save the newly created device to a file so that it can be recalled at a later time.



Select to connect to the device now.



Select the port the 373 is connected to. Notice that ControLog indicates which ports are in use. Refer to the 373 documentation for specific port settings:




The Connection Wizard dialog box is shown with the title bar "Connection Wizard". It contains two main sections:

- Select the port the device is connected to:** A dropdown menu labeled "Connect using" shows a list of ports: COM1, COM2 (in use), COM3, COM4, and COM5. COM1 is currently selected.
- Select the port settings for the device:** A series of dropdown menus for Baud Rate (9600), Data Bits (8), Parity (None), Stop Bits (1), and Handshake (None). Below these are two checkboxes: "Enable RTS" and "Enable DTR", both of which are unchecked.

At the bottom of the dialog are four buttons: "< Back", "Next >", "Finish", and "Cancel". The "Next >" button is highlighted in yellow.

The default access interval of 1.5 seconds is entered.




The Connection Wizard dialog box is shown with the title bar "Connection Wizard". It contains a single section:


- Access interval to the device:** A text input field contains the value "1.5", and a dropdown menu next to it is set to "seconds". Below this, a note states: "This is the rate at which ControLog will communicate with the device."

At the bottom of the dialog are four buttons: "< Back", "Next >", "Finish", and "Cancel". The "Finish" button is highlighted in yellow.

Once completed, ControlLog will attempt to establish communication with the 373.




Once communication is successfully established with the 373, a new parameter tab and data tab will be created. Note that both the parameter tab and data tab have the two data items we defined. Notice that since we defined what unit the data items were received in, ControlLog is able to convert the values into the selected units for easy reading.




Serial Connection Example 2

This example will demonstrate the creation of a serial connection to an MBW[®] DP-30 Precision Dew Point Hygrometer. This example will be working with a DP-30 that has the temperature and pressure option. The DP-30 does not require any request to receive data. Instead, it constantly outputs three data messages at a periodic rate.


Start by selecting “New” from the Connections menu.



Select “Serial” as the type of device connection.



Enter “DP 30 with Options” as the name for the device.




The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Enter a name for the device' and a text input field containing 'DP 30 with Options'. At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

For the DP-30 we do not need to send any setup commands.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Do you need to send a setup command or commands to configure this device?'. Below the heading are two radio button options: 'Yes' (unselected) and 'No' (selected). Below the options is a paragraph of text: 'Select 'Yes' if you need to send a setup command or commands to configure the device. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed. Select 'No' if you don't need to send any setup commands.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

We will receive three responses from the DP-30; one with the Dew Point Temperature value and Status, one with the Dry Temperature value and one with the Pressure value.



Connection Wizard


How many responses do you expect to receive?

3

Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains.

< Back **Next >** Finish Cancel

The DP-30 does not require any command to be sent to receive a response. The DP-30 automatically transmits the messages.



Connection Wizard

Do you need to send a command to receive a response from this device?

☐ Yes


☒ **No**

Select 'Yes' if you need to send a command to request data from the device.

Select 'No' if the device automatically transmits data periodically.

< Back **Next >** Finish Cancel

For the first message we enter the response identification of the Dew Point Temperature message.



The screenshot shows the 'Connection Wizard' window with the 'Response Identification' step. A text box contains 'DP'. Below it, explanatory text states: 'Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.' At the bottom are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Connection Wizard


Response Identification

DP

Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.

< Back Next > Finish Cancel

The DP-30 Dew Point Temperature message is of a fixed format and does not use any delimiter.



The screenshot shows the 'Connection Wizard' window with the question 'Are the response data items separated by a delimiter?'. There are two radio buttons: 'Yes' (unchecked) and 'No' (checked). Below the buttons, explanatory text states: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Connection Wizard

Are the response data items separated by a delimiter?

☐ Yes

☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

The Dew Point Temperature message contains three flags and one value that results in a total of four data items for this response.

How many data items does the response contain?

4

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

The first data item is a flag indicating whether the DP-30 is balanced. The DP-30 will send a “*” indicating that it is not balanced or a “B” indicating that it is balanced. Since we are not interested in the next two characters after this flag in the message we will tell ControLog to ignore them by adding two “X”s after the flag symbol. Next we define these flag definitions and give them a simple “Yes” or “No” description to allow the user to easily determine their state in the parameter tab.

Data Item Name
Balanced

Data Item Syntax
AXX

Data Item Unit
None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Flag	Numeric Value	Description
*	0	No
B	1	Yes

Legend
X = Ignore
N = Numeric
A = Flag

Response Syntax
 DPAXX

< Back Next > Finish Cancel

The second data item is a flag indicating whether the DP-30 is running. The DP-30 will send a “*” indicating that it is not running or an “R” indicating that it is running. We again define these flag definitions and give them a simple “Yes” or “No” description to allow the user to easily determine the state in the parameter tab.

Connection Wizard

Data Item Name: Run

Data Item Syntax: A

Data Item Unit: None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Flag Definitions

Flag	Numeric Value	Description
*	0	No
R	1	Yes

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax
DPAXXA

< Back Next > Finish Cancel

The third data item is a flag indicating whether the DP-30 is in standby. The DP-30 will send a “*” indicating that it is not in standby or an “S” indicating that it is in standby. We again define these flag definitions and give them a simple “Yes” or “No” description to allow the user to easily determine the state in the parameter tab.

Connection Wizard

Data Item Name: Standby

Data Item Syntax: A

Data Item Unit: None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Flag Definitions

Flag	Numeric Value	Description
*	0	No
S	1	Yes

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax
DPAXXA

< Back Next > Finish Cancel

The fourth data item is the dew point temperature value. The DP-30 will send a numeric value indicating the dew point temperature in degrees Celsius. We define the syntax exactly as the DP-30 documentation specifies. The response value is the actual value so there is no need to define any type of scaling.

The screenshot shows the 'Connection Wizard' dialog box. It has three input fields at the top: 'Data Item Name' with the text 'Dew Point', 'Data Item Syntax' with the text '+NNN.NN', and 'Data Item Unit' with a dropdown menu showing 'Temperature' and '°C'. Below these is a text area with instructions: 'Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.' To the right of this text is a 'Scaling' table with two columns: 'Signal Value' and 'Data Value'. The table is currently empty. To the right of the table is a 'Legend' box with the following text: 'X = Ignore', 'N = Numeric', and 'A = Flag'. Below the instructions and table is a 'Response Syntax' field containing the text 'DPAXXAA+NNN.NN'. At the bottom of the dialog are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Signal Value	Data Value

Response Syntax
DPAXXAA+NNN.NN

< Back Next > Finish Cancel

Notice that ControlLog shows the combined response syntax for this message at the bottom of the form. This syntax will closely resemble the syntax defined in the DP-30 documentation for the dew point temperature message. The bold portion of the syntax is the current data item's syntax within the response message.

This is a close-up of the 'Response Syntax' field from the previous screenshot. It shows the text 'DPAXXAA+NNN.NN' where 'DPAXXAA' is in a standard font and '+NNN.NN' is in a bold font.

Response Syntax
DPAXXAA+NNN.NN


The Dew Point Temperature message is terminated with a carriage return.

The screenshot shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'End of Transmit (EOT) character sent after the response'. Below this heading is a text input field containing '<CR>'. To the right of the input field is a 'Legend' box containing the following text: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. Below the input field and legend is a paragraph of text: 'Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.' At the bottom of the dialog are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Next we define the response identification for the second message. Referring to the DP-30 documentation we enter the identification for the Dry Temperature message.

The screenshot shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Response Identification'. Below this heading is a text input field containing 'DT'. Below the input field is a paragraph of text: 'Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.' At the bottom of the dialog are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

The DP-30 Dry Temperature message is of a fixed format and does not use any delimiter.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main text asks: 'Are the response data items separated by a delimiter?'. There are two radio buttons: 'Yes' (unchecked) and 'No' (checked). Below the radio buttons, there is explanatory text: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom, there are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Are the response data items separated by a delimiter?

☐ Yes


☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

The Dry Temperature message contains only one numeric value therefore there is only one data item for this response.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main text asks: 'How many data items does the response contain?'. There is a spin box with the value '1'. Below the spin box, there is explanatory text: 'Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.' At the bottom, there are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

How many data items does the response contain?

1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

The message is lead by a number of ASCII characters that we do not have any interest in. We will tell ControLog to ignore these leading characters by defining the syntax with five “X” characters then the numeric syntax. We define the numeric syntax exactly as the DP-30 documentation specifies. Again, the response value is the actual value so there is no need to define any type of scaling.

Connection Wizard

Data Item Name: Temperature

Data Item Syntax: XXXXX+NNN.NN

Data Item Unit: Temperature °C

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax

DTXXXXX+NNN.NN

< Back Next > Finish Cancel

The Dry Temperature message is terminated with a carriage return.

Connection Wizard

End of Transmit (EOT) character sent after the response

<CR>


Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend

- <CR> = Carriage Return
- <LF> = Line Feed
- <T> = Tab

< Back Next > Finish Cancel

We now define the response identification for the third and last message. Referring to the DP-30 documentation we enter the identification for the Pressure message.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main area is titled 'Response Identification'. It features a text input field containing 'PR'. Below the field is a paragraph of text: 'Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).


Response Identification

PR

Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.

< Back Next > Finish Cancel

The DP-30 Pressure message is of a fixed format and does not use any delimiter.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main area is titled 'Are the response data items separated by a delimiter?'. It contains two radio button options: 'Yes' (unchecked) and 'No' (checked). Below these is a paragraph of text: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Are the response data items separated by a delimiter?

☐ Yes

☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

The Pressure message contains only one numeric value therefore there is only one data item for this response.

Connection Wizard

How many data items does the response contain?

1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

Again the message is lead by a number of ASCII characters that we do not have any interest in. We will tell ControLog to ignore these leading characters by defining the syntax with five “X” characters then the numeric syntax. We define the numeric syntax exactly as the DP-30 documentation specifies. Note that we set the data item unit to the unit specified in the DP-30 documentation. This allows ControLog to convert the reading into a more desired unit if needed. Again, the response value is the actual value so there is no need to define any type of scaling.

Connection Wizard

Data Item Name
Pressure

Data Item Syntax
XXXXX+NN.NNN

Data Item Unit
Pressure Pa

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value


Legend

- X = Ignore
- N = Numeric
- A = Flag


Response Syntax
PRXXXXX+NN.NNN

< Back Next > Finish Cancel

The Pressure message is terminated with a carriage return.




Save the newly created device to a file so that it can be recalled at a later time.



Select to connect to the device now.




Select the port the DP-30 is connected to. Notice that ControLog indicates which ports are in use. Refer to the DP-30 documentation for the RS-232 port settings:




Use the default access rate of 1.5 seconds.



Once completed, ControLog will attempt to establish communication with the DP-30.




Once communication is successfully established with the DP-30, a new parameter tab and data tab will be created. Notice that ControLog automatically translates the ASCII flags based on the definitions described when creating the data items for the connection.




Analog Connection

An Analog Connection uses the Agilent® 34970A Data Acquisition/Switch Unit to acquire data from a single or multiple analog devices. The customizable interface provided by ControLog allows the user to define different analog types and scales to read various analog signals. Refer to the Agilent® documentation for more information on connecting analog devices to the Data Acquisition Unit.

To create a new analog connection, select “New” from the Connections menu. This will open a “Connection Wizard” dialog that will step the user through the connection definition process.



Select “Analog” as the type of connection the device uses.




Enter a unique name for the analog device or devices.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard' with a close button. The main area has a light gray background. In the center, there is a white rectangular box with the text 'Enter a name for the device' in bold. Below this text is a single-line text input field. At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Enter the number of analog signals you want to monitor. This is the combined amount of analog signals from each device you wish to monitor. Each analog signal connected to the 34970A will represent a single data item parameter.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard' with a close button. The main area has a light gray background. In the center, there is a white rectangular box with the text 'How many analog signals do you want to monitor?' in bold. Below this text is a spin box containing the number '1'. Underneath the spin box, there is a smaller text box containing the text: 'Enter the number of analog signals that you want to monitor. Each analog signal connected to the 34970A will represent a single data item parameter.' At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. This dialog will repeat for each monitored signal.

Connection Wizard

Data Item Name

Channel 101

Function Volts DC

Range 300 V

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit None

Scaling

Signal Value	Data Value

< Back Next > Finish Cancel

The function defines what type of analog signal the device uses.

Connection Wizard

Data Item Name

Channel 101

Function Volts DC

Range 300 V

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit None

Scaling

Signal Value	Data Value

< Back Next > Finish Cancel

The range specifies the range for which the analog signal operates within. The available range selection changes automatically to reflect what is available for the selected function. For example, the Resistance function has a Range from 100 Ohms to 100 Mega Ohms.

Connection Wizard

Data Item Name

Channel 101

Function Resistance

Range 100 MΩ

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit None

Scaling

Signal Value	Data Value

< Back Next > Finish Cancel

Selecting the unit for the data item will allow ControLog to convert the value to the selected system units for display in the parameter tab and record the value in the default SI units in the data tab. Remember this is the unit the device is sending the data item in, not the unit you wish to display the data item as. If “None” is selected then ControLog will treat the data item as a simple number and will display and record the value exactly as it is received.

Connection Wizard

Data Item Name

Channel 101

Function Volts DC

Range 10 V

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit

Temperature °C

Temperature °F

Pressure °C

Enthalpy K

Density

FlowRate

None

Scaling

Signal Value	Data Value

< Back Next > Finish Cancel

It is possible to scale an analog signal. The scaling consists of a two point definition for a linear scaling or a three to seven point definition for polynomial interpolation scaling. The number of points determines the degree of the polynomial used to scale the analog signal. Lagrangian Interpolation is used to determine the coefficients of the polynomial. Each point definition consists of a signal value and a data value. The signal value represents the “raw” analog signal. The data value represents the actual value or real world value at the given signal value.

Scaling allows the user to scale an analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0 %RH and 5 volts corresponds to 100 %RH. The user can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Connection Wizard

Data Item Name

Channel 101

Function Volts DC

Range 10 V

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit None

Scaling

Signal Value	Data Value
0	0
5	100

< Back **Next >** Finish Cancel

Select the name and location to save the new analog connection. Clicking the “Browse” button will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. All device connection files are saved in XML format with a (*.device) extension.



The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Save Device Configuration as". It contains a text input field and a "Browse" button. Below the input field, it says "Select the name and location to save this device configuration to." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Next the user can select whether to connect to the device now or to exit without connecting.


Note: The user can connect at any time by loading the device from the Connections menu.



The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Would you like to connect to this device now?". It contains two radio button options: "Yes" (which is selected) and "No". Below the options, it says "Select 'Yes' if you would like to connect to the device now." and "Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Select the communication port that the Agilent® 34970A Data Acquisition/Switch Unit is connected to and select the port settings for the 34970A.

Note: ControLog defaults to the default port settings for the Agilent® 34970A Data Acquisition/Switch Unit. Refer to the Agilent® documentation for instructions on viewing or setting the Data Acquisition RS-232 settings.



The image shows a software configuration window titled "Agilent 34970A Configuration". It is divided into two main sections. The left section, titled "Select the port the device is connected to", contains a label "Connect using" followed by a dropdown menu currently set to "COM1". The right section, titled "Select the port settings for the device", contains five dropdown menus: "Baud Rate" (57600), "Data Bits" (8), "Parity" (None), "Stop Bits" (1), and "Handshake" (None). Below these are two checkboxes, "Enable RTS" and "Enable DTR", both of which are unchecked. At the bottom of the window are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Section	Setting	Value
Select the port the device is connected to	Connect using	COM1
Select the port settings for the device	Baud Rate	57600
	Data Bits	8
	Parity	None
	Stop Bits	1
	Handshake	None
	Enable RTS	<input type="checkbox"/>
	Enable DTR	<input type="checkbox"/>

Navigation buttons: < Back, Next >, Finish, Cancel

Select the access rate to communicate with the Agilent® 34970A Data Acquisition/Switch Unit.

Note: Because of the nature of the 34970A and how it operates it is critical to specify an access interval that is long enough to allow the 34970A to complete its scan list within the specified interval. As a general rule of thumb, use 1.5 seconds per every 10 signals connected to the Data Acquisition Unit with a minimum access interval of 1.5 seconds.




The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. Inside the dialog, there is a section titled 'Access interval to the device'. Below this title, there is a text input field containing the value '1.5' and a dropdown menu set to 'seconds'. Below the input field, there is a descriptive text: 'This is the rate at which ControLog will communicate with the device.' At the bottom right of the dialog, there is a checkbox labeled 'Show Interface Console'. At the bottom of the dialog, there are four buttons: '< Back', 'Next >', 'Finish' (highlighted in yellow), and 'Cancel'.

Select “Show Interface Console” to automatically open a Device Interface Console tab once communication has been established with the device.


Analog Connection Example

This example will demonstrate the creation of an analog connection. This example will be working with a -10 to +10V input signal that we will scale to a Dew Point/Frost Point Temperature, a Temperature Thermistor and a 4-wire Temperature RTD.


Start by selecting “New” from the Connections menu.



Select “Analog” as the type of device connection.




Enter “Analog Dew Point Device” as the name for the device.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Enter a name for the device' and a text input field containing 'Analog Dew Point Device'. At the bottom, there are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

In this example we have three signals that we want to monitor; a voltage, thermistor and an RTD.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'How many analog signals do you want to monitor?' and a spin box containing the number '3'. Below the spin box, there is a text box with the following text: 'Enter the number of analog signals that you want to monitor. Each analog signal connected to the 34970A will represent a single data item parameter.' At the bottom, there are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

The first item is given the name “Dew/Frost Point”. We select the Channel number that the signal is connected to and since the signal is a voltage we select “Volts DC” as the function. The range is set to “10V” and we know the signal corresponds to a temperature in degrees Celsius. We also scale the value since we know -10V corresponds to a Frost Point temperature of -100 °C and +10V corresponds to a Dew Point temperature of +100 °C. By entering scaling ControLog will automatically scale the signal for display in the parameter tab and when recorded in the data tab.

The screenshot shows the 'Connection Wizard' dialog box. The 'Data Item Name' field contains 'Dew/Frost Point'. The 'Channel' dropdown is set to '109'. The 'Function' dropdown is set to 'Volts DC'. The 'Range' dropdown is set to '10 V'. The 'Data Item Unit' dropdown is set to 'Temperature' and the unit dropdown is set to '°C'. The 'Scaling' table is as follows:

Signal Value	Data Value
-10	-100
10	100

At the bottom, there are buttons for '< Back', 'Next >', 'Finish', and 'Cancel'.

The second item is given the name “Temperature 1”. We select the Channel number that the signal is connected to. Notice ControLog indicates which channels have already been configured and are in use.

The screenshot shows the 'Connection Wizard' dialog box for the second item. The 'Data Item Name' field contains 'Temperature 1'. The 'Channel' dropdown is set to '110'. The 'Function' dropdown is set to 'Volts DC'. The 'Range' dropdown is set to '10 V'. The 'Data Item Unit' dropdown is set to 'None'. The 'Scaling' table is empty. At the bottom, there are buttons for 'Next >', 'Finish', and 'Cancel'. A list of channels is visible on the left side of the dialog, with '109 (in use)' highlighted, indicating that channel 109 has already been configured.

The function is set to “Temperature Thermistor” and we set the Range to “10k” based on the type of Thermistor we are connecting. Notice that ControLog automatically selects temperature as the unit type but allows the user to select the desired temperature unit. Since the signal value is the actual temperature value, no scaling is required for this data item.

Connection Wizard

Data Item Name: Temperature 1

Channel: 102

Function: Temperature Thermistor

Range: 10k

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit: Temperature °C

Scaling:

Signal Value	Data Value

< Back Next > Finish Cancel

The third item is given the name “Temperature 2”. The function is set to “Temperature RTD (4 Wire)” and we set the Range to “100 ohms, 85” based on the type of RTD we connected. Next, select the Channel number the RTD is connected to. Notice the channel list is smaller because 4-wire signals are automatically paired with the selected channel plus 10 to provide the source and sense connections for an RTD. Given this, the RTD in this example will occupy both channel 110 and 120 to complete its 4-wire connection.

Connection Wizard

Data Item Name: Temperature 2

Channel: 110

Function: Temperature RTD (4 Wire)

Range: 100 ohms, 91

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit: Temperature °C

Scaling:

Signal Value	Data Value

Next > Finish Cancel

Again, notice that ControLog automatically selects temperature as the unit type but allows the user to select the desired temperature unit. Since the signal value is the actual temperature value, no scaling is required for this data item.

Connection Wizard

Data Item Name: Temperature 2
Channel: 110
Function: Temperature RTD (4 Wire)
Range: 100 ohms, 91

Enter a name or description of the data item and define the analog channel for it. The channel is the corresponding instrument channel that the device is connected to. For 4-wire measurements, the instrument automatically pairs channel n with channel n+10 to provide the source and sense connections. The function defines what type of analog signal the device uses. The range specifies the range for which the analog signal operates within. The scaling allows you to scale the analog signal into a given humidity value. For example, if you have an analog signal that ranges from 0 volts to 5 volts and it is known that 0 volts corresponds to 0% RH and 5 volts corresponds to 100%RH. You can then enter these scaling values and ControLog will automatically apply the scaling to the data item whenever it is displayed or logged.

Data Item Unit: Temperature °C

Scaling

Signal Value	Data Value

< Back Next > Finish Cancel

Save the newly created device to a file so it can be recalled at a later time.

Connection Wizard

Save Device Configuration as

Analog Dew Point Device Browse

Select the name and location to save this device configuration to.


< Back Next > Finish Cancel

Select to connect to the device now.



The image shows a 'Connection Wizard' dialog box. It has a title bar with a close button. The main content area contains the question 'Would you like to connect to this device now?' followed by two radio button options: 'Yes' (which is selected) and 'No'. Below these options, there is explanatory text: 'Select 'Yes' if you would like to connect to the device now.' and 'Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu.' At the bottom of the dialog, there are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Select the communication port that the Agilent® 34970A Data Acquisition/Switch Unit is connected to and select the port settings for the 34970A.




The image shows a 'Connection Wizard' dialog box with two main sections. The left section, titled 'Select the port the device is connected to', contains a label 'Connect using' and a dropdown menu currently showing 'COM1'. The right section, titled 'Select the port settings for the device', contains five dropdown menus: 'Baud Rate' (57600), 'Data Bits' (8), 'Parity' (None), 'Stop Bits' (1), and 'Handshake' (None). Below these are two checkboxes: 'Enable RTS' and 'Enable DTR', both of which are unchecked. At the bottom, there are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Enter a sufficient access interval for the amount of signals being monitored. In this example we have less than 10 signals so we can start with the minimum access interval of 1.5 seconds.



Once completed, ControLog will attempt to establish communication with the Agilent® 34970A Data Acquisition/Switch Unit.



Once communication is successfully established with the 34970A a new parameter tab and data tab will be created. Notice that ControLog automatically scales the voltage signal based on the definitions described when defining the data item.

3900 ControLog

File Mode Run Profile Graph Data Units Connections Help

Analog Dew Point Device Parameters

Actual

Temperature 1	25.391	°C
Dew/Frost Point	-86.616	°C
Temperature 2	23.508	°C

Analog Signal Data Items

New Analog Dew Point Device Parameter Tab

New Analog Dew Point Device Data Tab

Analog Dew Point Device Data 3900 Data

Date/Time Stamp	Temperature 1 [°C]	Dew/Frost Point [°C]	Temperature 2 [°C]
	25.372	-86.616266	23.475
	25.377	-86.616239	23.494
	25.378	-86.616266	23.497
	25.378	-86.616198	23.499

ControlLog automatically scales the signal

Status Log

- 11/4/2009 9:46:55 AM: Communication established with the Analog Dew Point Device device
- 11/4/2009 9:47:24 AM: Communication established with the 3900 generator
- 11/4/2009 9:47:24 AM: System is shutdown
- 11/4/2009 9:47:27 AM: 3900 Test Pressure transducer is disconnected
- 11/4/2009 9:47:28 AM: 3900 Test Temperature probe is disconnected

Analog Dew Point Device Parameters 3... Status Log


3900 Test Temperature probe is disconnected Point Time: 00:00:00

GPIB Connection


GPIB Connections use the National Instruments® GPIB-USB-HS USB-to-488.2 converter to acquire data from a given GPIB device. The customizable interface provided by ControLog allows the user to define the ASCII commands that are sent and/or received through the IEEE 488.2 interface to communicate with the GPIB device. The system supports both a request to receive type of communication as well as a receive only type of communication.

Note: The setup dialog for a GPIB device is virtually identical to the setup dialog for a Serial device. The only difference being the interface used to transmit and receive the data.

To create a new GPIB connection, select “New” from the Connections menu. This will open a “Connection Wizard” dialog that will step the user through the connection definition process.



Select “General Purpose Interface” as the type of connection the device uses.



Enter a unique name for the device.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard' with a close button. The main area has a light gray background. In the center, there is a white rectangular box with the text 'Enter a name for the device' above it. Below the text is a single-line text input field. At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select whether the device requires a setup command or commands. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard' with a close button. The main area has a light gray background. In the center, there is a white rectangular box. At the top of this box is the question 'Do you need to send a setup command or commands to configure this device?'. Below the question are two radio button options: 'Yes' (which is unselected) and 'No' (which is selected with a checkmark). Below these options is a block of text: 'Select 'Yes' if you need to send a setup command or commands to configure the device. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed.' Below this text is another line of text: 'Select 'No' if you don't need to send any setup commands.' At the bottom of the dialog, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

If setup commands are required then enter the ASCII setup command or commands that will be sent at the start of communication. Refer to the legend to enter special characters such as carriage returns and/or line feeds.


Note: All setup commands are case sensitive.

Note: End of Line (EOL) or End of Transmission (EOT) characters such as carriage returns and/or line feeds are very important and are the leading cause to failed communication. Refer to the manufacturer's documentation for the device to verify the required EOL or EOT characters.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area is titled "Setup Command or Commands to send" and contains a large text input field. Below the input field, there is instructional text: "Enter the ASCII setup command or commands that will be sent once at the start of communication. These are commands that are required to configure or setup the device. Refer to the legend to enter special characters such as carriage returns or line feeds." To the right of this text is a "Legend" box containing the following entries: "<CR> = Carriage Return", "<LF> = Line Feed", and "<T> = Tab". At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Enter the number of responses you expect to receive from the device. This is the amount of response messages you will receive from the device, not necessarily the number of data items. A device response message could contain multiple data items. Later you will select the number of data items that each response message contains.

For example, if you need to send one command to receive one response message then you would enter one. Or if your device automatically sends three messages periodically you would have three responses and you would enter three.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main content area has a question: 'How many responses do you expect to receive?'. Below the question is a text box containing the number '1'. Below the text box is a paragraph of explanatory text: 'Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains.' At the bottom of the dialog box are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.


How many responses do you expect to receive?

1

Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains.

< Back Next > Finish Cancel

Select whether a command needs to be sent to request a response. If the device requires a command to be sent to receive a response then select "Yes". If the device automatically outputs data without any request then select "No".



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main content area has a question: 'Do you need to send a command to receive a response from this device?'. Below the question are two radio button options: 'Yes' (which is selected) and 'No'. Below the radio buttons is a paragraph of explanatory text: 'Select 'Yes' if you need to send a command to request data from the device. Select 'No' if the device automatically transmits data periodically.' At the bottom of the dialog box are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Do you need to send a command to receive a response from this device?

☒ Yes
☐ No

Select 'Yes' if you need to send a command to request data from the device.
Select 'No' if the device automatically transmits data periodically.

< Back Next > Finish Cancel

If a command was required, then enter the ASCII Command or Commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns and/or line feeds.


Note: All commands are case sensitive.

The screenshot shows the 'Connection Wizard' dialog box with the title bar 'Connection Wizard' and a close button. The main area is titled 'Command or Commands to send' and contains a large text input field. Below the input field, there is instructional text: 'Enter the ASCII command or commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns or line feeds.' To the right of the input field is a 'Legend' box containing the following entries: '<CR> = Carriage Return', '<LF> = Line Feed', and '<T> = Tab'. At the bottom of the dialog are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Select the ASCII identification that will be used to identify the response if no command is required but the device automatically sends messages periodically. The identification is the unique characters at the beginning of the message that identifies the response received. This is often used on devices that automatically transmit multiple data messages without the use of a command to request the data.

The screenshot shows the 'Connection Wizard' dialog box with the title bar 'Connection Wizard' and a close button. The main area is titled 'Response Identification' and contains a text input field. Below the input field, there is instructional text: 'Select the ASCII identification that will be used to identify this response. The identification is the unique characters at the beginning of the response that identifies the data string received from the device. This is often used on devices that automatically transmit multiple data strings without the use of a command to request the data.' At the bottom of the dialog are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Select whether the response has a delimiter that is separating each data item. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then each item is separated by a comma delimiter.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'Are the response data items separated by a delimiter?' is displayed in bold. Below it are two radio button options: 'Yes' (unselected) and 'No' (selected). A text box below the options provides instructions: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Are the response data items separated by a delimiter?

☐ Yes


☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

Enter the number of data items the response message contains. A data item is the numeric value or flag portion for a single device parameter within the response message. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.



The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main content area has a light gray background. At the top, the question 'How many data items does the response contain?' is displayed in bold. Below it is a numeric spinner control showing the value '1'. A text box below the spinner provides instructions: 'Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

How many data items does the response contain?


1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

Enter a name or description for the data item and define the syntax, unit and scaling. The data item syntax is defined using the symbols in the Legend. Use the “X” symbol to indicate a character should be ignored, use the “N” symbol to represent a numeric ASCII character and the “A” symbol to indicate a flag or any ASCII character. This dialog will repeat for each data item in the response.

Note: The syntax cannot contain both a Flag and a Numeric syntax definition. If the user requires both, then create another data item to define them separately.



The Connection Wizard dialog box is used to define a data item. It includes fields for the Data Item Name, Data Item Syntax, and Data Item Unit. A Scaling table allows for defining signal and data values. A Legend explains the syntax symbols: X for Ignore, N for Numeric, and A for Flag. A large text area is provided for the Response Syntax. Navigation buttons at the bottom include < Back, Next >, Finish, and Cancel.

Connection Wizard

Data Item Name

Data Item Syntax

Data Item Unit
None ▼

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax

< Back **Next >** Finish Cancel

When defining a numeric syntax enter an “N” for each possible digit in the response. For example, if you know the device returns a six digit numeric value you would enter “NNN.NNN”. The decimal point is not required and its location is not important. Decimal points, plus signs and minus signs are treated the same as an “N” and are allowed merely to help make the syntax resemble a number value.

Note: It is important to have sufficient numeric definition to assure all possible numeric responses will be covered, especially when a device responses with scientific notation or varying precession.

Connection Wizard

Data Item Name
Data Point

Data Item Syntax
NNN.NNNN

Data Item Unit
None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Scaling

Signal Value	Data Value

Legend

X = Ignore
N = Numeric
A = Flag

Response Syntax
NNN.NNNN

< Back Next > Finish Cancel

It is possible to scale a numeric data item response. The scaling consists of a two point definition for a linear scaling or a three to seven point definition for polynomial interpolation scaling. The number of points determines the degree of the polynomial used to scale the data item response. Lagrangian Interpolation is used to determine the coefficients of the polynomial. Each point definition consists of a signal value and a data value. The signal value represents the “raw” output signal from the device. The data value represents the actual value or real world value at the given signal value.

Scaling allows the user to scale a numeric data item response into a given humidity value. For example, if you have a numeric data item response that ranges from -1 to 1 and it is known that 0 corresponds to -100 and 1 corresponds to 100. The user can then enter these scaling values and ControlLog will automatically apply the scaling to the data item whenever it is displayed or logged.

The screenshot shows the 'Connection Wizard' dialog box. The 'Data Item Name' is 'Data Point', 'Data Item Syntax' is 'NNN.NNNN', and 'Data Item Unit' is 'None'. The 'Scaling' section contains a table with 'Signal Value' and 'Data Value' columns. The 'Response Syntax' is 'NNN.NNNN'.

Signal Value	Data Value
-1	-100
1	100

Legend
X = Ignore
N = Numeric
A = Flag

Response Syntax
NNN.NNNN

Buttons: < Back, Next >, Finish, Cancel

Selecting the unit for the data item will allow ControlLog to convert the value to the selected system units for display in the parameter tab and record the value in default SI units in the data tab. Remember this is the unit that the device is sending the data item in, not the unit you wish to display the data item as. If “None” is selected then ControlLog will treat the data item as a simple number and will display and record the value exactly as it is received.

The screenshot shows the 'Connection Wizard' dialog box with the 'Data Item Unit' dropdown menu open. The 'Data Item Unit' is set to 'Temperature', and the unit is '°C'. The 'Scaling' section contains a table with 'Signal Value' and 'Data Value' columns. The 'Response Syntax' is 'NNN.NNNN'.

Signal Value	Data Value
-1	-1
1	100

Legend
X = Ignore
N = Numeric
A = Flag

Response Syntax
NNN.NNNN

Buttons: < Back, Next >, Finish, Cancel

When defining flag type syntax enter an “A” for each character in the response that represents the flag. The Flag Definitions define what each possible ASCII flag represents. The user must enter a numeric value for each flag definition which will be recorded in the data tab and a description for the flag that will be shown in the parameters tab.

Connection Wizard

Data Item Name
Stable

Data Item Syntax
XXXXA

Data Item Unit
None

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Flag Definitions

Flag	Numeric Value	Description
*	0	No
S	1	Yes

Legend
 X = Ignore
 N = Numeric
 A = Flag

Response Syntax
 XXXXA

< Back Next > Finish Cancel

Enter the End of Transmit (EOT) character that is sent after the response. This is the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns and/or line feeds.

Connection Wizard

End of Transmit (EOT) character sent after the response

Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend
 <CR> = Carriage Return
 <LF> = Line Feed
 <T> = Tab

< Back Next > Finish Cancel

Select the name and location to save the new GPIB connection. Clicking the “Browse” button will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. All device connection files are saved in XML format with a (*.device) extension.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Save Device Configuration as". It contains a text input field and a "Browse" button. Below the input field, the text reads: "Select the name and location to save this device configuration to." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Next, the user selects whether to connect to the device now or to exit without connecting.

Note: The user can connect at any time by loading the device from the Connections menu.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Would you like to connect to this device now?". It contains two radio button options: "Yes" (which is selected) and "No". Below these options, the text reads: "Select 'Yes' if you would like to connect to the device now." and "Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Select the GPIB Board ID and addresses for the device. Refer to the device's documentation for further information on these GPIB settings.



The image shows a Windows-style dialog box titled "Connection Wizard". Inside the dialog, there is a section titled "Select the GPIB settings for the device". This section contains three labels with corresponding input fields: "Board ID" with a numeric spinner set to 0, "Primary Address" with a numeric spinner set to 9, and "Secondary Address" with a dropdown menu set to "None". At the bottom of the dialog, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Connection Wizard

Select the GPIB settings for the device

Board ID 0

Primary Address 9

Secondary Address None

< Back Next > Finish Cancel

Select the access rate at which ControLog will communicate with the device.

Note: It is always recommended to start with the default 1.5 second access interval and to modify later as needed.

Warning: Do not set the access interval to small. If the device is not capable of communicating at the set interval then ControLog may inadvertently think communication has been lost when the device does not reply within the desired amount of time.




The screenshot shows a 'Connection Wizard' window. Inside, there's a section titled 'Access interval to the device'. It contains a text input field with '1.5' and a dropdown menu set to 'seconds'. Below this, a message states: 'This is the rate at which ControLog will communicate with the device.' At the bottom right of this section is a checkbox labeled 'Show Interface Console'. At the very bottom of the wizard are four buttons: '< Back', 'Next >', 'Finish' (highlighted in yellow), and 'Cancel'.

Select “Show Interface Console” to automatically open a Device Interface Console tab once communication has been established with the device.


GPIB Connection Example

This example will demonstrate the creation of a GPIB connection to an Agilent® 34401A Digital Multimeter (DMM) used to read the analog output of an MBW® K-1806. We will request the voltage reading from the DMM and scale it accordingly for the K-1806's Dew Point Temperature.

Start by selecting “New” from the Connections menu.



Select “General Purpose Interface Bus” as the type of device connection.



Enter “K-1806 (HP34401A)” as the name for the device.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Enter a name for the device' and a text input field containing 'K-1806 (HP34401A)'. At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

The DMM requires special setup commands before we can begin to take readings



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area has a heading 'Do you need to send a setup command or commands to configure this device?'. Below the heading are two radio button options: 'Yes' (selected) and 'No'. Below the options is a paragraph of text: 'Select 'Yes' if you need to send a setup command or commands to configure the device. Setup commands are only sent once at the start of communication. These commands are only required if you need to send special commands to configure the device before data request and response commands are processed. Select 'No' if you don't need to send any setup commands.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).


There are three setup commands to configure the DMM that we need to include. The first setup command resets the DMM, the second command clears the DMM and the third command configures the DMM to read DC voltage. Based on the DMM documentation each command must be terminated by a Line Feed.

The dialog box is titled "Connection Wizard" and has a close button in the top right corner. The main heading is "Setup Command or Commands to send". Below this is a text area containing the following commands: `*RST<LF>`, `*CLS<LF>`, and `CONFigure:VOLTage:DC<LF>`. To the right of the text area is a "Legend" box with the following entries: `<CR>` = Carriage Return, `<LF>` = Line Feed, and `<T>` = Tab. Below the text area is a paragraph of instructions: "Enter the ASCII setup command or commands that will be sent once at the start of communication. These are commands that are required to configure or setup the device. Refer to the legend to enter special characters such as carriage returns or line feeds." At the bottom of the dialog are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

We expect to only receive a single response message from the DMM.

The dialog box is titled "Connection Wizard" and has a close button in the top right corner. The main heading is "How many responses do you expect to receive?". Below this is a spin box containing the number "1". Below the spin box is a paragraph of instructions: "Enter the number of response messages you expect to receive from the device. For example if you need to send three commands to receive three response messages, then you would have three responses. Another example would be if your device automatically sends one message periodically, you would only have a single response. Remember this is only the amount of response messages that you receive from the device, not the amount of data items. Later you will select the number of data items that each response contains." At the bottom of the dialog are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

We need to request the DMM to take a voltage reading.



Connection Wizard

Do you need to send a command to receive a response from this device?

☒ **Yes**


☐ **No**

Select 'Yes' if you need to send a command to request data from the device.

Select 'No' if the device automatically transmits data periodically.

< Back Next > Finish Cancel

We will use the read command to request the voltage reading from the DMM. Again we must terminate the command with a Line Feed.



Connection Wizard

Command or Commands to send

Read?<LF>


Enter the ASCII command or commands to request the desired data from the device. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend

- <CR> = Carriage Return
- <LF> = Line Feed
- <T> = Tab

< Back Next > Finish Cancel

The DMM response is not separated by any delimiter.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main text asks: 'Are the response data items separated by a delimiter?'. There are two radio buttons: 'Yes' (unchecked) and 'No' (checked). Below the radio buttons, there is explanatory text: 'Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter. Select 'No' if the data items are not separated by a unique character.' At the bottom, there are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

Are the response data items separated by a delimiter?

☐ Yes


☒ No

Select 'Yes' if each data item is separated by a unique character. For example, if the response contains Dew Point, Frost Point and Pressure with a comma between each data item then you would enter a comma as the delimiter.

Select 'No' if the data items are not separated by a unique character.

< Back Next > Finish Cancel

The DMM response will contain a single voltage reading.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main text asks: 'How many data items does the response contain?'. There is a spin box with the value '1'. Below the spin box, there is explanatory text: 'Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.' At the bottom, there are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

How many data items does the response contain?

1

Enter the number of data items this response message contains. A data item is the numeric value or flag for a single device parameter. For example, if the response contains Dew Point, Frost Point and Pressure then you would have three data items. If the response only contains a single numeric value or flag then you would only have one data item for this response.

< Back Next > Finish Cancel

We will name the data item “Dew/Frost Point” based on the signal we are reading from the K-1806. The Data Item Syntax will be of numeric format. The Data Item’s Unit is defined as a Temperature in degrees Celsius based on the K-1806 documentation. We will apply a scaling to allow ControLog to scale the signal into terms of Dew/Frost Point. Minus one volt represents -100 °C Frost Point and plus one volt represents +100 °C Dew Point.

Connection Wizard

Data Item Name: Dew/Frost Point

Data Item Syntax: +N.NNNNNNNNNNNNN

Data Item Unit: Temperature °C

Enter a name or description of the data item and define the syntax for the response using the coded symbols defined in the legend. Note that the syntax cannot contain both a Flag and Numeric syntax definition. If you require both, then create another data item to define them separately.

Signal Value	Data Value

Legend

- X = Ignore
- N = Numeric
- A = Flag

Response Syntax

+N.NNNNNNNNNNNNN

< Back Next > Finish Cancel

The DMM documentation states the response will be terminated with a Line Feed.

Connection Wizard

End of Transmit (EOT) character sent after the response

<LF>

Select the ASCII character that is sent at the end of transmission of the response. Refer to the legend to enter special characters such as carriage returns or line feeds.

Legend

- <CR> = Carriage Return
- <LF> = Line Feed
- <T> = Tab

< Back Next > Finish Cancel

Save the newly created device to a file so that it can be recalled at a later time.



The screenshot shows a 'Connection Wizard' window with a title bar containing a close button. The main content area is titled 'Save Device Configuration as'. It features a text input field containing 'K-1806 (HP34401A)' and a 'Browse' button to its right. Below the input field, a message reads: 'Select the name and location to save this device configuration to.' At the bottom of the window, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select to connect to the device now.




The screenshot shows the same 'Connection Wizard' window, now at a different step. The title bar remains the same. The main content area is titled 'Would you like to connect to this device now?'. It contains two radio button options: 'Yes' (which is selected, indicated by a checkmark in the box) and 'No' (which is not selected). Below these options, a message reads: 'Select 'Yes' if you would like to connect to the device now. Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu.' At the bottom, the same four buttons are present: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select the GPIB Board ID and addresses for the DMM. Refer to the DMM's documentation for further information on these GPIB settings.




The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main area is titled 'Select the GPIB settings for the device'. It contains three settings: 'Board ID' with a numeric spinner set to 0, 'Primary Address' with a numeric spinner set to 22, and 'Secondary Address' with a dropdown menu set to 'None'. At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select the default access interval.




The image shows a 'Connection Wizard' dialog box with a title bar containing a close button. The main area is titled 'Access interval to the device'. It contains a numeric spinner set to 1.5 and a dropdown menu set to 'seconds'. Below these is the text: 'This is the rate at which ControLog will communicate with the device.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (disabled), 'Finish' (highlighted in yellow), and 'Cancel' (disabled).

Once completed, ControlLog will attempt to establish communication with the DMM.




Once communication is successfully established with the DMM a new parameter tab and data tab will be created. Notice ControlLog automatically scales the voltage signal into a Dew/Frost Point temperature value based on the definitions described when defining the data item.




Manual Connection

A Manual Connection allows the user to manually record data items for a device that either has no interface or has an interface that is not supported by ControLog. Manual devices still have their own parameter and data tab but the data values are manually entered by the user. When the user wants to record values they simply click on the actual field for the data item in the parameter tab and enter the value. Once all data items for the device have been entered, ControLog will record the values into the data tab for the device.

To create a manual connection, select “New” from the Connections menu. This will open a “Connection Wizard” dialog that will step the user through the connection definition process.



Select “Manual” as the type of connection the device uses.




Enter a unique name for the device.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area contains the text 'Enter a name for the device' above a single-line text input field. At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Select the number of manual data entries.



The image shows a 'Connection Wizard' dialog box. The title bar says 'Connection Wizard'. The main area contains the text 'How many data items do you want to manually enter data for?' above a spin box containing the number '1'. Below the spin box, there is a paragraph of text: 'Enter the number of data items that you want to manually enter data for. You will be able to manually enter point data for each data item and ControLog will record your entered data points just as it does for any automated device.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

Connection Wizard

Data Item Name

Data Item Unit

None

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

< Back Next > Finish Cancel

The user only selects the type of unit because all manual entries are entered in the currently selected system unit. For example, if the user creates a manual data item that is a temperature and has the system units set to degrees Celsius, then the user will enter manual values in degrees Celsius. If the system units are set to degrees Fahrenheit then the user will be required to enter manual values in degrees Fahrenheit.

Connection Wizard

Data Item Name

Data Item Unit

None

Temperature

Pressure

Enthalpy

Density

FlowRate

None

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

< Back Next > Finish Cancel

Select the name and location to save the new manual connection. Clicking the “Browse” button will open a save file dialog that will allow the user to specify the name and browse to the desired location to save the file. All device connection files are saved in XML format with a (*.device) extension.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Save Device Configuration as". It contains a text input field with the text "Manual Device" and a "Browse" button to its right. Below the input field, there is a line of text: "Select the name and location to save this device configuration to." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

Last, the user selects whether to connect to the device now or to exit without connecting.


Note: The user can connect at any time by loading the device from the Connections menu.

The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main content area is titled "Would you like to connect to this device now?". It contains two radio button options: "Yes" (which is selected with a checkmark) and "No". Below these options, there is a line of text: "Select 'Yes' if you would like to connect to the device now." and another line: "Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu." At the bottom of the window, there are four buttons: "< Back", "Next >", "Finish" (highlighted in yellow), and "Cancel".


Manual Connection Example

This example will demonstrate the creation of a Manual connection that will consist of three data items: Frost Point, Test Pressure and Test Temperature.

Start by selecting “New” from the Connections menu.



Select “Manual” as the type of connection the device uses.




Enter “Manual Device” as the name for the device.




The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area contains a text box with the label "Enter a name for the device" above it. The text box contains the text "Manual Device". At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

We have three data items that we want to manually enter for this device.



The screenshot shows a window titled "Connection Wizard" with a close button in the top right corner. The main area contains a text box with the label "How many data items do you want to manually enter data for?" above it. The text box contains the number "3". Below the text box, there is a paragraph of text: "Enter the number of data items that you want to manually enter data for. You will be able to manually enter point data for each data item and ControLog will record your entered data points just as it does for any automated device." At the bottom of the window, there are four buttons: "< Back", "Next >" (highlighted in yellow), "Finish", and "Cancel".

The first point is Frost Point and it will be of the temperature unit type.



The image shows a 'Connection Wizard' dialog box. It has a title bar with a close button. The main area contains two input fields: 'Data Item Name' with the text 'Frost Point' and 'Data Item Unit' with a dropdown menu showing 'Temperature'. Below these fields is a text instruction: 'Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Connection Wizard


Data Item Name
Frost Point

Data Item Unit
Temperature ▼

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

< Back **Next >** Finish Cancel

The second point is Test Pressure and it will be of the pressure unit type.



The image shows a 'Connection Wizard' dialog box. It has a title bar with a close button. The main area contains two input fields: 'Data Item Name' with the text 'Test Pressure' and 'Data Item Unit' with a dropdown menu showing 'Pressure'. Below these fields is a text instruction: 'Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.' At the bottom, there are four buttons: '< Back' (disabled), 'Next >' (highlighted in yellow), 'Finish' (disabled), and 'Cancel' (disabled).

Connection Wizard


Data Item Name
Test Pressure

Data Item Unit
Pressure ▼

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

< Back **Next >** Finish Cancel

The third point is Test Temperature and it will be of the temperature unit type.



The image shows a 'Connection Wizard' dialog box. It has a title bar with a close button. The main area contains two input fields: 'Data Item Name' with the text 'Test Temperature' and 'Data Item Unit' with a dropdown menu showing 'Temperature'. Below these fields is a descriptive text: 'Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.' At the bottom, there are four buttons: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Connection Wizard

Data Item Name
Test Temperature

Data Item Unit
Temperature

Enter a name or description of the manual device data item and specify the type of unit that you will be entering the data in.

< Back Next > Finish Cancel

Save the newly created device to a file so that it can be recalled at a later time.



The image shows the same 'Connection Wizard' dialog box, but at a different step. The title bar and close button are the same. The main area now features a section titled 'Save Device Configuration as' which contains a text input field with 'Manual Device' and a 'Browse' button. Below this is the text: 'Select the name and location to save this device configuration to.' At the bottom, the same four buttons are present: '< Back', 'Next >' (highlighted in yellow), 'Finish', and 'Cancel'.

Connection Wizard

Save Device Configuration as

Manual Device Browse

Select the name and location to save this device configuration to.

< Back Next > Finish Cancel

Select to connect to the device now.



The image shows a 'Connection Wizard' dialog box. It has a title bar with the text 'Connection Wizard' and a close button (X) in the top right corner. The main content area is light gray and contains a white rectangular box with a black border. Inside this box, the text 'Would you like to connect to this device now?' is displayed in bold. Below this text are two radio button options: 'Yes' (which is selected, indicated by a checkmark in the box) and 'No'. Below the radio buttons, there is explanatory text: 'Select 'Yes' if you would like to connect to the device now.' and 'Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu.' At the bottom of the dialog box, outside the white box, are four buttons: '< Back', 'Next >', 'Finish' (highlighted in yellow), and 'Cancel'.

Connection Wizard

Would you like to connect to this device now?

☒ **Yes**


☐ **No**

Select 'Yes' if you would like to connect to the device now.

Select 'No' if you would like to exit without connecting to the device. Note you can connect at any time using by loading the device file from the Connections menu.


< Back Next > **Finish** Cancel

Once completed a new parameter tab and data tab will be created.




To manually enter a value, click on the actual field you would like to enter. A Manual Entry box will appear and the title of the manual item being entered will be underlined. For example, to enter a manual value for the Frost Point item click on the Frost Point Actual field.

Note: The manual entries will not be recorded in the data tab until all data item values have been manually entered.




Enter the manual value into the Manual Entry box and select Ok.



A small dialog box with a text input field containing the value "-80.1235". Below the input field are two buttons: "Cancel" on the left and "Ok" on the right.


The manual entry will now be displayed in the actual field of the parameter tab.



A dialog box titled "Manual Device Parameters" with a dropdown arrow and a help icon. It contains a section labeled "Actual" with three rows of parameters:

Actual		
Frost Point	-80.124	°C
Test Pressure		psia
Test Temperature		°C


Repeat the process for the remaining data items.



A dialog box titled "Manual Device Parameters" with a dropdown arrow and a help icon. It contains a section labeled "Actual" with three rows of parameters:

Actual		
Frost Point	-80.124	°C
Test Pressure	14.700	psia
Test Temperature	23.500	°C

Once all data items have been manually entered, ControLog will record the values and clear out the Actual fields in preparation for the next set of manual entries.



Glossary of Terms

%RH

See Percent Relative Humidity

Absolute Humidity

Absolute Humidity is the weight of the water vapor per unit volume of humidified gas.

Access Interval

The access interval is the rate to which ControLog sends and receives commands to and from a specific device.

ASCII

American Standard Code for Information Interchange or ASCII is a set of codes that represent text characters.

Assurance Condition

An Assurance Condition forces an Auto Profile to wait until the measured values are within a specified tolerance and stability before the profile will start the Soak Phase.

Auto Profile

An Auto Profile is used as a road map to automatically control the 3900 generator. The profile defines which setpoint values to go to, at what rate to go from one setpoint to another, and how long to stay at a specific setpoint before moving to the next one.

Dew Point

Dew Point Temperature is the temperature to which a gas must be cooled in order to just begin condensing water vapor in the form of dew. Generally, Dew Point exists at temperatures above freezing. In many instances, Dew Point may actually exist at indicated values below freezing (super-cooled dew). However, it is important to note that Dew Point is not the same as Frost Point. Dew Point is independent of test chamber temperature.

Dew Point Control Mode

The Dew Point control mode, Td, is controlled at a constant value by varying saturation pressure, Ps, to compensate for any changes in either saturation temperature, Ts, or test pressure, Pt. While Dew Point is held constant, other humidity parameters may vary. While in Dew Point control mode, the saturation temperature setpoint is automatically

determined. Dew Point control mode is valid both above and below 0 °C, and Dew Point is independent of test temperature.

DMM

Digital Multimeter.

Dry Air Density

Dry Air Density is the partial density in weight per unit volume of only the dry air portion of a moist air sample. In other words, if the water vapor were removed from a fixed volume of air, the remaining dry air would exhibit this density.

Dry Air Mole Fraction

Dry Air Mole Fraction is the mole fraction of the dry air portion of a sample. The dry air portion is considered to be all constituents in a gas exclusive of the water vapor.

DTR

Data Terminal Ready (DTR) is a signal indicating that the serial device is ready for transmission.

Enhancement Factor

The enhancement factor corrects for the non-ideal behavior of a gas admixed with water vapor. The enhancement factor is a function of two independent variables, pressure and temperature.

Enthalpy

Enthalpy is a measure of the amount of energy required to change a gas from one temperature/humidity value to another. In application, enthalpy is not used as an absolute value, but rather it is the difference in enthalpy between two distinct points which are of interest. The datum point which results in zero enthalpy was therefore arbitrarily chosen at a test temperature of 0 °C and 0 %RH. Applying enthalpy is a matter of computing the difference in enthalpy between two or more distinct data points.

F@Td.Pt

Enhancement Factor at Dew/Frost Point Temperature and Test Pressure.

F@Ts.Ps

Enhancement Factor at Saturation Temperature and Pressure.

F@Tt.Pt

Enhancement Factor at Test Temperature and Pressure.

Frost Point

Frost Point Temperature is the temperature to which a gas must be cooled in order to just begin condensing water vapor in the form of frost or ice, and therefore only exists at values below 0.01 °C. When operating the system with indicated Frost Points above 0.01 °C, the values indicated are to be interpreted as Dew Points. However, Frost Point is not the same as Dew Point for values below freezing. Frost Point is independent of test chamber temperature.

Frost Point Control Mode

The Frost Point control mode, T_f , is controlled at a constant value by varying the saturation pressure, P_s , to compensate for changes in either saturation temperature, T_s , or test pressure, P_t . While Frost Point is held constant other humidity parameters may vary. While in Frost Point control mode, the saturation temperature setpoint is automatically determined. Frost Point is independent of test temperature.

GPIB

General Purpose Interface Bus or GPIB also known as IEEE-488 is a digital communications bus specification.

Grains/lb

Grains per pound is a ratio of the weight, in grains, of water vapor to the weight, in pounds, of the other constituents in the gas. (7000 grains = 1 pound). Once established, Grains/lb is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

HumiCalc

HumiCalc® is a form-based, high-accuracy humidity calculator by Thunder Scientific Corporation. The software application is designed using high accuracy equations to compute complex humidity conversions and uncertainty. The user may select known values to be used in solving for twenty-six humidity parameters. The user may also select the dimensional units, whether or not to use enhancement factors, vapor pressure over ice or water and ITS-90/IPTS-68 temperature scales.

Mixing Ratio by Volume

Mixing Ratio by Volume is a ratio of the partial pressure of the water vapor to the partial pressure of the remaining constituents in the sample. Mixing Ratio by Volume is independent of test chamber temperature.

Mixing Ratio by Weight

Mixing Ratio by Weight is a ratio of the weight of the water vapor to the weight of the remaining constituents in the sample. Mixing Ratio by Weight is independent of test chamber temperature.

Moist Air Density

Moist Air Density is the total weight per unit volume of a moist air sample. This density includes both the weight of the air and the weight of the water vapor.

Percent by Volume

Percent by Volume is a ratio (expressed as a percentage) of the partial pressure of the water vapor to the total pressure of the sample. Percent by Volume is independent of test chamber temperature.

Percent by Weight

Percent by Weight is a ratio (expressed as a percentage) of the weight of the water vapor to the total weight of the sample. Percent by Weight is independent of test chamber temperature.

Percent Relative Humidity

Percent Relative Humidity (%RH) is the ratio of the amount of water vapor in a given sample to the maximum amount possible at the same temperature and pressure.

Percent Relative Humidity Control Mode

The Percent Relative Humidity (%RH) control mode, %RH, is controlled at a constant value by varying saturation pressure, P_s , to compensate for any changes in saturation temperature, T_s , test temperature, T_t , or test pressure, P_t . While %RH is held constant, all other humidity parameters may vary. While in %RH control mode, the saturation temperature setpoint is automatically determined.

Permeation

Permeation refers to a continuous influx from or outgas to the humidity of the surrounding environment (such as the air within the laboratory) through small leaks or semi-permeable surfaces through the walls, fittings, valves, and dead spaces within the system.

PPMv

Parts Per Million by Volume is a ratio of the number of molecules of water vapor to the number of molecules of the other constituents in the gas. Once established, PPMv is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

PPMv Control Mode

The PPMv control mode is controlled at a constant value by varying saturation pressure, P_s , to compensate for any changes in saturation temperature, T_s . While PPMv is held constant, other humidity parameters may vary. While in PPMv control mode, the saturation temperature setpoint is automatically determined. PPMv is independent of test pressure and test temperature.

PPMw

Parts Per Million by Weight is a ratio of the weight of the water vapor in a sample to the weight of the remaining constituents in the gas. Once established, PPMw is pressure and temperature insensitive, and is therefore independent of test chamber temperature and test chamber pressure.

PPMw Control Mode

The PPMw control mode is controlled at a constant value by varying saturation pressure, P_s , to compensate for any changes in saturation temperature, T_s . While PPMw is held constant, other humidity parameters may vary. While in PPMw control mode, the saturation temperature setpoint is automatically determined. PPMw is independent of test pressure and test temperature.

Ramp Time

The Ramp Time is the desired amount of time that the 3900 should take to transition from one profile test point to another. Setting a ramp time of zero instructs the 3900 to make the transition as quickly as possible. Zero is the setting used for most applications.

RTD

Resistance thermal devices or resistance thermal detectors (RTD) are temperature sensors that are based on the predictable change in electrical resistance of a material (usually platinum) with changes in temperature. Platinum RTDs are often called platinum resistance thermometers (PRT).

RTS

Request to Send (RTS) is a serial signal sent to verify that the other device is ready for data.

Sample Size

The sample size is the time span that defines the sample points to use in the Standard Deviation calculation. The sample size should always be carefully considered based on the data storage interval. Too small of a sample size in relation to the data storage interval will result in a small number of points used to calculate the Standard Deviation.

Saturation Pressure Control Mode

The Saturation Pressure control mode, Ps, is controlled at a constant value independent of any other pressure, temperature, or humidity value. While saturation pressure is held constant, all humidity parameters may vary.

Saturation Vapor Pressure (SVP)

Saturation Vapor Pressure is a function of temperature. The function can best be described by a lab setup. Imagine a chamber whose temperature T can be controlled. The chamber is partially filled with water. Initially, the remaining space is a vacuum. The pressure P of the space over the water can be measured. At a fixed temperature, water molecules will leave the water and enter the space above at a fixed rate. As water molecules accumulate over the liquid water, the pressure there will increase, and molecules will re-enter the liquid at an increasing rate. Finally, water molecules will be entering and leaving the liquid at the same rate, giving equilibrium and a constant pressure P over the water. The equilibrium pressure P is the Saturation Vapor Pressure at temperature T.

Soak Time

The Soak Time is the desired amount of time to generate at a particular profile point. The soak time required depends on the application, but should be a significant amount of time based upon the humidity measurement devices being calibrated.

Specific Humidity

Specific Humidity is a ratio of the weight of the water vapor to the total weight of the humidified gas. Specific Humidity is independent of test chamber temperature.

Standard Deviation

Assurance condition Standard Deviation is a statistic used to measure the variation in the actual data and can be thought of as how spread out or stable the data is. ControLog calculates Standard Deviation from the device data tabs along the selected Sample Size. Once the actual standard deviation is less than the defined limit the standard deviation portion of the condition is considered meet.

Supply Pressure

Supply pressure is the regulated pressure reading of the gas supply. The supply pressure value is in gauge.

SVP@Td

Saturation Vapor Pressure (SVP) computed at the Dew/Frost Point Temperature.

SVP@Ts

Saturation Vapor Pressure (SVP) computed at the Saturation Temperature.

SVP@Tt

Saturation Vapor Pressure (SVP) computed at the Test Temperature.

Thermistor

A thermistor is a type of resistor whose resistance varies with temperature. Thermistors are generally made of a ceramic or polymer type material.

Tolerance

An assurance condition Tolerance is the allowable variation between the setpoint and the actual. This is best thought as a window based on a minimum and maximum value, the minimum being the setpoint minus the tolerance and the maximum being the setpoint plus the tolerance. Once the actual value is within the window the tolerance portion of the condition is considered met.

Vapor Mole Fraction

Vapor Mole Fraction is the mole fraction of water vapor in a sample.

Wet Bulb Temperature

Wet Bulb temperature is used in wet bulb / dry bulb aspirated Psychrometry, and is the temperature measured by a temperature probe whose tip is coated with water (typically by being covered with a wet sock). When aspirated at a constant air velocity, the wet bulb will cool due to evaporation of the water from the probe. The cool temperature, to which it equilibrates, is used in the calculation of humidity parameters.

Index

%

%RH 53

/

/ 40

3

3900 Interface Console 67

3900 Parameter Tab 49

3900 Reported Errors 66

A

About ControLog 1, 42

Absolute Humidity 53

Advancing to the Next Phase 132

Advancing to the Next Point 132

Analog Connection 174

Analog Connection Example 182

Assurance 114

Assurance Phase 129

Auto Profile Controls 72

Auto Profile Controls and Status 72

Auto Profile Status 72

Auto Profiling 113

Auto Scale 84

Auto Scale Graph 28

C

Calculated Humidity Parameters 53

Changing Setpoints 50

Clear Data 33

Clear Status Log 18

Close All Data Tabs 18

Close All Graphs 28

Close Connection 38

Close Data 33

Close Graph 28

Close Profile 26

Connect to 3900 38

Connection and Run Status 70

Connections Menu 38

Control Mode 114

Control Parameters 50

ControLog Minimum System Requirements 4

Copy All Data 33

Copyright 3

Crashes and Data Recovery 12

Create Data Summary 31

Creating a Data Summary 99

Creating a New Graph 76

Creating a New Profile 114

Current Point Time 73

Customizing a Graph 80

D

Data and Graph Tab Group 59

Data Backup 5

Data Menu 30

Data Storage Interval 33

Data Summary Tabs 98

Data Tabs 61

Density Unit 36

Device Data Tabs 94

Device Interface Console 68

Device Parameter Tabs 57

Device Parameters 58

Dew Point 53

Dew Point Control Mode 21

Disconnected 70

Dry Air Density 53

Dry Air Mole Fraction 53

E

Enthalpy 53

Enthalpy Unit 37

Example 1 130

Example 2 130

Example 3 130

Exit ControLog 20

F

F@Td.Pt 53

F@Ts.Ps 53

F@Tt.Pt 53

File Data Tabs 95

File Menu 17

Flow Unit 36

Frost Point 53

Frost Point Control Mode 21

G

- Generate Mode 23
- Generating 71
- GPIO Connection 190
- GPIO Connection Example 203
- Grains/lb 53
- Graph Menu 27
- Graph Properties 29, 85
- Graph Tabs 61

H

- Help Menu 42
- Help Topics 42
- Hold Profile Point 26
- Holding the Profile 131
- How to Record a Data Point 108

I

- Import Uncertainty Solution 18
- Information Tab Group 62
- Installing ControLog 5
- Interface Connections 6

L

- License 7
- License Agreement 2

M

- Manual Connection 212
- Manual Connection Example 216
- Manual Override of Profile 131
- Mass Flow Rate 114
- Menu Bar 16
- Mixing Ratio by Volume 53
- Mixing Ratio by Weight 53
- Mode Menu 21
- Moist Air Density 53

N

- New Connection 38
- New Event 71
- New Graph 27
- New Profile 26
- Next Profile Phase 26
- Next Profile Point 26
- Notice 1

O

- Open 3900 Interface Console 39
- Open Connection 38

- Open Data File 17, 31
- Open Device Interface Console 40
- Open Profile 26
- Opening a Profile 125

P

- Pan 80
- Parameters Tab Group 43
- Percent by Volume 53
- Percent by Weight 53
- Percent Relative Humidity Control Mode 21
- PPMv 53
- PPMv Control Mode 21
- PPMw 53
- PPMw Control Mode 22
- Pressure Unit 36
- Print Data 33
- Print Graph 28
- Print Profile 26
- Print Status Log 18
- Printing a Graph 90
- Printing a Profile 126
- Product Key 7
- Profile 113
- Profile Menu 25
- Profile Tab 61
- Profiling 113
- Purge Mode 23
- Purging 71

R

- Ramp Phase 129
- Ramp Time 114
- Record Data Point 31
- Recorded Data Points Tab 107
- Run Auto Profile 26
- Run Menu 23
- Running an Auto Profile 127

S

- Sample Size 114
- Saturation Pressure Control Mode 21
- Save All Data 17
- Save Data 33
- Save Graph 28
- Save Profile 26
- Save Status Log 18
- Saving a Graph 88
- Saving a Profile 124
- Scale 84
- Scale Graph 28
- Scale Graph's X Axis 29
- Scale Graph's Y Axis 29

- Scale X Axis 84
- Scale Y Axis 84
- Serial Connection 134
- Serial Connection Example 1 147
- Serial Connection Example 2 158
- Settings 31, 39
- Shutdown 24, 71
- Soak Phase 130
- Soak Time 114
- Specific Humidity 53
- Standard Deviation 114
- Starting ControLog 7
- Status Bar 70
- Status Log 66
- Stop Auto Profile 26
- Stopping the Auto Profile 132
- Supply Pressure 53
- SVP@Td 53
- SVP@Ts 53
- SVP@Tt 53
- System Requirements 4
- System timing 66

T

- Take Point Ctrl-P 32
- Technical Support 1
- Temperature Unit 35
- Test Pressure 114
- Test Temperature 114
- Thunder Scientific on the Web 42
- Tolerance 114
- Trademark 3

U

- Understanding Profile Phases 129
- Units Menu 35

V

- Vapor Mole Fraction 53

W

- Warranty 3
- Wet Bulb Temperature 53

Z

- Zoom 81
- Zoom Graph 28
- Zoom Graph's X Axis 28, 82
- Zoom Graph's Y Axis 28, 83