


**THUNDER SCIENTIFIC
CORPORATION**

2500 *ControLog* Software

FEATURES

- Powerful Graphing Capability Creates a Visual Picture of the Data
- Auto Profiling Feature Automates Humidity Generation
- Data Stored in a Familiar Spreadsheet Type Layout
- Customizable ASCII Interface Support for RS-232, GPIB and Analog Devices
- Uncertainty Calculated in Real-Time by HumiCalc with Uncertainty

Automation Software for the Model 2500 Humidity Generator



“Power by HumiCalc with Uncertainty”[®]

623 Wyoming Blvd. SE Albuquerque, NM 87123-3198
Copyright © 2010 - 2012, All Rights Reserved